



3M™ Petrifilm™ Plates

Guide to Dilution Preparation

Many foods require a dilution step before plating and the following is a guide to follow when preparing dilutions. Please refer to the 3M™ Petrifilm™ Plate Application Guide. Use with Dairy and Juice Products, for specific dilution recommendations.

Dilution	Sample size	Diluent
1:5	10.00 g	40 mL
	or	
	22.50 g	90 mL
1:10	1.00 g	9 mL
	or	
	10.00 g	90 mL
1:20	4.74 g	90 mL
	or	
	5.21 g	99 mL

Diluent Volume Calculations:

Volume of Buffer Needed = (Weight of Sample x Dilution) – Weight of Sample

Example: You have 7.5 mL of sample and want to prepare 1:10 dilution.


$$\left(\begin{array}{c} \text{7.5 mL} \\ \text{Weight of Sample} \end{array} \times \begin{array}{c} \text{10} \\ \text{Dilution} \end{array} \right) - \begin{array}{c} \text{7.5 mL} \\ \text{Weight of Sample} \end{array} = \begin{array}{c} \text{67.5 mL} \\ \text{Volume of Buffer Needed} \end{array}$$

Procedure for Determining Counts:

Use a multiplication factor to convert the plate count to the number of colony forming units (CFU) of bacteria per gram of sample. The multiplication factor is determined by dividing the dilution used by the volume plated.

Example 1: Single plate

1 mL plated of 1:10 dilution





count = 3

$$\frac{10 \text{ (1:10 dilution)}}{1 \text{ mL plated}} = 10 \text{ (dilution factor)}$$

$$3 \text{ (count)} \times 10 \text{ (dilution factor)} \div 1 \text{ (mL plated)} = 30 \text{ (CFU/g reported)}$$

Example 2: Multiple plates in order to achieve higher sensitivity

1 mL each plated of 1:10 dilution on two plates equals 2 mL total

total count = 3

$$\frac{10 \text{ (1:10 dilution)}}{2 \text{ mL plated}} = 5 \text{ (dilution factor)}$$

$$3 \text{ (count)} \times 10 \text{ (dilution factor)} \div 2 \text{ (mL plated)} = 15 \text{ (CFU/g reported)}$$

Sensitivity of 3M™ Petrifilm™ Plates:

If no colonies appear on the plate, do not report the count as zero. Report in CFU/g as shown at right. Most 3M Petrifilm Plates require a 1 mL volume of sample. 3M™ Petrifilm™ High-Sensitivity Coliform Count Plate requires a 5 mL volume of sample.

Dilution	3M Petrifilm Plate	CFU/g
1 : 10	1 mL plate	< 10.0
	5 mL plate	< 2.0
1 : 5	1 mL plate	< 5.0
	5 mL plate	< 1.0
1 : 2	1 mL plate	< 2.0
	5 mL plate	< 0.4
Undiluted	1 mL plate	< 1.0
	5 mL plate	< 0.0

For greater sensitivity, a lower dilution and /or more plates may be used.

For detailed WARNING, CAUTIONS, DISCLAIMER OF WARRANTIES / LIMITED REMEDY, LIMITATION OF 3M LIABILITY, STORAGE AND DISPOSAL information, and INSTRUCTIONS FOR USE see product's instructions online at 3M.com/foodsafety/Petrifilm



3M Food Safety
3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA

Phone 1-800-328-6553
Web 3M.com/foodsafety

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-800-364-3577

3M and Petrifilm are trademarks of 3M.
Printed in U.S.A. © 3M 2019.
70-2009-1806-1 (Rev-0219)