Nexcare™

ColdHot Therapy Pack Flexible

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en Reusable cold and hot gel pack

The 3M™ Nexcare™ ColdHot Therapy Pack Flexible can be used cold or hot and is a convenient solution for those who prefer drug-free therapy. A protective cover for the gel pack is included, allowing for easy and comfortable application. The gel is non-toxic and safe for the whole family, even accidental ingestion of a small amount of gel will not be harmful.

Cold therapy

Cold helps to bring relief in swelling, inflammation and pain.

When to use ColdHot Flexible pack

- Sprained ankle
- Strained muscles in leg/arm after sport activities
- Contusion swelling of knee
- Injured ligaments

Preparation - cold

Refrigerator or freezer
Keep the gel pack in the refrigerator or freezer for minimum 2h before use. The gel pack will be always ready to use if it is stored constantly in the refrigerator or freezer.

Always insert gel pack in the supplied cover or clean towel, when applying to the skin. Apply for no longer than 20 minutes. Then wait another 20 minutes before reapplying the cold pack again.

Hot therapy

Heat helps to reduce muscle spasm and to decrease pain, soothes and relaxes.

When to use ColdHot Flexible pack

- Repetitive strain injuries e.g. tennis, electronic games
- Knee or general leg muscle stiffness
- Wrist arthritis as an adjunct therapy for rheumatoid arthritis (without acute inflammation)

Gel pack care instructions

The ColdHot Flexible gel pack can be washed by hand in a mild detergent. Do not use alcohol or other solvents to clean the gel pack. Store in a dark place.

Disposal
Dispose of product in regular household waste.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Date of issue of the instruction of use 2019-10

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© 2022 3M. All rights reserved.
3M and Nexcare are trademarks of 3M.