

Tips for a trouble-free smile.

Forsus™ Class II Correctors



Making beautiful smiles easier.

We're harnessing the power of 3M Science to improve oral health – one smile at a time. Your Forsus™ Appliance will help efficiently move your teeth to your doctor's recommended positions, and can reduce or eliminate the need for headgear or extractions of permanent teeth. This appliance will help you reach your goal of a new smile, because it is always working for you.

Designed for maximum comfort and durability, the Forsus appliance can withstand normal force and constant use, but without appropriate care it can be damaged or broken. To get the most out of your Forsus appliance, follow these tips to ensure the easiest route to a beautiful smile.

It's important to keep all your scheduled appointments to monitor progress.



Contact your orthodontist if you have any questions.

Important Consumer Information – Forsus™ Class II Correctors:

Warning: Forsus Correctors contain nickel and chromium. A small number of people are known to be allergic to these metals. Inform your orthodontist if you have known allergies. Consult a physician immediately if an allergic reaction occurs or is suspected.

Important reminders:

Reviewed with Patient

- At 3M, we recognize the need to help your doctor create beautiful smiles with minimal impact on patients' lives during treatment. The Forsus appliance is carefully constructed to enable normal eating and speaking, however, opening your mouth extremely wide while wearing your Forsus appliance can damage it, and would require you to return to your orthodontist's office.
- The appliance is designed to stay clear of your teeth when you open or close your mouth. But if it does get between your teeth, avoid biting it. This could injure your mouth, or deform or break the appliance. Always call your orthodontist promptly if this problem occurs.
- As with any orthodontic appliance worn inside the mouth, some discomfort or irritation may occur while you adjust to the new appliance. Contact your orthodontist if a sore develops in your mouth, or discomfort and irritation continues.
- Ultimately, protecting your Forsus appliance from the foods you eat is up to you. Avoid hard or sticky foods, and reduce your food into smaller pieces. We recommend always brushing your teeth after each meal, and avoid dislodging or separating the appliance when brushing. To monitor your progress, it's important to keep all your scheduled appointments.