Wear it Right

Wearing Your Respirator

1



Place the respirator over your nose and mouth. Be sure the metal nose clip is on top. With models 8210 or 07048, pre-stretch the straps before wearing.

2



Pull the top strap over your head until it rests on the crown of your head above your ears.

3



Pull the bottom strap over your head until it rests just below your ears.

4



Using both hands starting at the top, mold the metal nose clip around your nose to achieve a secure seal.



Filtering Facepiece Non-Valved Respirator



Filtering Facepiece Valved Respirator

5

Check the seal of your filtering facepiece respirator each time you don the respirator.



Positive Pressure
User Seal Check
For Non-Valved
Respirators Place both
hands completely
over the respirator and
exhale. The respirator
should bulge slightly.
If air leaks between

the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.



Negative Pressure
User Seal Check
For Valved Respirators
Place both hands
over the respirator
and inhale sharply.
The respirator should
collapse slightly. If
air leaks between

the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.

