Wear it Right

Wearing Your Respirator



Place the respirator over your nose and mouth. Be sure the metal nose clip is on top. With models 8210 or 07048, pre-stretch the straps before wearing.



Pull the top strap over your head until it rests on the crown of your head above your ears.



Pull the bottom strap over your head until it rests just below your ears.



Using both hands starting at the top, mold the metal nose clip around your nose to achieve a secure seal.



Filtering Facepiece Non-Valved Respirator



Filtering Facepiece Valved Respirator

Check the seal of your filtering facepiece respirator each time you don the respirator.



Non-Valved Filtering Facepiece Respirator

Place both hands completely over the respirator. Be careful not to disturb the position of the respirator. Inhale and exhale sharply. If air leaks around nose, readjust the nosepiece. If air

leaks at the respirator edges, work the straps back along the sides of your head. If you **CANNOT** achieve a proper seal, DO NOT enter the contaminated area. See your supervisor.



Valved Filtering Facepiece Respirator

Place both hands over the respirator and inhale sharply. Be careful not to disturb the respirator. If air leaks around nose, readjust the nosepiece. If air leaks at the respirator edges, work the straps back

along the sides of your head. If you CANNOT achieve a proper seal, DO NOT enter the contaminated area. See your supervisor.

AWARNING