

# 3M Purification Inc.

## Zeta Plus™ Installation Qualification Test Procedure

### Purpose

The purpose of the Installation Qualification (IQ) test is to determine if Zeta Plus™ cartridges have been properly installed and sealed in a filter housing. This procedure is specifically intended for use with Zeta Plus cartridge grades 50, 60, 70, and 90. The procedure can be used for 8, 12, and 16 inch nominal diameter cartridge sizes. This document defines the general procedure and the necessary equipment required to determine if the Zeta Plus cartridge(s) have been properly installed in the filter housing. This procedure does not attempt to establish or define any correlation between the measured drop in pressure during the test and the ability of the filter to reduce or remove contaminants from the fluid being filtered.

### General

#### Equipment and Materials Required

- Stop watch
- 0 - 5 psi (0" - 137" H<sub>2</sub>O) Pressure Measuring Device with + 0.25% FS accuracy
- Programmable Automatic Test Instrument (Optional)
- Clean, Dry Regulated Gas Pressure Source Reagents
- Water

#### Applicable Documents

- Safety glasses with side shields must be worn at all times.
- Wear appropriate safety clothing and equipment as recommended in the MSDS for the fluid being filtered and operating conditions.
- Do not exceed the maximum pressure ratings for the test equipment and filter housings.
- Always vent the filter housing before opening.

### Procedure

#### Hardware Qualification Test (Ref. Figure 1)

1. Close and seal filter housing.
2. Close valves A, D, E, and F open valve C.
3. Attach a manual test assembly or an appropriate electronic test instrument to vent valve C.

**NOTE:** A manual test assembly consisting of an air regulator (capable of adjusting the out put pressure in increments of 0.05 psi) equipped with a 0 - 5 psi pressure measuring device (accuracy + .25% F.S.) followed by a shut-off valve (ball or needle type) (Ref. Figure 1A.)

4. Enter the following test parameters when a programmable electronic test instrument is used (table below):

5. Pressurize the filter housing to 2.9 psi (200 mbar).
6. When a manual test assembly is used, after three minutes (180 seconds), close valve B and start the stop watch.
7. After three minutes (180 seconds), record the pressure on the test gauge.
8. If the pressure is less than 2.83 psi (196 mbar) the system is leaking. Identify the source of the leak and stop it.
9. Repeat steps 1 through 8 until the pressure recorded at the end of the three minute test period is greater than 2.83 psi, then proceed to Cartridge Qualification Test.


Stabilization time	180 seconds (3 minutes)
Pressure hold test time	180 seconds (3 minutes)
Test Pressure	2.9 psi (200 mbar)
Maximum allowable pressure decay	0.07 psi (5 mbar)

#### Cartridge Qualification Test (Ref. Figure 1)

1. If attached, remove the manual test assembly or electronic test instrument from valve C.
2. Open the housing and install and seal the Zeta Plus filter cartridge(s) in the filter housing.
3. Close and seal filter housing.

4. Close valves A, D, E, and F; open valve C.
5. Connect the inlet of the filter housing to a filtered water source. (Note: The water should be filtered to a level equal to or better than the micron rating of the filter being tested.)
6. Fill the housing with water by slowly opening inlet valve A.
7. When fluid emerges from valve C, close the valve while simultaneously opening either outlet valve E or drain valve D.
8. Flow a volume of water through the filter that is equivalent to 54 L/m<sup>2</sup> of filter medium at a flow rate equal to or greater than the anticipated process flow rate.
9. Enter the following test parameters when a programmable electronic test instrument is used (table below):

Draining Time	300 seconds (5 minutes)
Stabilization time	180 seconds (3 minutes)
Pressure hold test time	180 seconds (3 minutes)
Low Drain Pressure	0.87 psi (60 mbar)
High Drain Pressure	1.16 psi (80 mbar)
Test Pressure	1.45 psi (100 mbar)
Maximum allowable pressure decay	0.73 psi (50 mbar)



**Figure 1. - Test Assembly**

10. Close valve A and open drain valve F.
11. Pressurize the filter housing to between 0.87 to 1.16 psi (60 - 80 mbar) using either the electronic test instrument or the manual test assembly.
12. When a manual test assembly is used, after five minutes (300 seconds), inspect drain lines D (or E) and F to determine if bulk water flow has stopped. If bulk water flow has not stopped, allow the housing to continue to drain until the water flow is reduced to a trickle. When bulk water flow stops close only drain valve F and proceed to Step 14. Valves D (or E) should remain open during the test.
13. When an electronic test instrument is used, after five minutes (300 seconds), inspect drain lines D (or E) and F to determine if bulk water flow has stopped. If bulk water flow has not stopped, discontinue the test and reprogram the instrument for a longer drain period. Prior to the start of the stabilization period, the water flow rate through valves D (or E) and F should be reduced to a trickle (Note: Do not reprogram instrument for a higher drain pressure setting.) When bulk water flow stops, immediately close only drain valve F. Valves D (or E) should remain open during the test.
14. Allow the electronic test instrument to complete the pre programmed test cycle and print out the test result.
15. When a manual test assembly is used, carefully increase the pressure to 1.45 psi and then close valve B and start the stop watch. (Note: Do not to exceed 1.45 psi.)
16. After three minutes (180 seconds), record the pressure on the test gauge.
17. If the pressure reading is greater than 0.72 psi, the installation has passed the IQ test.
18. If the pressure reading is less than 0.72 psi, either the cartridge(s) is damaged, improperly installed, or both. Check the installation by completing the following:
  - A. Disconnect the test assembly and relieve the pressure on the filter housing by opening vent valve C.
  - B. Open the housing and inspect the cartridge(s) for damage, replace if damaged.
  - C. Inspect seal nut and seal nut O-rings for signs of damage, replace if damaged.
  - D. Inspect cartridge(s) flat gaskets or O-ring for signs of damage, replace if damaged.
19. Repeat Steps 1 through 18 as required.

### Important Notice

The information described in this literature is accurate to the best of our knowledge. A variety of factors, however, can affect the performance of the Product(s) in a particular application, some of which are uniquely within your knowledge and control. **INFORMATION IS SUPPLIED UPON THE CONDITION THAT THE PERSONS RECEIVING THE SAME WILL MAKE THEIR OWN DETERMINATION AS TO ITS SUITABILITY FOR THEIR USE. IN NO EVENT WILL 3M PURIFICATION INC. BE RESPONSIBLE FOR DAMAGES OF ANY NATURE WHATSOEVER RESULTING FROM THE USE OF OR RELIANCE UPON INFORMATION.**

It is your responsibility to determine if additional testing or information is required and if this product is fit for a particular purpose and suitable in your specific application.

**3M PURIFICATION INC. MAKES NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER NATURE HEREUNDER WITH RESPECT TO INFORMATION OR THE PRODUCT TO WHICH INFORMATION REFERS.**

### Limitation of Liability

3M Purification Inc. will not be liable for any loss or damage arising from the use of the Product(s), whether direct, indirect, special, incidental, or consequential, regardless of the legal theory asserted, including warranty, contract, negligence or strict liability. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.



3M Purification Inc.  
400 Research Parkway  
Meriden, CT 06450, U.S.A.  
Tel (800) 243-6894  
(203) 237-5541  
Fax (203) 238-8977  
www.3Mpurification.com

3M is a trademark of 3M Company.  
Zeta Plus is a trademark of 3M Company used under license.  
© 2012 3M Company. All rights reserved.  
Please recycle. Printed in U.S.A.  
70-0201-8802-8 REV 0212b