

# Zeta Plus™ EXT Series Products

Featuring Dual Zone ZB Filter Media

The media used in Zeta Plus™ EXT filter products *extends* depth filtration capacity and enhances the protection of downstream membranes thereby *extending* filter life. Zeta Plus EXT series is a family of advanced dual zone depth filters designed to provide optimal clarification of bioprocess, biological and pharmaceutical fluids. Zeta Plus EXT Series Products contain two distinct layers, or “zones,” of filter media with the upstream zone more open than the downstream zone. This structure enhances the contaminant holding capacity of the filter media since larger particles are trapped in the upper zone of the filter media and smaller particles are trapped in the lower zone. This media structure helps reduce premature plugging and extends service life. The two filter zones can be independently selected and combined to optimize performance.

- Increased throughput with high turbidimetric efficiency
- Lower processing costs
- Smaller filter assemblies
- Consolidation of multiple filter stages
- Enhanced protection of downstream membrane filters

## Applications

Zeta Plus EXT Series Filters with ZB Media are optimized for clarification of mammalian cell harvest fluids and for protection of downstream sterilizing grade membrane filters. Their strong positive charge provides enhanced reduction of negatively charged contaminants compared to a mechanical filtration mechanism alone. In solution cells, cellular debris, endotoxins, and nucleic acids (DNA, RNA) carry a negative charge. 3M's application of positive zeta potential to ZB filter media creates a filter matrix that enhances the capture of negatively charged contaminants for maximum filtration performance. The strong positive charge over a broad pH range makes these filters ideal for cell broth clarification and for protection of downstream membrane filters.

## Regulatory Support

Zeta Plus™ EXT Series products with ZB Media are supported with a Regulatory Support File (3M literature 70-0203-0883-2) and Certificate of Quality to ease validation and regulatory submissions. The products are also documented in a Drug Master File (DMF) on record with the US Food and Drug Administration (FDA).



## Applications

Clarification of mammalian cell culture process fluids (cell separation)

Protection of downstream processes including membrane filters and chromatography columns

Clarification of bacteria and yeast cell lysates

Endotoxin and nucleic acid reduction

Clarification of colloidal and haze forming contaminants



Zeta Plus™ BC25 Capsule filters are ideal for laboratory evaluation and media selection

## Cartridge Construction

Individual cells of Zeta Plus EXT Series with ZB Media are assembled together with polypropylene separators under predetermined compression and unitized by three, 316 stainless steel bands.\* Each cell is constructed using thermoplastic elastomer or polypropylene molded edge seals and polypropylene separators for high performance. The standard gasket is composed of silicone, with other gasket materials available upon request. Filter cartridges are available in 8", 12", and 16" diameters, with surface area ranging from 0.16 m<sup>2</sup> (1.8 ft<sup>2</sup>) to 2.8 m<sup>2</sup> (29.8 ft<sup>2</sup>) per cartridge.

## Capsule Construction

Zeta Plus EXT series products are available in autoclavable and single-use BC25 and Zeta Plus Encapsulated System capsules. With 25 cm<sup>2</sup> (3.9 in<sup>2</sup>) of filtration area, the BC filter capsule is ideal for laboratory feasibility and initial product sizing evaluations. The Zeta Plus Encapsulated System includes scale-up capsule filters with effective filtration areas at 170 cm<sup>2</sup> (0.18 ft<sup>2</sup>), 340 cm<sup>2</sup> (0.37 ft<sup>2</sup>), 1020 cm<sup>2</sup> (1.10 ft<sup>2</sup>); and large capsule filters with 0.23 m<sup>2</sup> (2.4 ft<sup>2</sup>), and 1.6 m<sup>2</sup> (17.1 ft<sup>2</sup>). They are ideal for both scale-up studies as well as large scale production. For more information on the Zeta Plus Encapsulated System, please refer to 3M literature 70-0202-3967-2.

## Grade Selection

The Zeta Plus EXT Series Filters with ZB Media are available in six grades. Operating conditions and the fluid being filtered impact retention performance. 3M Purification recommends small-scale pilot runs to confirm grade selection prior to full scale-up. 3M Purification's Technical Services staff can assist in grade selection as well as perform on-site evaluations.




Table 2: Effective Surface Area

| Cartridge Configuration                        | Surface Area                                 |
|--|--|
| BC25 Capsule                                   | 25 cm <sup>2</sup> (3.9 in <sup>2</sup> )    |
| E0170FSA Scale-up Capsule                      | 170 cm <sup>2</sup> (0.18 ft <sup>2</sup> )  |
| E0340FSA Scale-up Capsule                      | 340 cm <sup>2</sup> (0.37 ft <sup>2</sup> )  |
| E1020FSA Scale-up Capsule                      | 1020 cm <sup>2</sup> (1.10 ft <sup>2</sup> ) |
| Z08 E05 (8 diameter cartridge plug-in, 5-cell) | 0.16 m <sup>2</sup> (1.8 ft <sup>2</sup> )   |
| Z08 E06 (8 diameter cartridge DOE, 6-cell)     | 0.20 m <sup>2</sup> (2.1 ft <sup>2</sup> )   |
| Z08 E07 (8 diameter cartridge DOE, 7-cell)     | 0.23 m <sup>2</sup> (2.4 ft <sup>2</sup> )   |
| Z16 E01 (16 diameter cartridge, 1-cell)        | 0.23 m <sup>2</sup> (2.4 ft <sup>2</sup> )   |
| Z16 E02 (16 diameter cartridge, 2-cell)        | 0.46 m <sup>2</sup> (5.0 ft <sup>2</sup> )   |
| Z16 E08 (16 diameter cartridge, 8-cell)        | 1.8 m <sup>2</sup> (19.8 ft <sup>2</sup> )   |
| Z16 E12 (16 diameter cartridge, 12-cell)*      | 2.8 m <sup>2</sup> (29.8 ft <sup>2</sup> )   |

\*Z16E12 (16" diameter cartridge, 12-cell) wet weight (29 kg (65 lb)) may require lifting assistance.

Table 3: Recommended Operating Parameters

|                               |                                |  |
|-------------------------------|--------------------------------|--|
| Maximum Operating Pressure    | BC25 Capsules                  | 2.8 bar (40 psig) maximum inlet pressure, 2.4 bar (35 psid) maximum capsule pressure drop. |
|                               | Scale-Up Capsules              |  |
|                               | Standard Zeta Plus™ Cartridges | 2.4 bar (35 psid) maximum cartridge pressure drop.   |
| Maximum Operating Temperature | BC25 Capsules                  | 40 °C (104 °F)   |
|                               | Scale-Up Capsules              |  |
|                               | Zeta Plus™ Cartridges          | 82 °C (180 °F)   |
| Recommended Pre-use Rinse     |                                | 54 L/m <sup>2</sup> (1.33 gal/ft <sup>2</sup> )  |
| Sterilization Parameters      | BC25                           | Autoclave 30 min. at 121 °C (250 °F) (1 cycle).  |
|                               | Scale-Up Capsules              | Autoclave 30 min @ 126 °C (259 °F)   |
|                               | Standard Zeta Plus™ Cartridges | Autoclave or <i>In situ</i> steam sterilize 30 minutes @ 126 °C (259 °F) (3 cycles)        |

## Sanitary Housings

3M Purification provides a wide array of standard and specially designed sanitary filter housings to accommodate Zeta Plus™ EXT series filter cartridges. All housings are designed to the requirements of the pharmaceutical and bioprocessing industries. The housings feature 316L polished surfaces and easy-to-clean components, along with accessories such as CIP spray-balls and the innovative Zeta Plus Lifting Device. Custom design skid systems incorporating Zeta Plus and membrane cartridge filter housings, together with process piping, flow control and complete validation documentation are also available.

Zeta Plus EXT series cartridges can be retrofitted into existing standard Zeta Plus housing offerings. For more information, please call your local 3M Purification Sales Representative or Distributor.

Figure 1: Nominal Retention Ratings





Figure 2: Cartridge Flow vs. Differential Pressure



# Zeta Plus™ EXT Series Products Featuring Dual Zone ZB Filter Media Ordering Guide

## 8" Diameter Cartridges

| Catalog Number | Configuration | Number of Cells | Cartridge Construction                            | Gasket Material | Grade  |
|----------------|---------------|-----------------|---|-----------------|--|
| Z08            | E - EXT       | 05 - 5 cell     | P - Polypropylene Plug In                         | A - Silicone    | 60ZB05A 120ZB05A<br>90ZB05A 120ZB08A<br>90ZB08A 120ZB10A |
|                |               | 07 - 7 cell     | A - Stainless Steel Bands<br>B - Hastelloy® Bands |                 |  |

Part Number Example: Z08E05PA60ZB05A. All 8" cartridges are packaged two units per box.

## BC Series for Feasibility & Sizing Trials

| Catalog Number                       | Grade   |
|--------------------------------------|---|
| BC0025L (Luer)<br>BC0025S (Sanitary) | 60ZB05A<br>90ZB05A<br>90ZB08A<br>120ZB05A<br>120ZB08A<br>120ZB10A |

Part Number Example:  
BC0025L60ZB05A.

## 12" Diameter Cartridges

| Catalog Number | Configuration | Number of Cells | Cartridge Construction                            | Gasket Material | Grade  |
|----------------|---------------|-----------------|---|-----------------|--|
| Z12            | E - EXT       | 11 - 11 cell    | A - Stainless Steel Bands<br>B - Hastelloy® Bands | A - Silicone    | 60ZB05A 120ZB05A<br>90ZB05A 120ZB08A<br>90ZB08A 120ZB10A |

Part Number Example: Z12E11AA60ZB05A. All 12" cartridges are packaged two units per box.

## 16" Diameter Cartridges

| Catalog Number | Configuration | Number of Cells | Cartridge Construction                            | Gasket Material | Grade  |
|----------------|---------------|-----------------|---|-----------------|--|
| Z16            | E - EXT       | 01 - 1 cell     | A - Stainless Steel Bands<br>B - Hastelloy® Bands | A - Silicon     | 60ZB05A 120ZB05A<br>90ZB05A 120ZB08A<br>90ZB08A<br>120ZB10A* |
|                |               | 02 - 2 cell     |   |                 |  |
|                |               | 08 - 8 cell     |   |                 |  |
|                |               | 12 - 12 cell    |   |                 |  |

Part Number Example: Z16E08AA120ZB08A.  
8-cell cartridges are packaged one unit per box.  
\*Not available in 12 cell.

### Important Notice

The information described in this literature is accurate to the best of our knowledge. A variety of factors, however, can affect the performance of the Product(s) in a particular application, some of which are uniquely within your knowledge and control. **INFORMATION IS SUPPLIED UPON THE CONDITION THAT THE PERSONS RECEIVING THE SAME WILL MAKE THEIR OWN DETERMINATION AS TO ITS SUITABILITY FOR THEIR USE. IN NO EVENT WILL 3M PURIFICATION INC. BE RESPONSIBLE FOR DAMAGES OF ANY NATURE WHATSOEVER RESULTING FROM THE USE OF OR RELIANCE UPON INFORMATION.**

It is your responsibility to determine if additional testing or information is required and if this product is fit for a particular purpose and suitable in your specific application.

**3M PURIFICATION INC. MAKES NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED INCLUDING WITHOUT LIMITATION ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER NATURE HEREUNDER WITH RESPECT TO INFORMATION OR THE PRODUCT TO WHICH INFORMATION REFERS.**

### Limitation of Liability

3M Purification Inc. will not be liable for any loss or damage arising from the use of the Product(s), whether direct, indirect, special, incidental, or consequential, regardless of the legal theory asserted, including warranty, contract, negligence or strict liability. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.



Your Local Distributor:

### 3M Purification Inc.

400 Research Parkway  
Meriden, CT 06450  
U.S.A.  
(800) 243-6894  
(203) 237-5541  
Fax (203) 630-4530  
www.3Mpurification.com

3M and Zeta Plus are trademarks of 3M Company.  
All other trademarks are the property of their respective owners.  
© 2014 3M Company. All rights reserved.  
70-0201-8862-2  
REV 0114