Intended Use and Product Description

The information provided in this insert applies to the 3M™ Red Dot™ Resting EKG Electrodes 2330, 2360 and 9650. The 3M Red Dot resting electrodes are disposable with a silver/silver chloride sensing element, are designed for a patient undergoing a short-term EKG diagnostic procedure while resting and are radiolucent.

Caution

- U.S.A. Federal Law restricts this device to sale by or on the order of the physician.
- The 3M™ Red Dot™ Resting EKG Electrodes 2330, 2360 and 9650 were determined to be MR-Conditional via non-clinical testing. For information on the use of a 3M Red Dot resting electrode during an MRI procedure please see below.

Precautions

- To prevent dryout, electrodes should be used immediately after opening the package, except the 2360 resting tab electrode which should be used within 7 days of opening the pouch. The unused electrodes should remain stored in the opened pouch.
- Do not apply electrodes over pantyhose.
- For proper skin management and to minimize skin irritation:
 - Avoid placing an electrode on an irritated skin site.
 - Do not abrade a skin site more than one time.
 - Avoid placing electrodes on skin still wet from an alcohol wipe (dry thoroughly).
 - Do not use for long-term monitoring.

Instructions for Use**1. Before Electrode Application:**

- Select and prepare an electrode site according to your health care facility's protocol for resting EKG monitoring.
- The electrode site should be clean, dry and free of skin oil or lotion before application of the electrode to allow optimal electrode adhesion to skin. Application sites with heavy or curly hair that cannot be parted should be clipped.
- For best trace quality, skin abrasion should be part of the skin preparation. Use 3M™ Red Dot™ Trace Prep as part of the skin preparation.

2. Electrode Application:

- Remove electrode from the liner. Apply to skin at standard placement sites.
- With tabs pointing toward the patient's waist, firmly press electrodes to assure uniform contact.
- For difficult areas such as hairy sites or contours (chest), use more downward pressure to increase adhesion.
- Place all electrodes on patient before any clips are connected.
- Recommended placement order is limb lead sites, followed by chest lead sites.

3. Leadwire/Clip Connection:

- Attach clip-on connector to tab portion of electrode without clipping into adhesive.
- Clip electrodes in the order of electrode application, limb leads first and then chest leads.
- Ensure that metal of clip is in contact with conductor side (underside) of electrode tab.
- Avoid leadwire tension on electrode tab. Position cable yoke to provide adequate leadwire slack. Rotate clip/adaptor on leadwire terminal pin to reduce tension.

Note: This electrode can be repositioned upon initial application. Select a new skin site for reapplication.

Shelf Life

For shelf life, refer to the expiration date that is printed on each package.

MRI Safety Information

MR Conditional

Non-clinical testing demonstrated that the 3M™ Red Dot™ Resting EKG Electrodes 2330, 2360 and 9650 are MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)(extrapolated) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined for the 3M Red Dot resting electrodes are expected to produce a maximum temperature rise of 1.3°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

Artifact Information

In non-clinical testing, the image artifact caused by the 3M Red Dot resting electrodes extends approximately 1-mm from the 3M Red Dot resting electrodes when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

IMPORTANT NOTE: To ensure patient safety and proper recording of the electrocardiogram during an MRI procedure, use only electrocardiograph (EKG/ECG) electrodes tested specifically with regard to the MRI environment and recommended by the manufacturer of the MR Conditional monitoring equipment.

Thank you for choosing 3M™ Red Dot™ Resting Electrodes. If you are not completely satisfied or if you have questions or comments, please call the 3M Health Care Customer Helpline at 1 800 228-3957. Outside of the U.S.A., please contact your 3M Distributor or your 3M Sales Representative.

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3M and Red Dot are trademarks of 3M.
© 2016, 3M. All rights reserved.