3M™ Attest™ 13911/13913/13915 Instrument Protector en

Product Description: 3M™ Attest™ Instrument Protectors are thick, disposable paper holders designed to protect delicate surgical instruments during the sterilization process. Instrument tips are secured in a 2 mil clear film pocket to help prevent penetration of a peel pouch. It also allows for viewing and easy identification of each instrument. The anti-lock flap keeps jawed instruments open to permit steam or gas penetration.

Instrument protectors are available in three sizes:

13911: 5.16 cm x 12.9 cm (2 in. x 5 in.) for reusable needles and small surgical uninged instruments up to 12.9 cm (5 in.).

13913: 9.03 cm x 16.77 cm (3.5 in. x 6.5 in.) for instruments up to 18,06 cm (7 in.).

13915: 14.19 cm x 24.51 cm (5.5 in. x 9.5 in.) for instruments up to 28.38 cm (11 in.).

Indications for Use:

Use Attest™ Instrument Protectors to secure delicate surgical or medical instruments during any steam or EO sterilization process.

Contraindications: None.






Precautions:

Do not use in dry heat sterilization processes. Product contains cellulose. Do not use in vaporized hydrogen peroxide sterilization processes.

Directions for Use:

- Select proper size protector for instrument to be sterilized (see Product Description section above for details.)
- Insert instrument so that tip is completely covered by the clear plastic flap.
- Open hinged instruments and fold the anti-lock flap into place (for 13913 and 13915).
- Position the instrument protector into a peel pouch with the instrument facing the clear side of the peel pouch.
- Insert an appropriate internal chemical indicator per your facility policy and procedures.
- Seal peel pouch and sterilize.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

| Symbol Title | Symbol | Description and Reference |
|--------------------------|---|--|
| Manufacturer |  | Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1 |
| Date of Manufacture |  | Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3 |
| Use-by date |  | Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4 |
| Batch code |  | Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5 |
| Caution |  | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4 |
| Does not contain lead |  | Indicates product does not contain lead |
| Steam indicator |  | Indicates product is designed for use with steam sterilization processes. Source: ISO 11140-1, 5.6 |
| Ethylene oxide indicator |  | Indicates product is designed for use with ethylene oxide sterilization processes. Source: ISO 11140-1, 5.6 |

For more information see, HCBRegulatory.3M.com

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