

## **3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad**

### **Description:**

The 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad is a waterproof, bacterial and viral\* barrier dressing. The dressing consists of a non-adherent, absorbent pad bonded to a larger thin film backing with a non-latex, hypoallergenic adhesive.

\**In vitro* testing shows that the transparent film dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.

### **Indications:**

The 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad is designed for covering acute wounds such as cuts, superficial and partial- thickness burns, abrasions, IV catheter sites, and surgical incisions. It is also designed for covering low exuding superficial and partial- thickness chronic wounds including pressure ulcers, venous ulcers and diabetic foot ulcers. Follow your “gauze and tape” protocol for use. This product is not designed, sold or intended for use except as indicated.

### **Contraindications:**

None known.

### **Warnings:**

Do not use the dressing as a replacement for sutures and other primary wound closure methods.

### **Precautions:**

1. Stop any bleeding at the site before applying the dressing.
2. Do not stretch the dressing during application as tension can cause skin trauma.
3. Make sure the skin is dry and free of soap residue and lotion to prevent skin irritation and to ensure good adhesion.
4. The dressing may be used on an infected site only when under the care of a health care professional.

### **Instructions for use:**

#### **Dressing selection:**

Choose a dressing size large enough to provide a margin that adheres to dry healthy skin around the catheter or wound site.

#### **Site preparation:**

1. Prepare the site according to institution protocol.
2. Clipping of hair at the site may improve adhesion. Shaving is not recommended.
3. Allow all prep liquids to completely dry before applying the dressing to prevent skin irritation and ensure good adhesion.

#### **Application:**

1. Open package and remove sterile dressing.
2. Peel the paper liner from the paper-framed dressing, exposing the adhesive surface.
3. Position the framed window over the wound site or catheter insertion site and apply dressing.
4. Press the dressing into place.
5. Remove the paper frame from the dressing while smoothing down the dressing edges. Seal securely around catheter or wound site. Firmly smooth adhesive border to the skin.

#### **Site Care:**

1. The site should be observed for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or an unusual odor or discharge.
2. Change the dressing according to institution protocol, or when the barrier properties have been compromised.

#### **Removal:**

Gently grasp an edge and slowly peel the dressing from the skin in the direction of hair growth. Avoid skin trauma by peeling the dressing back, rather than pulling it up from the skin.

#### **Storage and Shelf Life:**

For best results, store in a cool, dry place. For shelf life, refer to the expiration date which is printed on each package.

#### **How Supplied and Ordering Information:**

Supplied in boxes of individually packaged sterile dressings. Sterility of dressing guaranteed unless individual package is damaged or opened.

For further information, please contact your local 3M representative or contact us at 3M.com and select your country.

### **Explanation of Symbols**

 Do not use if package is damaged or open

 Not Made With Natural Rubber Latex

 Sterilized using irradiation

 Do not reuse

 Batch code

 Use by date

 Caution, see instructions for use

 Manufacturer

Made in U.S.A. for

 **3M Health Care**

2510 Conway Ave.

St. Paul, MN 55144 U.S.A.

1-800-228-3957

3M.com/Tegaderm

3M and Tegaderm are trademarks of 3M.

Used under license in Canada.

© 2015, 3M. All rights reserved.

Issue Date: 2015-08

34-8716-6451-1

3M Canada Company

P.O. Box/C.P. Box 5757

London, Ontario, N6A 4T1

1-800-364-3577

3M.com/Tegaderm