












3M Advanced Wound Care Product Guide

Stage I Pressure Ulcer	Stage II Pressure Ulcer	Stage III/IV Pressure Ulcer	Unstageable Pressure Ulcer	Suspected Deep Tissue Injury
<p>Intact skin with non-blanchable redness of a localized area usually over a bony prominence. Darkly pigmented skin may not have visible blanching; its color may differ from the surrounding area.*</p>	<p>Partial thickness loss of dermis presenting as a shallow open ulcer with a red pink wound bed, without slough. May also present as an intact or open/ruptured serum-filled blister.*</p>	<p>Stage III: Full thickness tissue loss. Subcutaneous fat may be visible but bone, tendon or muscle are <i>not</i> exposed. Slough may be present but does not obscure the depth of the tissue loss. May include undermining and tunneling.* Stage IV: Full thickness tissue loss with exposed bone, tendon, or muscle. Slough or eschar may be present on some parts of the wound bed. Often includes undermining and tunneling.*</p>	<p>Full thickness tissue loss in which the base of the ulcer is covered by slough (yellow, tan, gray, green or brown) and/or eschar (tan, brown or black) in the wound bed.*</p>	<p>Purple or maroon localized area of discolored intact skin or blood-filled blister due to damage of underlying soft tissue from pressure and/or shear. The area may be preceded by tissue that is painful, firm, mushy, boggy, warmer or cooler as compared to adjacent tissue.*</p>
		 		
<p>CLEAN: 3M™ Cavilon™ No-Rinse Skin Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>COVER: 3M™ Tegaderm™ Transparent Film Dressing</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>COVER: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ Transparent Film Dressing</p> <p><u>Minimal to Moderate Drainage</u> 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing 3M™ Tegaderm™ Hydrocolloid Dressing 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing 3M™ Tegaderm™ Superabsorber Dressing</p> <p><u>Sacral Ulcer</u> 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing, Sacral Design 3M™ Tegaderm™ Hydrocolloid Dressing, Sacral Design</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ Hydrogel Wound Filler</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ Alginate Dressing</p> <p>COVER: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing</p> <p><u>Minimal to Moderate Drainage with Minimal Depth</u> 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing 3M™ Tegaderm™ Hydrocolloid Dressing 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing 3M™ Tegaderm™ Superabsorber Dressing</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>COVER: See cover dressing choices noted under Stage III/IV for management of exudate level</p> <p>NOTE: Debridement of ulcers with devitalized tissue is necessary. Notify MD/Wound Consultant for debridement options</p> <p>Stable (dry, adherent, intact without erythema or fluctuance) eschar on the heels serves as “the body’s natural (biological) cover” and should not be removed</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>COVER: 3M™ Tegaderm™ Transparent Film Dressing</p>

*National Pressure Ulcer Advisory Panel. Updated Pressure Ulcer Stages, <http://www.npuap.org/pr2.htm>. Accessed July, 2007. NPUAP images: NPAUP copyright and used with permission.



Skin Tears	Lower Extremity Venous Ulcer	Lower Extremity Arterial Ulcer	Neuropathic Ulcer/ Diabetic Ulcers	Stalled/Hard to Heal/ Bioburden Suspected
<p>A skin tear is a traumatic wound that results in separation of the epidermis from the dermis. It can be partial thickness or full thickness, usually due to friction and/or shear forces.</p>	<p>Lower extremity venous leg ulcers are caused by chronic venous hypertension. Failure of valves in the veins or faulty calf muscle pump action in the legs results in inadequate fluid return from the legs. Fluid extravasation, tissue ischemia and, eventually, ulceration is the common course of this condition.</p>	<p>Lower extremity arterial ulcers are caused by lower extremity vascular arterial disease.</p>	<p>Neuropathic ulcers result from neurologic and musculoskeletal changes leading to a lack of protective sensation and altered weight-bearing. Tissue damage most often results from trauma and/or repetitive pressure. Inadequate arterial perfusion may also be a factor.</p>	<p>A chronic wound does not heal in a timely, orderly process. Stalled in a phase of healing. Failed to heal in 4-6 weeks with standard wound care protocol.</p>
				
<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: 3M™ Tegaderm™ Alginate Dressing</p> <p>COVER: <u>Light to Moderate Drainage</u> 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ Silicone Foam Border 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing</p> <p>3M™ Coban™ LF (Latex Free) Self-Adherent Wrap [secural for Tegaderm™ High Performance Foam Non-Adhesive Dressing]</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: 3M™ Tegaderm™ Non-Adherent Contact Layer 3M™ Tegaderm™ Alginate Dressing</p> <p>COVER: 3M™ Tegaderm High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing 3M™ Tegaderm™ Superabsorber Dressing</p> <p>COMPRESS 3M™ Coban™ 2 Layer Compression System (ABPI ≥ 0.8) 3M™ Coban™ 2 Layer Lite Compression System (ABPI ≥ 0.5)</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ Non-Adherent Contact Layer 3M™ Tegaderm™ Hydrogel Wound Filler</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ Alginate Dressing</p> <p>COVER: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ Non-Adherent Contact Layer 3M™ Tegaderm™ Hydrogel Wound Filler</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ Alginate Dressing</p> <p>COVER: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing 3M™ Tegaderm™ Superabsorber Dressing</p>	<p>CLEAN: 3M™ Wound Cleanser</p> <p>PROTECT: 3M™ Cavilon™ No Sting Barrier Film</p> <p>FILL: <u>Stalled/Bioburden Suspected</u> 3M™ Tegaderm™ Alginate Ag Silver Dressing 3M™ Tegaderm™ Ag Mesh Dressing with Silver</p> <p>COVER: <u>Non-Draining to Minimal Drainage</u> 3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad 3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing</p> <p><u>Minimal to Heavy Drainage</u> 3M™ Tegaderm™ High Performance Foam Non-Adhesive Dressing 3M™ Tegaderm™ High Performance Foam Adhesive Dressing 3M™ Tegaderm™ Silicone Foam Border Dressing</p> <p><u>Moderate to Heavy Drainage</u> 3M™ Tegaderm™ Superabsorber Dressing</p>