1	3M™ Littmann® StethAssist™ Heart and Lung Sound Visualization Software User Manual - Version 4
13	3M™ Littmann® StethAssist™ Logiciel de visualisation des bruits du cœur et des poumons Manuel de l'utilisateur - Version 4
27	3M™ Littmann® StethAssist™ Software zur Visualisierung von Herz- und Lungentönen Benutzerhandbuch - Version 4
41	3M™ Littmann® StethAssist™ Software per la visualizzazione dei rumori cardiaci e polmonari Manuale d'uso - Versione 4
55	3M™ Littmann® StethAssist™ Software para la visualización de los sonidos del corazón y los pulmones Manual del usuario – Versión 4
69	3M™ Littmann® StethAssist™ Software voor de visualisatie van hart- en longgeluiden Gebruikershandleiding - Versie 4
83	StethAssist™ Littmann® da 3M™ Software de Visualização do Som do Coração e do Pulmão Manual do Utilizador - Versão 4
97	3M™ リットマン® StethAssist™ 心音・肺音可視化ソフトウェア ユーザーマニュアル - 第4版
111	3M™ Littmann® StethAssist™ 心肺音視覺軟體 軟體使用者操作手冊 - 第 4 版
125	3M™ Littmann® StethAssist™ 心肺音可视化软件 用戶手冊 - 第 4 版
139	3M™ Littmann® StethAssist™ 심음 및 폐음 시각화 소프트웨어 사용 설명서 - 버전 4

User Manual - Version 4

 SGS 9300 Winnetka Ave. N Brooklyn Park, MN 55445 (763) 488-5700		CUST ID#		PRINT METHOD		SGS# 4664696		ROUND 6				
		DATE		CHANGE LOG				INT				
Production Art	Description: 3M™ Littmann® StethAssist™ Heart and Lung Sound Visualization Software User Manual BOOK - Version 4 Base File Name: PA: 34872176853.indb (BOOK) GA: Software: InDesign CC 2014 SGS Contact: Email: Printer: Printer Location: Supplier:	Item Spec#: 34-8721-7685-3 Supersedes#: Structure#: Cat/Product#: Reference: GS Series Requester: Katelyn Early Die # / Doc. Size: 8.5" x 11"	10.16.17	New Build				kmh				
			10.24.17	Remove empty page and update all CopyRight Dates.				kmh				
INKS	 PROCESS BLACK	 PROCESS CYAN	 PROCESS MAGENTA	 PROCESS YELLOW	MATCH					SPECIAL		
LIMITATION OF LIABILITY – The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proofread all supplied materials before producing image carriers and printing final product. Should the supplied materials be found to be incorrect please contact an SGS representative immediately. The liability of SGS shall be limited to correcting its own product. SGS will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost. UPC NOTE – Because of the differences in scanning equipment, SGS cannot guarantee the scannability of UPC barcodes. The Uniform Code Council recommends a UPC symbol be printed black on white at 100% with no truncation. All UPC codes are test scanned and verified before leaving our facility.												

TABLE OF CONTENTS

1. General Information.....	1
2. Getting Started.....	2
3. Create a New Recording or Import Recording(s) from the Stethoscope.....	5
4. Review Recordings.....	7
5. Save, Open, and Export Recordings.....	9
6. Navigating Menus	9
7. Improving Recording Quality.....	11
8. Troubleshooting	12

1. General Information

1.1 Intended Use

3M™ Littmann® StethAssist™ Heart and Lung Sound Visualization Software is intended to be used with the 3M™ Littmann® Electronic Stethoscope Model 3200 for the recording and display of heart sounds and other physiological acoustic signals. The device is intended for use in a physician office setting. It is not intended as a means of diagnosis. This product does not report any clinical findings or analysis.

1.2 Safety Information

Consult Instructions for use

Please read, understand, and follow all safety information contained in these instructions prior to using this software. Retain these instructions for future reference.

A CAUTION or NOTE used in this user manual will be boxed and will have the following significance:

CAUTION: Indicates a situation in which the software may present inaccurate or incomplete information to the user.

NOTE: Information that is considered important to emphasize.

1.3 System Requirements and Specifications

Computer Requirements:

- Windows 7 32-bit, Windows 7 64-bit, Windows 8 32-bit, Windows 8 64-bit, Windows 10 32-bit, Windows 10 64-bit.
- Mac OS 10.12
- 1024 MB of RAM or greater
- Intel Core 2 Duo compatible or later
- 1024 × 768 graphic display resolution or better
- USB 1.1 port or better to support the USB Bluetooth wireless adaptor

1.4 Patient Privacy

The privacy and security of certain individually identifiable health information may be protected by state, federal, international, and/or foreign laws that regulate and restrict how such information can be used, disclosed, stored and transmitted. In some instances, you may even be prohibited from permitting any third party to access such information without first obtaining prior written authorization to do so.

You are solely and fully responsible for determining, understanding, and complying with any and all laws that affect how you may store and transmit any electronic individually identifiable health information through use of the software. If you determine that you are not able to comply with a law or restriction that applies to your use and disclosure (including transmission) of individually identifiable health information through the software, then you should not proceed to save such information. By proceeding to save any information, you hereby assume all risks, responsibilities and liabilities associated with using, storing, and/or transmitting such information, including sound files, electronically or otherwise.

This program does not require you to enter individually identifiable health information in order to function. Records can be stored and recalled through the use of a timestamp or an arbitrarily assigned identifier. If you choose to enter a meaningful Patient Name or Patient Identifier, then you agree to assume any and all risks of and liabilities incurred with transmitting and storing such individually identifiable health information through the software, including risks to the privacy and security of such individually identifiable health information.

2. Getting Started

2.1 Install and Register the 3M™ Littmann® StethAssist™ Software

Follow these steps to install 3M™ Littmann® StethAssist™ Software on your computer:

1. Go to www.Littmann.com/StethAssist
2. Click on the Download StethAssist link
3. Click Run on the message box that appears
4. Follow the installation instructions
5. Complete online registration

NOTE: To complete registration offline, please contact 3M™ Littmann® Customer Support at 1-800-228-3957 or littmannsoftwaresupport@mmm.com.

2.2 Insert the USB Bluetooth Wireless Adaptor (Dongle)

To install the USB Bluetooth wireless adaptor, insert the device into a USB port on the computer. The USB Bluetooth wireless adaptor allows communication between the computer on which the software is installed and the 3M™ Littmann® Electronic Stethoscope Model 3200. The USB Bluetooth wireless adaptor receives and transmits data via the commonly accepted Bluetooth® short range wireless connectivity standard.

To connect to a stethoscope, **the USB Bluetooth wireless adaptor must be inserted into the computer prior to launching the 3M™ Littmann® StethAssist™ Software**. If the USB Bluetooth wireless adaptor is not detected by the software, the following message will be displayed. If this message is displayed, insert the USB Bluetooth wireless adaptor, and then close and re-launch the StethAssist software.

2.3 Pair & Connect the Stethoscope with the Software

The following screen will help you pair and/or connect a stethoscope to the software.


2.3.1 Pair New Stethoscope

To connect a Littmann® Model 3200 Stethoscope to the 3M™ Littmann® StethAssist™ Software for the first time, you must first pair the stethoscope to the software. To pair a new stethoscope, select “Add New Stethoscope” from the dropdown menu and click Pair (or, select Add New Stethoscope from the Tools menu on the main screen). Follow the onscreen instructions. After the stethoscope is paired, the software will automatically connect to it. Every time you launch the stethoscope thereafter, simply connect to the stethoscope.

2.3.2 Connect Stethoscope

To connect to a previously paired stethoscope, select the stethoscope name from the dropdown menu and click Connect. On the stethoscope, press the (M) button twice to select “CONNECT”. The screen below will appear each time you attempt to connect a stethoscope to the software. To review recordings without connecting a stethoscope, click Cancel.

If your stethoscope is connected, the Stethoscope Connection button will display a green check mark . If not connected, the Stethoscope Connection button will display a red X .

NOTE: Follow these tips to optimize wireless connections

1. Auscultate the patient as close to the USB wireless adaptor as possible
2. Avoid having objects or people between the stethoscope and USB Bluetooth wireless adaptor
3. Avoid having your hand cover the entire chestpiece of the stethoscope (location of the stethoscope's Bluetooth® antenna)

NOTE: You can pair multiple stethoscopes to the 3M™ Littmann® StethAssist™ Software, but only one stethoscope can be actively connected at a time. The stethoscope can only be actively connected to one computer.

3. Create a New Recording or Import Recording(s) from the Stethoscope

Recording Set-Up Window

The screenshot shows the 'Recording Set-Up' window with the following fields and controls:

- Patient ID:** Text box containing '01234'. Annotation: 'Enter patient information (optional)'.
- Patient Name:** Three text boxes for 'First' (Jane), 'Middle' (empty), and 'Last' (Doe). Annotation: 'Enter patient information (optional)'.
- Recording Length (Seconds):** A slider and a text box showing '15'. Annotation: 'Select the recording length (5-60 seconds)'.
- Select Recording Site:** A dropdown menu showing 'Mitral site (Apex)'. Annotation: 'Select a recording site from the dropdown or click on the blue torso to manually designate the site. To switch to a different image, click an alternate thumbnail view.'
- Patient Posture:** A dropdown menu showing 'Supine'. Annotation: 'Select the patient posture'.
- Torso Image:** A blue-tinted image of a human torso with a white dot on the left chest. Annotation: 'Select a recording site from the dropdown or click on the blue torso to manually designate the site. To switch to a different image, click an alternate thumbnail view.'
- Thumbnails:** A row of six small thumbnail images below the torso image.
- Recording 1 of 12:** Text at the bottom of the window.
- Buttons:** 'OK' and 'Cancel' buttons at the bottom.

3.1 Create a New Recording

- 3.1.1. To start a new recording, click on the New Encounter button. In the Recording Set-up Window, enter the Patient ID and Patient Name (optional, please see Section 1.4 for patient privacy information). Patient IDs may contain any combination of letters and numbers up to 20 characters long. The Patient ID will only allow the following characters: a-z, 0-9, -, and _
- 3.1.2. Select the recording length, recording site and patient posture. Then click OK.
- 3.1.3. Obtain good sound quality at the desired patient site, and then press the (M) button on the stethoscope to begin recording.
- 3.1.4. To restart the recording, press the (M) button on the stethoscope to stop the recording, and then press the (M) button again to restart the recording.
- 3.1.5. After the recording is completed, you may make another recording with the same patient by clicking the “Add Recording” button. To make a recording on a different patient, select the “New Encounter” button. The 3M™ Littmann® StethAssist™ Software will automatically ask if you would like to save the current set of recordings.
- 3.1.6. Click on the Summary view to enter patient comments (optional)

NOTE: The volume level of the stethoscope does not affect the volume level of recordings. Also, the recording can be played back in all three filter settings (Bell, Diaphragm, Extended Range), regardless of what filter setting was used during the recording.

3.2 Import Recording(s) from the Stethoscope Memory

To transfer a saved recording(s) from the on-board memory of Littmann® Model 3200 Stethoscope to the software, click on the Import button or select Import From Stethoscope from the File menu. Follow the steps below:

1. Select Recorded Track(s) to Import

2. Enter Patient Information (optional).

3. Click on the Play icon to listen to the track Comment (if present). Note: Track Comment is played through stethoscope eartips.

4. Enter Patient Recording Site (optional).

5. Enter Patient Posture (optional).

Listen to Recorded Track (optional)


Lock or Unlock the Recorded Track in the stethoscope memory (optional)

Delete Recorded Track from the stethoscope memory (optional)

6. After entering the information for all of the recorded tracks, click Import Selected.

The screenshot shows a software window titled "Import Stethoscope Recordings". At the top, it says "Listen through the stethoscope to a track comment, then enter patient info, recording site and posture." Below this is a table with 12 rows, each representing a track. The columns are: Track (with checkboxes), Patient Info (with a play button and a dropdown), Recording Site (with a dropdown), Posture (with a dropdown and a play button), Lock/Unlock (with a lock icon), Delete (with a red X icon), and Imported (with a checkmark icon). Step 1 points to the checkboxes in the Track column. Step 2 points to the Patient Info dropdown. Step 3 points to the play button in the Patient Info column. Step 4 points to the Recording Site dropdown. Step 5 points to the Posture dropdown. A callout box "Listen to Recorded Track (optional)" points to the play button in the Posture column. Another callout box "Lock or Unlock the Recorded Track in the stethoscope memory (optional)" points to the Lock/Unlock column. A third callout box "Delete Recorded Track from the stethoscope memory (optional)" points to the Delete column. At the bottom, there are "Import Selected" and "Close" buttons. Step 6 points to the "Import Selected" button.

4. Review Recordings

4.1 Signal Pan & Zoom

To zoom in or out within a sound display, first select the recording by clicking on it, then use either the Horizontal Zoom or Vertical Zoom “+” or “-” buttons. Zooming in provides more detail. Zooming out provides less detail. The Vertical Zoom will also increase or decrease the volume of the recording as it affects both the visual and the audible amplitude of the recording.

If zoomed in horizontally, you can scroll through the recording using the Pan buttons. To move the recording to the left, click the Pan Left button. To move the recording to the right, click the Pan Right button. It is also possible to pan left and right through the sound displays by using the Drag button located on the toolbar. Click this hand-shaped button and you have the ability to pan left or right through any of the recordings by placing the mouse cursor on a sound display and then dragging the cursor in the direction you wish to scroll. Deactivate this feature by clicking the Drag button again.

4.2 Spectrogram Display Option

The spectrogram is a time-frequency distribution presented as a two-dimensional color plot of the sound energy distributed over time and frequency. Black represents minimum energy while red is maximum energy. The x-axis represents the time index while the y-axis is a linear frequency range from 0 – 500 Hz.

The spectrogram provides an additional perspective on the recorded sound. To view a sound display in spectrogram mode, click the Spectrogram **Sp** button appearing to the left of the sound display. To return to the standard, default sound display view, click the Phonocardiogram **Ph** button that appears to the left of the sound display.

4.3 *Playing a Recording*

Press the Play button to begin playback. A red cursor will scroll across the sound display in synchrony with the recording. Click the Stop button to end the playback.

By default, the software will play the sound track through the Littmann® Model 3200 Stethoscope (if connected). To play the sound track through the computer speakers, external speakers, or headphones, click the Play to Headphones checkbox.

CAUTION: The audio quality of the sounds played through the computer speakers or computer audio jack may not be a faithful reproduction of the actual sounds due to the acoustical performance of the computer’s sound card, sound card drivers, and/or the frequency response of the headphones/speakers, potentially resulting in not being able to hear clinically relevant sounds. It is strongly recommended that a 3M™ Littmann® Stethoscope Model 3200 is used to listen to the recorded sounds.

NOTE: The amplitude or volume of the recorded sounds may be less than what was heard during the initial recording. Use the “+” button on the vertical axis of the waveform graph or on the stethoscope keypad to increase the volume.

4.4 *½ Speed*

The ½ Speed option allows you to play back recordings at half speed with no significant change in the spectral properties of the signal. Click the ½ Speed checkbox, and then click the Play button to play back the recording at half speed. Once clicked, the ½ Speed checkbox stays activated until it is clicked again.

CAUTION: The ½ Speed option is available to assist in identifying physiological events, but should not be solely used to make a physiological diagnosis.

4.5 *Playing Back a Specific Portion of a Recording*

To play back a specific portion of a recording, simply click with your mouse on the area of the sound display where you want playback to begin and drag your mouse to the right to highlight the area of interest. Once an area of the sound display has been highlighted in this way, click the Play button to listen to the specific section. Note that you must select a segment of sound of at least one second in duration in order for this feature to function. Clear the highlighted area by left-clicking anywhere on the sound display.

4.6 *Looped Playback*

The Loop checkbox appearing at the top of the screen allows you to play back an entire recording, or a specific portion of a recording, in a loop that will continue until the Stop button is clicked. To play a full recording in loop mode, select the corresponding sound display by clicking on it, click the Loop checkbox, and then click the Play button. To play a specific portion of a recording in loop mode, use the cursor to highlight the area of interest on a sound display, click the Loop checkbox, and then click the Play button. Once clicked, the Loop checkbox stays activated until it is clicked again.

4.7 *Global Scaling*

The global scaling feature normalizes all of the recordings contained within a single recording set. This allows the user to accurately compare the amplitudes of multiple recordings.

4.8 *Printing Feature*

To print a copy of the “Patient Encounter Report” for a given recording or recording set, click the Print button. The resulting print-out will display the patient ID, patient name (if entered), and the ID of the stethoscope used for the recordings. The report will also contain a listing of the posture, auscultation site and time-length for the recording(s) in question, along with any corresponding user comments that were entered.

4.9 Adding Comments to a Recording

While in the “Detailed View”, you can add comments to the sound file by switching to the Summary view or by clicking on the “Comments” button located to the left of the Patient ID. A green check mark on the Comments button indicates that comments exist for the recording(s).

5. Save, Open, and Export Recordings

5.1 Save Recording

Click the “Save” button in the upper left-hand corner of the screen or select Save from the File menu. The StethAssist software will automatically prompt you to save the recording if starting a new recording or exiting the software program.

5.2 Open a Previously Recorded Sound File

Click the “Open” button in the upper left-hand corner of the screen or select Open from the File menu. Click on the heading to sort by that field.

5.3 Export Recordings to .wav file

To export a recording in a .wav file format, select Export recording as .wav from the File menu. The .wav file format is very common and can be played using most computer media players. Be sure to select the filter that is applied to the recording i.e. Bell, Diaphragm or Extended.

6. Navigating Menus

3M™ Littmann® StethAssist™ Software contains the following menu selections:
File Tools View Help

6.1 File

File > Save
Save the current recording set.

File > Open
Open a previous recording set.

File > Close
Close the current recording set.

File > New Encounter
Begin a new patient recording.

File > Comments
Open the comments window.

File > Import From Stethoscope
Import a recording from the stethoscope’s on-board memory.

File > Import Saved Recording(s) to Current Directory
Transfer previously saved recording files from another location on the computer to the current recording directory. Moving files out of the current directory will improve program startup speed as well as reduce the amount of time it takes to open a previously saved file.

File > Export Saved Recording(s)
Transfer previously saved recording files from the current recording directory to another location on the computer. Moving files out of the current directory will improve program startup speed as well as reduce the amount of time it takes to open a previously saved file.

File > Import e4k File

Import a file in the .e4k format. The .e4k file format was used in the previous 3M™ Littmann® Electronic Stethoscope Models 4000 and 4100.

File > Export recording as .wav

Open the File Export window to save the recording in the .wav file format.

File > Print Report

Print the Patient Encounter Report for a given recording set.

File > Recently Viewed

Open a recently viewed recording set.

File > Exit

Exit the StethAssist application.

6.2 Tools

Tools > Add New Stethoscope

Open the “Add New Stethoscope” window.

Tools > Stethoscope Settings

Open the Stethoscope Settings window, where you can:

- Delete a paired stethoscope.
- Change the stethoscope name.
- Adjust the length of the “lead time”, or waiting period, before the StethAssist™ Software begins the recording after pressing the (M) button on the stethoscope. Lead time options: 0, 1, 2 or 4 seconds.
- Select the default patient posture.
- Select the default recording length.
- Select the default view (summary or detailed).

Tools > Options

Change the default directory where the recordings are saved and/or change the default language.

6.3 View

Select Detailed or Summary View

6.4 Help

Help > Online Help

Open an internet browser and connect to the StethAssist webpage at www.littmann.com/StethAssist

Help > User Manual

Open the software user manual.

Help > Learn About Littmann® Products

Open an internet browser and connect to Littmann.com

Help > About

View the product version number along with phone numbers and addresses for requesting further information or making service inquiries.

Help > Check for Latest Version

Check for newer versions of StethAssist.

7. Improving Recording Quality

The quality of the recordings depends on the environment, the patient, and user technique:

The acoustic environment:

Conduct the auscultation procedure in an environment that is as quiet as possible. Listen carefully for background sounds, such as that caused by traffic noises entering through open windows or by fans and air conditioners, and eliminate these noises to the greatest extent possible. Make sure that there are no background conversations and prevent interruptions during recording.

The patient:

If possible, the patient should be properly prepared so that the stethoscope can make direct contact with the patient's skin. The patient should be instructed to remain as still as possible throughout the recording.

The placement and handling of the stethoscope:

When using the Littmann® Model 3200 stethoscope, there are several techniques that are important to consider that can substantially improve your auscultation experience.

The sound sensor, where all sound acquisition occurs, is quite small and located in the very center of the chestpiece (for demonstration purposes, you can peel off the diaphragm and look for the indentation under the white sticker). The area around the sound sensor does not contribute to sound acquisition. Therefore, it is important to focus the point of contact with the patient on the very center of the chestpiece.

Only light contact is necessary with the Littmann® Model 3200 stethoscope. This includes both the handling of the stethoscope and the contact on the patient. Excessive pressure on the patient may result in artifact sounds.

Contact with the patient must remain steady and consistent for optimal sound quality. Movement with the chestpiece will induce frictional noise.

It is a best practice to start your auscultation exam at a lower amplification level. If you are trying to listen to a difficult-to-hear sound, adjust the amplification upward. Using the amplification on an as-needed basis will help reduce distracting noises during your exam.

Finally, it is important to recognize the difference in technique between non-electronic and electronic stethoscopes. Holding the Littmann® Model 3200 stethoscope lightly, making light contact with the patient, focusing the point of contact on the center of the chestpiece, and remaining steady will ensure the highest sound quality auscultation experience.

8. Troubleshooting

Problem	Potential Cause	Solution
Cannot connect to stethoscope	No USB Bluetooth wireless adaptor installed/ USB Bluetooth wireless adaptor not properly installed	Insert USB Bluetooth wireless adaptor. Disable the internal computer Bluetooth radio. Update the Bluetooth driver to the Generic Bluetooth Adaptor
	Stethoscope not paired	Pair stethoscope. Select Add New Stethoscope from the Tools menu on the main screen
	Each unique Stethoscope ID can only represent one stethoscope	If your office has multiple stethoscopes, make sure you are trying to connect to the correct Stethoscope ID.
	Undetermined	Turn off the stethoscope, wait a few seconds, and then turn the stethoscope back on and try again.
Noisy recorded signal	Excessive background noise	Take actions to reduce the background noises heard in the auscultation environment.
Signal shown in the display plot is clipped	Noisy recording	See the Section 7 of this user manual titled, "Improving Recording Quality".
Poor data quality reported	Excessive background noise	Take actions to reduce the background noises heard in the auscultation environment.
	Excessive stethoscope motion	Hold the stethoscope with a steady hand, lightly pressing against the patient's skin while recording.

For questions or concerns regarding to the 3M™ Littmann® StethAssist™ Software or the Littmann® Model 3200 Stethoscope, please contact our customer service hotline at 1-800-228-3957 or find your local country 3M™ Littmann® Stethoscope website at www.Littmann.com for Frequently Asked Questions and the Contact Us page.

3M Health Care
Littmann® Brand Products
3M Center, Building 275-4E-01
2510 Conway Ave.
St. Paul, MN 55144
1-800-228-3957
www.littmann.com

3M, Littmann, and StethAssist are registered trademarks of 3M, used under license in Canada.
Windows is a registered trademark of Microsoft. Bluetooth is a registered trademark of the Bluetooth SIG, Inc.
Mac is a registered trademark of Apple Inc.
© 2017 3M. All rights reserved.
34-8721-7685-3

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The screenshot shows a software window titled 'Import Stethoscope Recordings'. It contains a table with 12 rows, each representing a recording. The columns are: Track (checkbox), Patient Info (dropdown), Recording Site (dropdown), Posture (dropdown), Lock/Unlock (lock icon), Delete (X icon), and Imported (checkbox). Arrows point from numbered instructions to specific elements in the window:

- 1. Points to the 'Track' column checkboxes.
- 2. Points to the 'Patient Info' dropdowns.
- 3. Points to the 'Afspeel' (play) button next to a track.
- 4. Points to the 'Recording Site' dropdowns.
- 5. Points to the 'Posture' dropdowns.
- 6. Points to the 'Import Selected' button at the bottom.

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
The screenshot shows a software window titled "Import Stethoscope Recordings". At the top, it says "Listen through the stethoscope to a track comment, then enter patient info, recording site and posture." The window contains a table with 12 rows, each representing a recording track. The columns are: "Track" (with checkboxes), "Patient Info" (with a play button and a "< Select >" dropdown), "Recording Site" (with a "< Select >" dropdown), "Posture" (with a play button and a "< Select >" dropdown), "Lock/Unlock" (with a lock icon), "Delete" (with a red X icon), and "Imported" (with a checkmark icon). At the bottom are "Import Selected" and "Close" buttons. Annotations 1 through 6 point to specific elements: 1 points to the "Track" checkbox; 2 points to the "Patient Info" dropdown; 3 points to the play button in the "Patient Info" column; 4 points to the "Recording Site" dropdown; 5 points to the "Posture" dropdown; 6 points to the "Import Selected" button. Three callout boxes provide additional instructions: one points to the play button in the "Recording Site" column with the text "録音トラックを再生します (任意)"; another points to the lock icon in the "Lock/Unlock" column with the text "聴診器のメモリーの録音トラックをロックするかまたはロックを解除します (任意)"; and a third points to the red X icon in the "Delete" column with the text "録音トラックを聴診器のメモリーから削除します (任意)".

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The screenshot shows a software window titled "Import Stethoscope Recordings" with the instruction "Listen through the stethoscope to a track comment, then enter patient info, recording site and posture." The window contains a table with 12 rows, each representing a recording track. The columns are: Track (checkbox), Patient Info (dropdown), Recording Site (dropdown), Posture (dropdown with a play button), Lock/Unlock (lock icon), Delete (red X), and Imported (checkbox). Annotations point to various elements: 1. Points to the checkbox in the first row. 2. Points to the Patient Info dropdown in the first row. 3. Points to the play button in the first row. 4. Points to the Recording Site dropdown in the first row. 5. Points to the Posture dropdown in the first row. 6. Points to the "Import Selected" button at the bottom. Additional callouts include: "收聽錄音 (可選)" pointing to the play button, "在聽診器記憶體中鎖定或解開鎖定錄音 (可選)" pointing to the lock icon, and "從聽診器記憶體刪除錄音 (可選)" pointing to the red X icon.

Track	Patient Info	Recording Site	Posture	Lock/Unlock	Delete	Imported
☒	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	
☐	< Select >	< Select >	< Select >	🔒	✖	

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