## 3M™ Steri-Strip™ C1544 Compound Benzoin Tincture

(EN)

### Applicator Exterior Sterile

#### Product Description

3M™ Steri-Strip™ Compound Benzoin Tincture is an alcoholic solution of compound benzoin.

#### Indications for Use

Steri-Strip™ Compound Benzoin Tincture is intended for use as an adhesive adjunct, when enhancement of surgical tape and/or adhesive wound closure adhesion is desired. This device is to be used by medical professionals.

#### Contraindications

Steri-Strip™ Compound Benzoin Tincture should not be used on patients with previous history of sensitivity to compound benzoin tincture.

#### Warnings

- WARNING flammable liquid and vapor.
- Compound Benzoin Tincture is flammable until it has completely dried on the skin.
- Compound Benzoin Tincture should only be applied when no ignition sources, flames, or heat-producing devices are in use.
- This product is not intended for reuse. Reuse may result in patient infection.

#### Precautions

- Allergic contact dermatitis to compound benzoin tincture has been reported.
- Spillage of compound benzoin tincture into the wound should be avoided, as it may impair the wound healing process.
- Avoid contact with eyes and/or mucous membranes.
- For external use only.

#### Directions for Use

- Clean and dry skin for at least 2 ½ inches (6.4 cm) around the wound making sure to remove skin oil or exudate. For best results, use a sponge moistened with alcohol or saline followed by a firm wipe with a dry sponge.
- Using aseptic technique, peel back lid and deliver vial.
- Gently squeeze vial to crush inner ampoule, thereby releasing compound benzoin tincture. Do not roll crushed ampoule between fingers as this may cause glass to penetrate outer vial.
- Carefully apply a thin film of compound benzoin tincture to the skin on each side of the wound with strokes parallel to the wound. Begin at the wound edges and work outward. Avoid contact with the wound.
- Allow tincture to dry before applying surgical tapes and/or adhesive wound closure bandages (approximately 30-45 seconds).

#### Storage/Shelf Life/Disposal


This product should be stored at room temperature. Avoid excessive heat or humidity. Do not use if expired or if expiration date is missing or illegible. Dispose of product per facility protocol.

#### How supplied

Supplied in individual packages. Applicator sterile unless package is opened or damaged. Please report a serious incident occurring in relation to the device to 3M.

**For further information, please contact your local 3M representative or contact us at 3M.com and select your country.**

#### Explanation of Symbols

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Unique device identifier		Indicates bar code to scan product information into patient electronic health record
Flammable		Product is classified as flammable. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals

**For more information see, HCBRegulatory.3M.com**

Made in USA with globally sourced materials

#### 3M Company

2510 Conway Ave.  
St. Paul, MN 55144 USA  
1-800-228-3957 (USA Only)  
3M.com

Issue Date: 2020-02

3M and Steri-Strip are trademarks of 3M. Used under license in Canada. © 2020, 3M. Unauthorized use prohibited. All rights reserved.

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