

- (en) Warming Unit Model 775 Operator's Manual
- (fr) Manuel d'utilisation de l'unité de réchauffement modèle 775
- (de) Benutzerhandbuch zur Wärmeeinheit, Modell 775
- (it) Manuale d'uso dell'Unità riscaldante, modello 775
- (es) Manual del operador de la unidad de calentamiento modelo 775
- (nl) Bedieningshandleiding verwarmingsunit model 775
- (sv) Värmeenhet modell 775 Bruksanvisning
- (da) Betjeningsmanual til Warming Unit Model 775
- (no) Brukerhåndbok for varmeenhet modell 775
- (fi) Lämmitysyksikkö malli 775 käyttöopas
- (pt) Manual do Utilizador da Unidade de Aquecimento, Modelo 775
- (el) Εγχειρίδιο Χειριστή Μονάδας Θέρμανσης, Μοντέλο 775
- (pl) Aparat do terapii grzewczej Model 775 — instrukcja obsługi
- (hu) Melegítőegység, 775-ös modell – Kezelői kézikönyv
- (cs) Ohřívací jednotka, model 775 – Návod k obsluze
- (sk) Návod na obsluhu ohrievacej jednotky model 775
- (sl) Navodila za uporabo grelne enote model 775
- (et) Soojendusseadme mudel 775 kasutaja käsiraamat
- (lv) Sildītāja modelis 775, lietošanas rokasgrāmata
- (lt) Šildymo įrenginio modelio 775 naudotojo vadovas
- (ro) Manualul operatorului pentru unitatea de încălzire model 775
- (ru) Инструкция по эксплуатации устройства для обогрева пациентов, модели 775
- (hr) Korisnički priručnik za uređaj za grijanje, Model 775
- (bg) Ръководство за оператора на затоплящо устройство модел 775
- (tr) Isıtma Ünitesi Model 775 Kullanım Kılavuzu
- (zh) 775 型加温器操作手册
- (ar) دليل المشغل لوحدة الاحترار طراز 775

(en) English.....	1	(cs) Čeština.....	164
(fr) Français.....	12	(sk) Slovensky	175
(de) Deutsch.....	24	(sl) Slovenski	186
(it) Italiano	36	(et) EESTI.....	197
(es) Español.....	48	(lv) Latviesu	208
(nl) Nederlands	60	(lt) Lietuviskai.....	219
(sv) Svenska	72	(ro) Romana	230
(da) Dansk	83	(ru) РУССКИЙ	241
(no) Norsk.....	94	(hr) HRVATSKI	255
(fi) Suomi	105	(bg) БЪЛГАРСКИ	266
(pt) Português	116	(tr) Turkce.....	279
(el) Ελληνικά.....	128	(zh) 简体中文.....	290
(pl) Polski	141	(ar) اللغة العربية	300
(hu) Magyar.....	153		

Introduction	2
Description of the 3M™ Bair Hugger™ Temperature Management System	2
Indications for use	2
Patient Population and Settings	2
Contraindications, Warnings, Cautions and Notices	2
Explanation of Signal Word Consequences	2
Notices	3
Proper Use and Maintenance	3
Read Before Servicing Equipment	3
Warming Unit Power-On-Reset	4
Overview of Control Panel	4
Mounting the Warming Unit to an IV Pole	5
Instructions for Use	5
Viewing the Temperature Mode Timer	6
What to Do in Case of an Over-Temperature Condition	6
What to Do If a Fault Condition occurs	6
Cleaning Instructions	6
Notice:	6
Storage	6
When You Call for Technical Support	6
In-Warranty Repair and Exchange	6
Symbol Glossary	7
Specifications	8
Physical Characteristics	8
Temperature Characteristics	8
Average temperatures at the end of the hose:	8
Classification	8
Diagnostics	8
Safety System	9
Electrical Characteristics	9

Operator's Manual

Introduction

Description of the 3M™ Bair Hugger™ Temperature Management System

The Bair Hugger temperature management system consists of a Model 775 warming unit (with optional rolling cart and other accessories) along with a 3M™ Bair Hugger™ warming blanket or 3M™ Bair Hugger™ warming gown.

The warming unit is attached to the warming blanket or warming gown by means of a flexible hose. Warm air is generated in the warming unit and flows through the warming unit hose and into the blanket or gown. Depending on the model, the 3M blanket or gown is placed either around, over, or underneath the patient. Small perforations on the blanket or gown allow the warm air to be dispersed over the patient.

This manual includes operating instructions and warming unit specifications for the Model 775 warming unit. You can use the Model 775 warming unit in clinical settings where the patient may feel too warm or too cold, including the operating room to provide patient temperature management. For information about using 3M blankets or gowns with the Model 775 warming unit, refer to the respective Instructions for Use.

Indications for use

The Bair Hugger family of temperature management systems are indicated for hypothermic patients or normothermic patients for whom induced hypothermia or localized temperature therapy is clinically indicated. In addition, the temperature management systems can be used to provide patient thermal comfort when conditions exist that may cause patients to become too warm or too cold. The temperature management systems can be used with adult and pediatric patients.

The Bair Hugger temperature management system should only be used by trained medical professionals.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency departments and other departments in the hospital setting where patient temperature management is required.

Contraindications, Warnings, Cautions and Notices

Explanation of Signal Word Consequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE: Indicates a situation which, if not avoided, could result in property damage only.

CONTRAINDICATION: To reduce the risk of thermal injury:

- Do not apply heat to lower extremities during aortic cross-clamping. Thermal injury may occur if heat is applied to ischemic limbs.

WARNING: To reduce the risk of thermal injury:

- The Bair Hugger Model 775 warming unit has been designed to operate safely ONLY with 3M disposable warming products. Use with other products may cause thermal injury. To the full extent permitted by law, the manufacturer and/or importer declines all responsibility for thermal injury resulting from the warming unit being used in conjunction with non-3M products.

- Do not treat patients with the warming unit hose alone. Always attach the hose to a 3M warming blanket/gown before providing warming therapy. Do not allow the patient to lie on the warming unit hose. Do not allow the warming hose to directly contact the patient's skin during warming therapy.
- Do not leave neonates, infants, children and other vulnerable patient populations unattended during warming therapy.
- Do not leave patients with poor perfusion unmonitored during prolonged warming therapy.
- Do not place the non-perforated side of the warming blanket/gown on the patient. Always place the perforated side (with the small holes) directly on top of the patient in contact with the patient's skin.
- Do not connect a torn or damaged warming blanket/gown to the warming unit.
- For over-body warming blankets and gowns; do not place patient securement device (i.e. safety strap or tape) over the warming blanket/gown.
- For underbody or side channel warming blankets; if a securement device (i.e. safety strap, tape) is used, ensure the warming channels are not occluded.
- Do not place the warming blanket/gown directly over a dispersive electrode pad.
- Do not continue warming therapy if the red Over-temp indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.
- For the Bair Hugger flex warming gown; ensure that the blood pressure cuff, ECG, IV or other lines or cords are not between the upper sleeve insert and the gown prior to deploying the upper sleeve warming insert(s), as this could result in tearing of the insert during deployment.
- Do not perform the over-temperature detection system test while the warming unit is being used for warming therapy.

WARNING: To reduce the risk of patient injury or death due to altered drug delivery:

- Do not use a warming blanket/gown over transdermal medication patches.

WARNING: To reduce the risk of injury due to interference with ventilation:

- Do not allow the warming blanket/gown or head drape to cover the patient's head or airway when the patient is not mechanically ventilated.

WARNING: To reduce the risk of injury due to patient falls:

- Do not use a warming blanket/gown to transfer or move the patient.

WARNING: To reduce the risks associated with hazardous voltage and fire:

- Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device.
- Only connect to outlets marked "Hospital only," "Hospital Grade," or a reliable grounded outlet.
- Use only the 3M-supplied power cord intended for this product and certified for the country of use.
- Do not allow the power cord to get wet.

- Do not use the warming unit when it appears the warming unit, power cord or any component is damaged. Replace the warming unit. Contact 3M Technical Support.
- Do not disassemble the warming unit unless you are a qualified service technician. There are electrically live parts within the warming unit when it is connected to a power source.
- Connect each warming unit being tested to a separate power source.

CAUTION: To reduce the risk of cross-contamination:

- Except for specific Bair Hugger warming blanket models, 3M warming blankets/gowns are not sterile. Each warming blanket/gown is intended for single patient use ONLY. Placing a sheet between the warming blanket/gown and the patient does not prevent contamination of the product.
- Clean the warming unit and the outside of the warming unit hose after each patient use. See “Cleaning Instructions” section.
- Follow applicable regulations when disposing of this warming unit or any of its electrical components.
- Do not attempt to clean the air filter as it may be contaminated from use. Discard the filter in a manner consistent with institutional protocol.
- Do not operate the Bair Hugger warming unit with the hose detached from the 3M warming blanket/gown. The Bair Hugger warming unit is intended to be used with the hose properly attached to its corresponding warming blanket/gown and in accordance with good practices for operating room sterile technique.

CAUTION: To reduce the risk of patient or caregiver injury:

- If mounted on an IV pole, the distance from the bottom of the warming unit to the floor must be less than 44" (112 cm) and the IV pole wheelbase diameter must be at least 28" (71 cm) to prevent tipping.

CAUTION: To reduce the risk of fire:

- 3M warming blankets and gowns are classified as Class I Normal Flammability as defined by the Consumer Product Safety Commission's flammable fabric regulation, 16 CFR 1610. Follow standard safety protocols when using high intensity heat sources.

CAUTION: To reduce the risk of thermal injury, hyperthermia or hypothermia:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.
- Do not place the warming unit on a soft uneven surface, such as a bed, or a visibly wet surface as the air intake may become blocked and cause the warming unit to overheat, compromising the warming unit's performance.
- Perform all temperature testing of the warming unit with a 3M Model 22110 temperature test unit. The model 22110 temperature test unit itself requires no calibration.

CAUTION: To reduce the risk associated with electromagnetic interference (EMI) due to portable and mobile radio frequency (RF) communications equipment:

- The 3M Bair Hugger temperature management system has been tested to be resistant to both EMI and electrostatic discharge (ESD). Install and put into service the 3M Bair Hugger temperature management system according to the electromagnetic compatibility (EMC) information provided in the Guidance and Manufacturer's Declaration.
- Should interference occur, move away from the portable or mobile RF communications equipment.

Notices

1. The Bair Hugger warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, connect the warming unit to a different power source.
2. To avoid warming unit damage:
 - Use proper Electrostatic Discharge (ESD) procedures when performing maintenance.
 - Do not modify this equipment without authorization from the manufacturer.
 - Do not immerse the warming unit, warming unit parts or accessories in any liquid or subject them to any sterilization process.
 - Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.

Proper Use and Maintenance

3M assumes no responsibility for the reliability, performance, or safety of the warming unit if the following events occur:

- Modifications or repairs are performed by unqualified personnel.
- The warming unit is used in a manner other than that described in the Operator's Preventive Maintenance Manual.
- The warming unit is installed in an environment that does not meet the appropriate electrical and grounding requirements.

The warming unit not maintained in accordance with the procedures described in the Preventive Maintenance Manual.

Read Before Servicing Equipment

All repair, calibration and servicing of the warming unit require the skill of a qualified, medical equipment service technician who is familiar with good practice for medical device repair. If service does not require the manufacturer's attention, the Model 775 Preventive Maintenance Manual provides the technical information needed to service the warming unit. Perform all repairs and maintenance in accordance with the instructions in the Preventive Maintenance Manual. For additional service information please contact 3M.

Overview and Operation


Figure A. Model 775 Warming Unit – front view


Figure B. Model 775 Warming Unit – rear view

Warming Unit Power-On-Reset

The Model 775 warming unit automatically performs the following power-on-reset sequence after you connect the warming unit to a properly grounded power source:

- Performs all self-test functions.
- Illuminates indicator lights and all pixels in the alphanumeric display momentarily.
- Displays the text “MD 700’s” and the software revision level in the alphanumeric display.
- Sounds the alarm (three low-level clicks).
- Enters the **Standby** mode with the high fan setting selected by default.

Note: The low fan setting indicator light will not illuminate during power up.

If the warming unit loses power for less than 1 second, the warming unit software will return the warming unit to the operating mode you selected prior to the power loss. However, if the warming unit loses power for longer than 1 second, the warming unit software will reset when you restore power. The warming unit will then enter the **Standby** mode with the *high fan setting selected by default*.

Overview of Control Panel


Figure C. Model 775 Warming Unit Control Panel

ALPHANUMERIC DISPLAY

The alphanumeric display shows the temperature at the end of the Bair Hugger warming unit hose in degrees Celsius.

TEMPERATURE MODES

- Press the 32°C (Low), 38°C (Med), or 43°C (High) button to select the desired temperature.
- Press the Ambient button to supply room temperature air.

When you select a temperature mode, the following events occur:

- Corresponding temperature indicator and fan setting lights illuminate.
- Blower operates at selected fan speed.
- Temperature at the end of the hose appears in the alphanumeric display.
- Heater activates except in ambient mode.
- Temperature mode timer and hour meter activate.
- Temp in range indicator light illuminates when the temperature at the blanket end of the hose is within +/- 1.5°C of the selected setting; this indicator light does not illuminate in the ambient mode.

AIRFLOW MODES

The Model 775 warming unit has two preset fan speed settings: a default or high fan setting (🌀) and a reduced or low fan setting (🌀). The Model 775 warming unit will reset/default to the high fan setting each time the warming unit is powered up or *Standby* is selected. The low fan setting (🌀) may be pre-selected while in *Standby* prior to selecting the desired temperature.

STANDBY

Upon power up, the warming unit defaults to *Standby* and the high fan setting. When the warming unit is in 32°C (Low), 38°C (Med), or 43°C (High) or *Ambient* mode, press the **Standby** button to place the warming unit in *Standby* mode. When you select the *Standby* mode, the following events occur:

- **Standby** indicator light illuminates.
- Blower and heater are turned OFF.
- Alphanumeric display deactivates.
- Alarm and fault detection functions remain active.
- Operating timer pauses.
- Airflow mode resets/defaults to the high fan setting.

TEMP IN RANGE INDICATOR LIGHT

The *Temp in Range* indicator light illuminates when the temperature at the blanket or gown end of the hose is $\pm 1.5^{\circ}\text{C}$ of the selected setting; this indicator light does not illuminate in the *ambient* mode.

FAULT INDICATOR LIGHT

When a system fault occurs, the amber *Fault* indicator light flashes and an audible alarm sounds. Please refer to *What to Do in Case of a Fault Condition* on page 7 for additional information.

UNDER-TEMP FAULT

If the warming unit senses an under-temperature condition, the amber *Fault* indicator light flashes and an audible alarm sounds. Please refer to *What to Do in Case of a Fault Condition* on page 7 for additional information.

OVER-TEMP INDICATOR LIGHT

If the warming unit senses an over-temperature condition, the red *Over-temp* indicator light flashes and an alarm sounds. Please refer to *What to Do in Case of an Over-Temperature Condition* on page 7 for additional information.

Mounting the Warming Unit to an IV Pole

To prevent tipping, clamp the Model 775 warming unit to an IV pole at a height that provides stability. Clamp the warming unit no higher than 44" (112 cm) from the floor on an IV pole with a minimum 28" (71 cm) diameter wheelbase.

Caution: If mounted on an IV pole, the distance from the bottom of the warming unit to the floor must be less than 44" (112 cm) and the IV pole wheelbase diameter must be at least 28" (71 cm) to prevent tipping.

METHOD

To mount the Model 775 warming unit to an IV pole, simply turn the clamp handle clockwise to tighten the clamp on the IV pole and counterclockwise to release it. Ensure the IV pole bumper rests against the pole. See Figure D.


Figure D. Model 775 warming unit on an IV Pole

Instructions for Use

The instructions below describe how to operate the Model 775 warming unit. For information about using 3M blankets, and/or gowns with the Model 775 warming unit, refer to the respective Instructions for Use.

1. If the Model 775 warming unit is not mounted on an IV pole or Bair Hugger rolling cart place the warming unit on a flat, hard, dry surface, such as a table, before beginning warming therapy.
Caution: Do not place the warming unit on a soft uneven surface, such as a bed, or a visibly wet surface as the air intake may become blocked and cause the warming unit to overheat, compromising the warming unit's performance.
2. Insert the end of the Bair Hugger warming unit hose in the warming blanket or gown hose port. Use a twisting motion to ensure a snug fit. A visual marker is located around the mid-section of the hose end to guide the depth of hose insertion. (Figures E). Support the hose to ensure secure attachment.
Warning: Do not treat patients with the Bair Hugger warming unit hose alone. Always attach the hose to a Bair Hugger warming blanket/gown before providing warming therapy.


Figure E.

3. Connect the warming unit to a properly grounded power source. The warming unit will be in Standby mode and the Standby indicator light will illuminate. The high fan setting (fan icon) will be pre-selected by default with its indicator light illuminated. The low fan setting (fan icon) may be pre-selected prior to selecting the desired temperature mode.
4. Press the appropriate button (i.e., 32°C, 38°C, 43°C, or Ambient) to select the desired temperature. This will turn the blower and heater ON. When the warming unit reaches the selected temperature, the Temp in Range indicator light will illuminate; this indicator light does not illuminate in the Ambient mode.

Caution: Patient Monitoring Recommendations:

- 3M recommends continuously monitoring core temperature. In the absence of continuous monitoring, monitor the temperature of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Monitor cutaneous responses of patients who are incapable of reacting, communicating and/or who cannot sense temperature a minimum of every 15 minutes or according to institutional protocol.
- Adjust air temperature or discontinue therapy when the therapeutic goal is reached, if elevated temperatures are recorded or if there is an adverse cutaneous response in the warmed area.

5. Press the *Standby* button when warming therapy is complete. Disconnect the warming unit hose from the warming blanket/gown and discard the blanket/gown per hospital policy.

6. Unplug the warming unit to disconnect it from the power source.

Note: At any time, the fan speed setting can be switched between high fan setting (🌀) and low fan setting (🌀) by pressing the corresponding button. However, these buttons will not place the warming unit into *Standby*.

Viewing the Temperature Mode Timer

The temperature mode timer records the amount of time the warming unit has been operating in the selected temperature mode. The timer resets itself every time you select a different temperature mode. To view the temperature mode timer, press the button and hold for 3 seconds for the current temperature mode. The alphanumeric display will show the temperature mode time, then go back to showing the temperature mode setting.

What to Do in Case of an Over-Temperature Condition

If an over-temperature condition occurs, the red *Over-temp* indicator light flashes and an audible alarm sounds and the warming unit automatically turns OFF the heater, blower, and operating mode indicator lights. A fault code appears in the alphanumeric display.

Note: Pressing any button will silence the alarm.

IF AN OVER-TEMPERATURE CONDITION OCCURS:

Warning: Do not continue warming therapy if the red *Over-temp* indicator light illuminates and the alarm sounds. Unplug the warming unit and contact a qualified service technician.

1. Discontinue all warming therapy. If you are using the 241 blood/fluid warming set, immediately stop fluid flow and discard the blood/fluid warming set.
2. Press any button on the control panel to silence the alarm.
3. Unplug the warming unit.
4. Replace warming unit. Do not attempt to restart or reset the warming unit.
5. Contact a qualified service technician.

What to Do If a Fault Condition occurs

The Model 775 warming unit software recognizes several non-hazardous conditions and reports those conditions as faults. When a system fault occurs, the software stores the fault code in memory, the amber *Fault* indicator light flashes, and an alarm sounds. The warming unit automatically turns the heater, blower, and operating mode indicator lights OFF. The control panel will not respond to commands.

Note: Pressing any button will silence the alarm.

IF A FAULT CONDITION OCCURS:

1. Unplug the warming unit and wait for 30-60 seconds.
2. Reconnect the warming unit to a grounded power source. The warming unit will perform the normal power-on-reset sequence and then enter the *Standby* mode.
3. Reselect the temperature setting.
4. If the warming unit does not return to normal operation, contact a qualified service technician.

General Maintenance and Storage

Cleaning Instructions

1. Disconnect the warming unit from the power source before cleaning.
 - Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use; wipe the warming unit, the outside of the warming unit hose, and any other surfaces that may have been touched. Use a damp, soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray.

The following active ingredients are acceptable for use in cleaning the warming unit:

- Oxidizers (e.g. 10% Bleach)
- Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
- Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
- Alcohols (e.g. 70% Isopropyl Alcohol)

2. Let air dry.

Notice:

- Do not use cleaning solutions with greater than 80% alcohol or solvents, including acetone and thinner, to clean the warming unit or hose. Solvents may damage the labels and other plastic parts.
- Do not immerse the Bair Hugger warming unit, warming unit parts, or accessories in any liquid or subject them to any sterilization process.

Storage

Store all components in a cool, dry place when not in use.

Technical Service and Order Placement

USA

1-800-228-3957

OUTSIDE USA

Contact your local service provider or distributor.

When You Call for Technical Support

We will need to know the serial number of your Bair Hugger warming unit. The serial number label is located on the back of the warming unit.

In-Warranty Repair and Exchange

USA

Call 3M customer service if your Model 775 warming unit requires factory service. A customer service representative will give you a Return Authorization (RA) number. Please use this RA number on all correspondence concerning your warming unit. Your customer service representative will also send a shipping carton to you at no charge, if needed. Call your local supplier or sales representative to inquire about borrowing a warming unit while we service your warming unit. For more detailed instructions on returning warming units for service, please refer to the Model 775 Service Manual.

OUTSIDE OF THE USA

Contact your local distributor concerning in-warranty repair and exchange.

PLEASE REPORT A SERIOUS INCIDENT OCCURRING IN RELATION TO THE DEVICE TO 3M AND THE LOCAL COMPETENT AUTHORITY (EU) OR LOCAL REGULATORY AUTHORITY.

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified. Source: ISO 15223, 5.1.7
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
No Free Hosing		Safety sign indicating no free hosing
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.

Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Maximum safe working load		Indicates the maximum safe working load in less than the number reported.
UL Classified		Indicates product was evaluated and Listed by UL for the USA and Canada.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Refer to instruction manual/booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002
Protective earth (ground)		To identify any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode. Source: IEC 60417, 5019
Equipotentiality		To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential. Source: IEC 60417-5021
Fuse		Indicates a replaceable fuse
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
Type BF applied part		To identify a type BF applied part complying with IEC 60601-1. Source: IEC 60417-5333

For more information see, HCBGRegulatory.3M.com

Specifications

Physical Characteristics

DIMENSIONS

13" high x 13" deep x 14" wide (33 cm high x 33 cm deep x 36 cm wide)

WEIGHT

16 lb (7.3kg)

RELATIVE NOISE LEVEL

53 dBA (High fan setting)

48 dBA (Low fan setting)

HOSE

Flexible; Length: 80" (203 cm)

Length: 80" (203 cm)

FILTRATION SYSTEM

MERV 14*

*See Frequently Asked Questions on Filtration (p. 347)

RECOMMENDED FILTER CHANGE

Every 12 months or 500 hours of use.

MOUNTING

Can be clamped to an IV pole, placed on a hard surface, or mounted to the rolling cart accessory.

Temperature Characteristics

RECOMMENDED OPERATING ENVIRONMENT TEMPERATURE

15°C-25°C

HUMIDITY

Max relative humidity 80% (up to 31°C) decreasing linearly to 50% relative humidity at 40°C.

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OPERATING TEMPERATURES

Average temperatures at the end of the hose:

HIGH: 43° ± 1.5°C (109.4° ± 2.7°F)

MED: 38° ± 1.5°C (100.4° ± 2.7°F)

LOW: 32° ± 1.5°C (89.6° ± 2.7°F)

TIME TO REACH OPERATING TEMPERATURE

2 - 5 minutes (dependent on blanket model)

Time required for the contact surface temperature to heat up from 23 ± 2°C to 37°C (73 ± 2°F to 99°F).

STORAGE/TRANSPORT TEMPERATURE

-20°C to 45°C (-4°F to 113°F)

Store all components in a cool, dry place when not in use.

Classification

- Protection against electric shock: Class I Medical Electrical Equipment with Type BF applied part.
- Protection against ingress of water: Ordinary equipment (IPX0)
- Mode of operation: Continuous operation.

Diagnostics

- A qualified service technician can perform over-temperature detection system testing, temperature output testing, operating temperature calibration, and fault code troubleshooting.

Not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

MEDICAL—GENERAL MEDICAL EQUIPMENT AS TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD (2012), CAN/CSA-C22.2 No. 60601-1 (2008) + (2014), IEC 60601-2-35 (2009) + AMD (2016), IEC 60601-1-6:2010 (Third Edition) + A1:2013, IEC 60601-1-8:2008 (Second Edition) + A1:2012; Control No. 4HZ8

Safety System

THERMOSTAT

Independent electronic circuit; thermal cutoff shuts the heater OFF to ensure hose end air remains below 56°C (53°C ± 3°C typical); back-up over-temperature detection at warming unit hose inlet.

ALARM SYSTEM

Over-temperature (≤56°C, 53°C ± 3°C typical): red *Over-Temp* indicator light flashes, alarm sounds, heater and blower shut down, operating indicator lights turn OFF, and control panel becomes unresponsive.

FAULT

Amber *Fault* indicator light flashes, alarm sounds.

OVERCURRENT PROTECTION

Dual input fused lines.

Electrical Characteristics

HEATING ELEMENT

1400 W Resistive

LEAKAGE CURRENT

Meets IEC 60601-1 requirements.

BLOWER MOTOR

Operating speed: 4,700 rpm (high fan setting)

4,100 rpm (low fan setting)

Airflow: up to 48 cfm or 23 L/s

POWER CONSUMPTION

Peak: 1550 W

Average: 800 W

POWER CORD

15 ft., SJT, 3 cond., 13 A

15 ft., SJT, 3 cond., 15 A

4.6 m, HAR, 3 cond., 10 A

DEVICE RATINGS

110-120 VAC, 50/60 Hz, 11.7 A, or

220-240 VAC, 50/60 Hz, 7.2 A, or

100 VAC, 50/60 Hz, 15 A

FUSES

TYPE

Fast acting ceramic fuses, 250 VAC

AMP RATING

12 A (110 - 120 VAC)

8 A (220 - 240 VAC)

15 A (100 VAC)

OPERATING SPEED

F (Fast Acting)

BREAKING CAPACITY

15 A, 12 A: 750 A @ 250 Vac

15 A, 12 A: 10,000 A @ 125 Vac

8 A: 200 A @ 250 Vac

8 A: 10,000 A @ 125 Vac

Guidance and manufacturer's declaration – electromagnetic emissions		
The Model 775 warming unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 775 warming unit should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Model 775 warming unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Model 775 warming unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
The Model 775 warming unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 775 warming unit should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line(s) to earth	±1 kV line to line ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Model 775 warming unit requires continued operation during power mains interruptions, it is recommended that the Model 775 warming unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The Model 775 warming unit is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 775 warming should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	3 Vrms 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the Model 775 warming unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3\sqrt{P} \text{ 800 MHz to 2,5 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 775 warming unit is used exceeds the applicable RF compliance level above, the Model 775 warming unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model 775 warming unit.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the Model 775 warming unit			
The Model 775 warming unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Model 775 warming unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 775 warming unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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3M, Bair Hugger, and the Bair Hugger logo are trademarks of 3M, used under license in Canada.
All rights reserved.
3M Deutschland GmbH