

3M™ Red Dot™ Pediatric Micropore™ Monitoring Electrodes

Commonly Asked Questions

Question:

Why are there so many different kinds of backing that are available in the Red Dot product line?

Answer:

Some are for use in long-term wear situations and others are for short-term situations and still others are for problem adhesion conditions.

- Micropore™.
- Foam.
- Clear Plastic.
- Soft Cloth

Question:

What is the difference(s) between a short-term and long-term electrode?

Answer:

Short-term electrodes:

- Used for a shorter period.
- Less porous/occlusive.
- Good initial tack
- Used for demanding applications.

Long-term electrodes:

- Used for a longer period.
- More porous/breathable.
- Clinically tested to be safe and effective for long-term wear (3 days).

Please note that there have been cases where long-term electrodes have worked well in place of short-term electrodes. The preference to use long- or short-term electrodes depends on customer preference and their clinical knowledge and experience.

Question:

Why is skin abrading recommended prior to attaching electrodes?

Answer:

Skin abrading is not recommended for neonatal and pediatric applications. Skin abrading does reduce skin impedance and expose the inner conductive layer of the epidermis to enhance trace quality. 3M™ Red Dot™ Trace Prep is a "tool" to solve trace quality issues and artifacts (if skin impedance related). Virtually all electrodes face trace quality and skin impedance issues. 3M™ Red Dot™ Trace Prep is designed to minimize those issues. Skin abrading is not only meant for 3M™ Red Dot™ electrodes, but is a good practice to enhance trace quality. 3M has designed this tool to solve customer issues related to trace quality.

**Health Care**

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
USA
1 800 228-3957
www.3M.com/healthcare

For More Information

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.