

# 3M™ Steri-Drape™ Angiography Drapes

## Commonly Asked Questions

**Question:**

Do any of the Steri-Drape angiography drapes contain latex in the drape or packaging?

**Answer:**

3M Steri-Drape angiography drapes do not contain either natural rubber latex or dry natural rubber as components in the product or its packaging.

**Question:**

Do Steri-Drape angiography drapes come sterile and non-sterile?

**Answer:**

Yes. Steri-Drape angiography drapes are packaged sterile for single use and are also sold non-sterile to custom tray vendors for sterilization within a tray.

**Question:**

What is the shelf life? or What is the expiration date?

**Answer:**

The shelf life of the product can be determined by the expiration date which is included on the product's labeling. The word LOT in a box and the hourglass are symbols that represent the lot number and expiration date.

The hourglass is followed by a year and month which represent the expiration date (year and month: 2005-10). When the alpha digits are added, the entire line after the hourglass represents the lot number (2005-10 AZ).

**Question:**

Can an opened but unused drape be resterilized?

**Answer:**

3M does not recommend resterilization of any of our specialty drapes. However, 3M provides a drape replacement program that will exchange opened and unused 3M drapes for sterile drape product. Please contact your 3M Sales Representative for additional information on the drape replacement program.

**Question:**

What is the flammability rating of 3M surgical drapes?

**Answer:**

Steri-Drape surgical drapes classify as having normal flammability (Class 1) using the Consumer Product Safety Commission standard CFR 1610 for the Flammability of Clothing Textiles. The test is the current industry standard, and is used to differentiate materials having normal flammability from those materials which exhibit rapid and intense burning.



**Health Care**

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA  
1 800 228-3957  
[www.3M.com/healthcare](http://www.3M.com/healthcare)

**For More Information**

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.