


3M™ True Definition Scanner

Quick Start Guide for Scanning and TS150™ Mill Workflow


Steps to Complete a Case

 <p>1. Video review of the preparation and contacts.</p>	 <p>2. Scan prepared tooth and contacts.</p>	 <p>3. Choose prep number and review.</p>	 <p>4. Scan remainder of operative arch.</p>	 <p>5. Scan opposing arch.</p>	 <p>6. Scan bite.</p>	 <p>7. Submit prescription.</p>
---	---	--	---	--	--	--

Scan Path


Posterior

1. Scan the occlusal surface from the first premolar to the last molar.
2. Rotate 90° and scan the lingual surface.
3. Rotate over the occlusal surface to scan the buccal.
4. Finish on occlusal.



Anterior

1. Start on the occlusal surface of #28.
2. Scan lingual to #21, concentrating on incisal edges.
3. From #21, rotate occlusal to the buccal.
4. Scan anterior facial surfaces as shown, past midline.
5. Remove wand, return to occlusal of #28, rotate to buccal, scan anterior facial surfaces as shown, past midline. **NOTE:** this may be done in reverse, beginning on #21.



Lab Prescription Form

Enter the following required information:

1. Date
2. Select **TS150™ Chairside Design Station**.
3. Select **Doctor**.
4. Press “**Touch here to select a restoration**”.
5. Select **Restoration** type.
6. Select each prepped tooth by touching tooth once (green).
7. Select **Brand**.
8. Verify Rx for accuracy, then **Sign** (enter PIN).
9. The case will automatically arrive in the “Cases” tab of FastDesign™ CAD/CAM Design Software.*



The screenshot displays the 3M True Definition Scanner's software interface. At the top, the 3M logo and 'True Definition Scanner' are visible, along with a 'Home' button. The main form is titled 'Lab Prescription Form' and contains several sections:

- Date Requested:** A text field containing 'October 17, 2012' with an asterisk indicating it is a required field. A blue circle with the number '1' is next to it.
- Destination:** A dropdown menu showing 'TS150 Chairside Mill' with an asterisk. A blue circle with the number '2' is next to it.
- Doctor:** A text field containing 'Dr. Deckard' with an asterisk. A blue circle with the number '3' is next to it.
- Shipping Address:** A section with an asterisk and the text '* Required field'. It lists: 'Dr. Deckard's Office', '3M Center', 'Building 275', and 'St. Paul, MN 55144'.
- Billing Address:** A section with the same address information as the shipping address.
- Notes:** A large empty text area.
- Buttons:** 'Save as Defaults' and 'Return to Scan'.

At the bottom of the screen, there is a 'Restoration:' section with the instruction: 'Restoration: Touch top line in restoration to select indication.' The device model 'HP Compaq L2105tm' is printed at the very bottom.



Patient: C


Scan Date: 2012-10-24

Status: Scheduled awaiting Scan

ID:

Restoration 1 of 1

- Restoration 4
- 5
- 6
- 7
- 8



Add Restoration
Delete

Sign

Data Plan: **Advanced Plan**

Trusted Workflows

Indications



Materials



Connections



Open Workflows



3M and TS150™ Customer Support



3M™ True Definition Scanner, Cart and Mobile Edition

- Hardware
- Software
- Wireless Connection
- Scan Data/Case Upload
- Scan Protocol

FastDesign™ CAD Station

- Software
- Laptop
- Networking
- Design Protocol
- Indications
- Training

TS150™ Mill and Job Server

- Hardware
- Software
- Networking
- Training
- Restoration Fit

3M True Definition Scanner,
3M Connection Center and
3M Software Questions
1-800-634-2249, Option 3

TS150 Hardware, Software, Application
and Restoration Questions
1-866-791-9538 | CTS@glidewell dental.com

Customer Care Center: 1-800-634-2249
www.3M.com/TrueDef/TS150 for more information



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web www.3M.com/TrueDef

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685



3978 Sorrento Valley Blvd
Suite 200
San Diego, CA 92121
USA
858-202-3360

3M is a trademark of 3M. Used under license in Canada. Please recycle. Printed in U.S.A.

TS150 and FastDesign are registered trademarks of IOS Technologies, Inc.

© 3M 2015. All rights reserved.
70-2013-0601-9 (Rev. 5/2017)