



Evaluation of Four Membrane Filter Materials for Use with 3M™ Petrifilm™ *E. coli* Coliform Count Plates to Enumerate *Escherichia coli* in Water Samples

Robert S. Donofrio¹, Amy Harrison¹, Robin Bechanko¹, DeAnn L. Benesh² and Cynthia Zook²
 1. NSF International, Ann Arbor, MI 48105, 2. 3M, St. Paul, MN 55144



ABSTRACT

Testing was conducted to evaluate the performance of 3M™ Petrifilm™ *E. coli* Coliform Count [EC] Plates in recovering varying concentrations of *Escherichia coli* from water samples. The membrane filtration method specified by the American Public Health Association [SM 9222B] was adapted to incorporate 3M Petrifilm EC Plates as the growth medium. Four commercially available membrane filter types were assessed. All filters showed statistically similar recoveries at a low level of inoculation. A statistically significant difference ($P < 0.05$) was observed for inoculum level in the 100-150 CFU/100 mL range. This demonstrates that although there may be variability amongst filter manufacturers, the 3M Petrifilm EC Plate platform is an adequate culture media for the recovery and enumeration of *E. coli* from water samples.

INTRODUCTION

The American Public Health Association Standard Method 9222B for enumeration of *E. coli* in drinking-water requires filtration of 100 mL water through a membrane filter, followed by plating the filter onto a selective and differential growth medium. 3M Petrifilm EC Plates have global method approvals for use within the food industry, but no method validations for use in drinking or bottled water.

OBJECTIVES

The goals of this project were to:

- Investigate the feasibility of using 3M Petrifilm EC Plates coupled with membrane filters to detect and enumerate *E. coli* in drinking water samples.
- Compare performance of commercially available membrane filters in the above plating application

METHODS

Membrane filters from four manufacturers were evaluated as candidates for use with 3M Petrifilm EC Plates. The following filters (all 0.45 µm pore size and 47 mm in diameter) were utilized for this study:

- Pall GN-6 Metrical Grid
 - Hydrophilic mixed cellulose ester membrane
- Cuno NM04708 BNA045
 - Polyethersulfone membrane
- Whatman Nucleopore
 - Polycarbonate track etch membrane
- Millipore S-PAK Membrane
 - Mixed cellulose

A total of 264 corresponding samples of water were inoculated with *E. coli* ATCC 25922. Sterile phosphate buffered saline (pH 7.0) was utilized as the test water. Half of the samples were inoculated at a “low” concentration of 10 to 50 CFU/100 mL and the other half at a “high” concentration of 100 to 150 CFU/100 mL. Inoculum levels were confirmed via spread plating appropriate dilutions on Tryptic Soy Agar (with incubation under the same conditions as specified in SM 9222B)




Figure 1. Illustration of the work flow process for coupling membrane filtration with 3M™ Petrifilm™ EC Plates. All plates were incubated at 35 ± 0.5 °C for 24 ± 2 hr.




Figure 3. Examples of 3M Petrifilm EC Plates with membrane filters of each of the filter types evaluated. All images were acquired via the 3M™ Petrifilm™ Plate Reader. Difference in media coloration and colonial morphologies can be observed.




Figure 2. Comparison of *E. coli* recovery using membrane filtration coupled to 3M Petrifilm EC Plates. The target for the low inoculum spike was 10 to 50 CFU/100 mL. The target for the high inoculum spike was 100 to 150 CFU/100 mL. Four commercially available membrane filters were evaluated. All filters possessed a pore size of 0.45 µm.

LOW CONCENTRATION <i>E. COLI</i> COMPARISON						
Anova: Single Factor						
SUMMARY						
Groups	Count	Sum	Average	Variance		
Millipore	33	1641	31.54545	28.38068		
Cuno	33	1083	32.81818	38.65591		
Whatman	33	1060	31.81818	31.82841		
Pall	33	1097	33.24242	20.93939		
ANOVA						
Source of Variation	SS	df	MS	F	P-value	F crit
Between Groups	64.20485	3	21.40162	0.769048	0.513386	2.675385
Within Groups	3662.061	128	27.8266			
Total	3626.266	131				
HIGH CONCENTRATION <i>E. COLI</i> COMPARISON						
Anova: Single Factor						
SUMMARY						
Groups	Count	Sum	Average	Variance		
Millipore	33	8349	162.0909	190.2162		
Cuno	33	4706	142.8661	169.1937		
Whatman	33	5231	158.5152	167.7576		
Pall	33	4843	146.7576	113.0644		
ANOVA						
Source of Variation	SS	df	MS	F	P-value	F crit
Between Groups	8548.083	3	2849.361	17.8025	1.02E-09	2.675385
Within Groups	20486.91	128	160.054			
Total	29034.99	131				

Table 1. ANOVA results for filter comparison study. Results from low and high *E. coli* concentration spikes are presented. $\alpha = 0.05$. A significant difference was only observed at the high concentration, and may be attributed to these targets sitting at the edge of the upper limit of the method.

Filter Comparison	T Critical (one tailed)	T Stat Value	P Value ($\alpha = 0.05$)
Millipore vs. Whatman	1.669	1.041	0.15
Pall vs. Cuno	1.669	-1.478	0.07
Cuno vs. Whatman	1.669	-5.038	2.05×10^{-6}
Millipore vs. Cuno	1.669	5.787	1.17×10^{-7}
Pall vs. Whatman	1.669	-4.031	7.51×10^{-5}
Millipore vs. Pall	1.669	4.886	3.61×10^{-6}

Table 2. Student T-test results for filter comparison study. Results from the high *E. coli* concentration spikes are presented. $\alpha = 0.05$.

METHODS

Samples were filtered through each of the four membranes, plated onto hydrated 3M Petrifilm EC plates, incubated and enumerated both manually and using the 3M™ Petrifilm™ Plate Reader (Figure 1). Comparison of the filter recovery was performed using Student *t*-tests and ANOVA ($\alpha = 0.05$).

RESULTS

The results are presented in Figure 2 and 3 as well as Tables 1 and 2. Analysis of data (both pre and post logarithmic transformation) showed counts between the four membranes were not statistically different for low count samples ($P > 0.05$). When comparing all 4 filter recoveries at the high range, a small, but significant difference was observed ($P = 1.40 \times 10^{-9}$). Student *T*-tests revealed that the four filters could be grouped into two statistically similar pairs. Review of images captured by the 3M Petrifilm Plate Reader indicate other factors may potentially affect results: membrane appearance (transparent versus opaque); incubation conditions such as height of plate stacks, location of plate within the stack, and plate moisture loss. The impact of these factors will be investigated in future studies.

SIGNIFICANCE

- The differences observed between filters at the high inoculum spike level may be due to the upper limitations of the method and not necessarily to the filters themselves.
- Preferred counting range for 3M Petrifilm EC Plates is 15-150 CFU per plate³
- These preliminary studies indicate that the 3M Petrifilm EC Plate coupled with membrane filtration shows promise for consideration as an acceptable alternative quantitative *E. coli* method to SM 9222B.

REFERENCES

- American Public Health Association. 1998. Standard Methods for the Examination of Water and Wastewater. 20th edition, Washington, DC.
- USEPA. 2004. EPA 821-B-03-004 - EPA Microbiological Alternate Test Procedure (ATP) Protocol for Drinking Water, Ambient Water, and Wastewater Monitoring Methods.
- American Public Health Association. 2001. Compendium of Methods for the Microbiological Examination of Foods. 4th edition. Washington, D.C.

CONTACT INFORMATION

Robert S. Donofrio, M.S., Ph.D.
 Director of Microbiology
 NSF International
 789 Dixboro Rd., Ann Arbor, MI
 donofrio@nsf.org
 Office - 734-827-6894
 www.nsf.org

DeAnn Benesh
 Regulatory Affairs
 3M Food Safety Department
 3M Center; Bldg 260-8B-01
 St. Paul, MN 55144-1000
 Office: 651-736-3594
 dbenesh1@mm.com