- (EN) ***Salmonella* Enrichment Supplement**
- (FR) ***Salmonella* Supplément pour Enrichissement**
- (DE) ***Salmonella* Anreicherung - Zusatz**
- (IT) **Supplemento per l'arricchimento di *Salmonella***
- (ES) **Suplemento para el Enriquecimiento de *Salmonella***
- (NL) ***Salmonella* Verrijking Supplement**
- (SV) ***Salmonella* anrikningsbuljong tillägg**
- (DA) ***Salmonella* Opformerings supplement**
- (NO) ***Salmonella* supplement til oppformeringsbuljong**
- (FI) ***Salmonella* Rikastussupplementti**
- (PT) **Suplemento para Enriquecimento de *Salmonella***
- (EL) **Σαλμονέλλα, Συμπλήρωμα Εμπλουτισμού**
- (PL) ***Salmonella* Dodatek wzbogacający**
- (HU) ***Salmonelladúsító* kiegészítő**
- (CS) **Obohacující supplement pro stanovení *Salmonelly***
- (RO) **Supliment de îmbogățire pentru *Salmonella***
- (RU) **Добавка к среде обогащения для бактерий рода *Salmonella***
- (TR) ***Salmonella* Zenginleştirme Saplementi**
- (JA) **サルモネラエンリッチメントサプリメント**
- (ZH) **沙门氏菌增菌补充物**
- (TH) **ซัพพลายเม้นท์สำหรับอาหารเลี้ยงเชื้อซาลโมเนลล่า**
- (KO) **살모넬라 증균 보충제**

Product Instructions

Salmonella Enrichment Supplement

DESCRIPTION

For use with 3M™ *Salmonella* Enrichment Base in conjunction with the 3M™ Petrifilm™ *Salmonella* Express (SALX) System.

SAFETY

The user should read, understand, and follow all safety information in the instructions for the 3M *Salmonella* Enrichment Supplement. Retain the safety instructions for future reference.

⚠ WARNING Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

NOTICE Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

⚠ WARNING

To reduce the risks associated with false-negative results leading to release of contaminated food or beverage product:

- Always follow the proper enrichment scheme (time, temperature, and dilution).
- Always use care when pipetting or transferring samples.
- Follow the protocol and perform the tests exactly as stated in the product instructions for use.

To reduce the risks associated with false negative or false positive results leading to the release of contaminated product or the destruction of good product:

- Always follow product instructions for preparation, handling, and storage.

To reduce the risks associated with exposure to chemicals and biohazards:

- Dispose of enriched samples according to all applicable government regulations and applicable laboratory procedures.
- Read, understand, and follow the safety information in the 3M Petrifilm *Salmonella* Express System Instructions for Use and the associated MSDS.
- Always follow standard good laboratory safety practices (GLP)¹, including proper containment procedures, wearing appropriate protective apparel and eye protection while handling testing materials and test samples.

To reduce the risks associated with exposure to chemicals which, if not avoided, may result in allergic skin reaction, allergic respiratory reaction, target organ effects, genotoxic or mutagenic effects and may interact with genetic material and possibly alter gene expression:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using this product.
- If eye contact does occur, flush eyes with large amounts of water. If signs/symptoms persist, get medical attention.
- If skin contact occurs, remove contaminated clothing and shoes. Immediately flush skin with large amounts of water. Get medical attention, wash contaminated clothing, and clean shoes before reuse.
- Do not eat, drink, or smoke when using this product.
- Always wash exposed areas thoroughly with soap and water.
- If inhaled, remove person to fresh air. If signs/symptoms develop, get medical attention.
- If swallowed, do not induce vomiting unless instructed to do so by medical personnel. Give victim two glasses of water. Never give anything by mouth to an unconscious person. Get medical attention.

To reduce the risks associated with environmental contamination:

- Follow current industry standards and local regulations for disposal of contaminated waste.

NOTICE: To reduce false-positive results due to cross contamination that result in re-testing or the rejection of food or beverage product:

- Always wear gloves.
- Use a new pipette tip for each sample transfer.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

LIMITATION OF 3M LIABILITY

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE AND DISPOSAL

Storage: 2-8°C. Material is hygroscopic. Keep container tightly closed. Protect from light, heat and moisture.

Disposal: After use, enriched samples may contain microorganisms that may be a potential biohazard.

Follow current industry standards and local regulations for disposal of contaminated waste. Consult the Material Safety Data Sheet for additional information.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory safety practices (GLP)¹.

Use an analytical balance capable of reading at least three to four decimal places to weigh the supplement.

Preparation of Supplement: Allow vial to come to ambient temperature. Gently invert the vial 3-4 times to mix the supplement powder. Aseptically add 50 mg (\pm 2.5 mg) of the 3M™ *Salmonella* Enrichment Supplement to 1 L of sterilized 3M™ *Salmonella* Enrichment Base tempered to 20 to 42.5°C. Mix thoroughly.

Opened Vial Shelf Life: Once the vial has been opened, it is acceptable to use for 5 months as long as the expiration date has not been exceeded.

Composition: (typical formula g/L):
Selective Agents 0.05 g

REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.

EXPLANATION OF SYMBOLS

Consult instructions for use

The lot in a box and the hourglass symbols are symbols that represent lot number and expiration date. The hourglass is followed by a year, month, and day which represents the expiration date (year, month, and day: YYYY- MM-DD). The entire line after the LOT represents the lot number.

Storage conditions

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz - Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Asia Pacific

No 1, Yishun Avenue 7
Singapore, 768923
65-64508869

3M Japan

3M Health Care Limited
33-1, Tamagawadai 2-chrome
Setagaya-ku, Tokyo
158-8583, Japan
81-570-011-321

3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Health Care

2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

© 2013, 3M. All rights reserved.
3M and Petrifilm are trademarks of 3M. Used under license in Canada.
34-8711-0869-1

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