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Meta-analysis and trial sequential analysis of prophylactic negative pressure therapy for groin wounds in vascular surgery.¹⁷

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
Meta-analysis of negative pressure wound therapy of closed groin incisions in arterial surgery.¹⁸

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
Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study.¹⁹


Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Boning A. *Int Wound J*. 2018;15(1):75-83.

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Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy.²⁰

Grauhan O, Navasardyan A, Hofmann M, et al. *J Thorac Cardiovasc Surg.* 2013;145:1387–1392.

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

Sternal wound debridement following a three-vessel CABG.

Tanna N, MD, FACS, Northwell Health, New York, USA.

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
Niveau van klinisch bewijs.

Geciteerde bron	Type wond/operatie	Niveau van klinisch bewijs*
Engelhardt M, Rashad NA, Willy C, Müller C, Bauer C, Debus S, Beck T. Closed-incision negative pressure therapy to reduce groin wound infections in vascular surgery: a randomised controlled trial. <i>Int Wound J.</i> 2018;doi: 10.1111/iwj.12848.	Lieswonden	1b ●
Gombert A, Babilon M, Barbati ME, Keskei A, et al. Closed incision negative pressure therapy reduces surgical site infections in vascular surgery: a prospective randomised trial (AIMS Trial). <i>Eur J Vasc Endovasc Surg.</i> 2018;doi:10.1016/j.ejvs.2018.05.018.	Lieswonden	1b ●
Kwon J, Staley C, McCullough M, Goss, S, et al. A randomized clinical trial evaluating negative pressure therapy to decrease vascular groin incision complications. <i>J Vasc Surg.</i> 2018;1-9.	Lieswonden	1b ●
Lee AJ, Sheppard CE, Kent WD, Mewhort H, Sikdar KC, Fedak PW. Safety and efficacy of prophylactic negative pressure wound therapy following open saphenous vein harvest in cardiac surgery: a feasibility study. <i>Interact Cardiovasc Thorac Surg.</i> 2017;24(3):324–328.	Open oogsten van vena saphena bij hartchirurgie	1b ●
Lee K, Murphy PB, Ingves MV, Duncan A, DeRose G, Dubois L, Forbes TL, Power A. Randomized clinical trial of negative pressure wound therapy for high-risk groin wounds in lower extremity revascularization. <i>J Vasc Surg.</i> 2017;66(6):1814–19.	Lieswonden	1b ●
Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Böning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. <i>Int Wound J.</i> 2018;15(1):75–83. doi: 10.1111/iwj.12836.	Incisies in de lies	1b ●
Colli A. First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients. <i>J Cardiothoracic Surg.</i> 2011;6(1):160.	Sternotomie	2 ●
Weir G. The use of a surgical incision management system on vascular surgery incisions: a pilot study. <i>Int Wound J.</i> 2014;11 Suppl 1:10–2.	Vasculaire bypass	2 ●
Grauhan O, et al. Effect of surgical incision management on wound infections in a poststernotomy patient population. <i>Int Wound J.</i> 2014;11 Suppl 1:6–9.	Sternotomie	2b ●
Grauhan O, et al. Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. <i>J Thorac Cardiovasc Surg.</i> 2013;145(5):1387–92.	Sternotomie	2b ●
Matatov T, et al. Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients. <i>J Vasc Surg.</i> 2013;57(3):791–5.	Vasculaire bypass	3 ●
Santarpino G, Gazdag L, Sirch J, Vogt F, Ledwon M, Fischlein T, Pfeiffer S. A Retrospective Study to Evaluate Use of Negative Pressure Wound Therapy in Patients Undergoing Bilateral Internal Thoracic Artery Grafting. <i>Ostomy Wound Manage.</i> 2015 ec;61(12):26–30.	Thoracale-arterietransplantatie	3 ●
Simon K, et al. [Use of Negative Pressure Wound Therapy on Surgical Incisions (Prevana™) after Surgery of Pectus Deformities Reduces Wound Complications.]. <i>Zentralblatt für Chirurgie.</i> 2014. [Duitstalig]	Sternotomie	3 ●
Atkins BZ, Tetterton JK, Petersen RP, et al. Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications? <i>Surg Innov.</i> 2009;16(2):140–6.	Sternotomie	4 ●
Atkins BZ, et al. Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy. <i>Int Wound J.</i> 2011;8(1):56–62.	Sternotomie	4 ●
Reddy VS. Use of Closed Incision Management with Negative Pressure Therapy for Complex Cardiac Patients. <i>Cureus.</i> 2016;8(2):e506.	Sternotomie	4 ●

Geciteerde bron	Type wond/operatie	Niveau van klinisch bewijs*	
Simon K, <i>et al.</i> [Use of Negative Pressure Wound Therapy on Surgical Incisions (Prevena™) after Surgery of Pectus Deformities Reduces Wound Complications.]. <i>Zentralblatt für Chirurgie</i> . 2014. [Duitstalig]	Sternotomie	3	●
Chopra K, Tadisina KK, Singh DP. The 'French Fry' VAC Technique: Hybridization of Traditional Open Wound NPWT with Closed Incision NPWT. <i>Int Wound J</i> . 2016;13(2):216-9.	Massief gelokaliseerd lymfoedeem	5	●
Dohmen PM, <i>et al.</i> Can post-sternotomy mediastinitis be prevented by a closed incision management system? <i>GMS Hyg Infect Control</i> . 2014;9(3):Doc19.	Sternotomie	5	●
Dohmen PM, <i>et al.</i> Use of incisional negative pressure wound therapy on closed median sternal incisions after cardiothoracic surgery: clinical evidence and consensus recommendations. <i>Med Sci Monit</i> . 2014;20:1814-25.	Sternotomie	5	●
Haghshenas Kashani A, Varcoe RL. A New Negative Pressure Dressing (Prevena) to Prevent Wound Complications Following Lower Limb Distal Arterial Bypass. <i>Br J Diabetes Vasc Dis</i> . 2011;11(1):21-4.	Vasculaire bypass	5	●
Wu RT, Sumpio BJ, Miller S, Sumpio BE. Use of closed-incision negative pressure therapy: cardiothoracic and vascular surgery. <i>Plast Reconstr Surg</i> . 2019;143:31S.	Sternotomie	5	●

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Literatuur

- McGillicuddy E, Ozaki K, *et al.* The impact of vascular surgery wound complications on quality of life. *Journal of Vascular Surgery*. Volume 64, nr. 6, december 2016, pag. 1780-1788.
- Fowler VG Jr, O'Brien SM, Muhlbaier LH, *et al.* Clinical predictors of major infections after cardiac surgery. *Circulation*. 2005;112 (9 Suppl.):1358-65.
- Willy C, Engelhardt M, *et al.* The impact of surgical site occurrences and the role of closed incision negative pressure therapy. *International Wound Journal*. 2016, pag. 35-46.
- Giles KA, Hamdan AD, Pomposelli FB, Wyers MC, Siracuse JJ, Schermerhorn ML. Body mass index: surgical site infections and mortality after lower extremity bypass from the national surgical quality improvement program 2005-2007. *Ann Vasc Surg*. 2007;24:48-56.
- Davenport DL, Zwischenberger BA, Xenos ES. Analysis of 30-day readmission after aortoiliac and infrainguinal revascularization using the American College of Surgeons national surgical quality improvement program data set. *J Vasc Surg*. 2011;60:1266-74.
- Kuy S, Dua A, Desai S, Dua A, Patel B, Tondravi N, *et al.* Surgical site infections after lower extremity revascularization procedures involving groin incisions. *Ann Vasc Surg*. 2014;28:53-8.
- Ozaki CK, Hamdan AD, Barshes NR, Wyers M, Hevelone ND, Belkin M, *et al.* Prospective, randomized, multi-institutional clinical trial of a silver alginate dressing to reduce lower extremity vascular surgery wound complications. *J Vasc Surg*. 2015;61:419-27.e1.
- Wiseman JT, Fernandes-Taylor S, Barnes ML, Saunders RS, Saha S, Havlena J, *et al.* Predictors of surgical site infection after hospital discharge in patients undergoing major vascular surgery. *J Vasc Surg*. 2012;62:1023-31.e5.
- Greenblatt D, Rajamanickam V, *et al.* Predictors of surgical site infection after open lower extremity revascularization. *Journal of Vascular Surgery*. 2011; 54(2): 433-439.
- Graf K, Ott E, Vonberg RP, *et al.* Surgical site infections e economic consequences for the health care system [Review] *Langenbecks Archs Surg* 2011;396:453e9.
- Coello R, Charlett A, Wilson J, *et al.* Adverse impact of surgical site infections in English hospitals. *J Hosp Infect* 2005;60:93e103.
- Jenks PJ, Laurent M, McQuarry S, Watkins R. Clinical and economic burden of surgical site infection (SSI) and predicted financial consequences of elimination of SSI from an English hospital. *J Hosp Infect* 2014;86:24e33.
- Wilkes RP, Kilpadi DV, Zhao Y, *et al.* Closed incision management with negative pressure wound therapy (CIM): Biomechanics. *Surg Innov*. 2012;19(1):67-75.
- Kilpadi DV, Cunningham MR. Evaluation of Closed Incision Management with Negative Pressure Wound Therapy (CIM): Hematoma/Seroma and Involvement of the Lymphatic System. *Wound Repair Regen*. 2011;19:588-596.
- Glaser DA, Farnsworth CL, Varley ES, *et al.* Negative pressure therapy for closed spine incisions: A pilot study. *Wounds*. 2012;24(11):308-316.
- Federal Drug Administration. De Novo Classification Request for PREVENA 125 and PREVENA PLUS 125 Therapy Units. De Novo Summary (DEN180013), 2019. https://www.accessdata.fda.gov/cdrh_docs/reviews/DEN180013.pdf
- Antoniou GA, Onwuka CC, Antoniou SA, Russell D. Meta-analysis and trial sequential analysis of prophylactic negative pressure therapy for groin wounds in vascular surgery. *J Vasc Surg*. 2019; doi 10.016/jvs.2019.01.083.
- Svensson-Björk R, Zarrouk M, Ascitto G, Hasselmann J, Acosta S. Meta-analysis of negative pressure wound therapy of closed groin incisions in arterial surgery. *BJS*. 2018; doi:10.1002/bjs.11100.
- Pleger SP, Nink N, Elzien M, Kunold A, Koshty A, Boning A. Reduction of groin wound complications in vascular surgery patients using closed incision negative pressure therapy (ciNPT): a prospective, randomised, single-institution study. *Int Wound J*. 2018;15(1):75-83.
- Grauhan O, Navasardyan A, Hofmann M, *et al.* Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy. *J Thorac Cardiovasc Surg*. 2013;145:1387-1392.
- Willy, C., Agarwal, A., Andersen, C.A., Santis, G.D., Gabriel, A., Grauhan, O., Guerra, O.M., Lipsky, B.A., Malas, M.B., Mathiesen, L.L., Singh, D.P. and Reddy, V.S. (2017), Closed incision negative pressure therapy: international multidisciplinary consensus recommendations. *Int Wound J*, 14: 385-398. doi:10.1111/iwj.12612.

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© 2020 3M. Alle rechten voorbehouden. 3M en de andere weergegeven merken zijn merken of geregistreerde merken. Ongeoorloofd gebruik verboden. PRA-PM-NL-00062 (07/20). OMG102183.