



Littmann®
Stethoscopes

Classic

Love your stethoscope.

The compact and sensitive **3M™ Littmann® Classic III™ Monitoring Stethoscope** is the perfect instrument for clinicians who monitor and assess children or adults in non-critical care environments. Manufactured with strong yet lightweight materials for hours of comfortable use.

Monitor and assess a wide range of patients.

- Useful in non-critical care environments such as a medical office, general care unit, OB/GYN, ambulatory clinic or urgent care
- Fun and vibrant; match color and finish to your personality
- Designed for use with adult and pediatric patients
- Tunable, dual-sided stainless steel chestpiece with open or closed bell
- Non-chill bell sleeve for greater patient comfort
- Soft-sealing eartips provide optimal comfort and excellent sound occlusion
- Made in the USA without natural rubber latex



Product Features

Acoustic Performance	8 (1–10 scale)
Chestpiece Design/Technology	Machined stainless steel, two-sided
Headset Material	Lightweight, aircraft-grade, aluminum alloy
Approximate Length	27"/69 cm
Approximate Weight	150 grams
Diaphragm Details/Diameter	Two tunable diaphragms 1.7"/4.3 cm, 1.4"/3.4cm
Warranty	5 years



3M Health Care
3M Medical Solutions Division
2510 Conway Ave., St. Paul, MN 55144 USA
1-800-228-3957 | Littmann.com

Explore Littmann stethoscopes at Littmann.com or call the 3M Health Care Helpline at 1-800-228-3957.

3M, Littmann, the L Littmann logo and Classic III are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A.
© 3M 2019. All rights reserved. 70-2011-7198-3. Rev. 05/2019

3M™ Littmann® Classic III™ Monitoring Stethoscope



Standard stainless steel



5620 Black 5621 Gray 5622 Navy 5623 Caribbean 5630 Ceil Blue



5633 Pearl Pink 5627 Burgundy 5831 Plum 5832 Lavender 5839 Lemon-Lime



5835 Turquoise

Smoke/smoke



5811 Black 5812 Dark Olive 5872 Turquoise 5873 Gray 5875 Lime Green




5871 Raspberry

Mirror/smoke



5862 Raspberry 5863 Navy 5865 Lavender

Rainbow/brass finish



5806 Raspberry 5807 Caribbean

Copper finish



5809 Chocolate

Rainbow/black



5870 Black

Champagne/smoke



5861 Black 5864 Burgundy

Mirror/standard



5874 Caribbean Blue

Black/black



5803 Black 5867 Royal Blue 5868 Burgundy 5869 Caribbean



3M Health Care
 3M Medical Solutions Division
 2510 Conway Ave., St. Paul, MN 55144 USA
 1-800-228-3957 | Littmann.com

Product Features	
Acoustic Performance	7 (1-10 scale)
Chestpiece Design/Technology	Machined stainless steel, two-sided
Headset Material	Lightweight, aircraft-grade aluminum alloy
Length/Total Weight	27"/69 cm, 150 grams
Warranty	5 years
Diaphragm Details/Diameter	Two, single-piece tunable diaphragms 1.7"/4,3 cm, 1.4"/3,4 cm

Colors shown are approximate.