

3M CASTING & SPLINTING PRODUCTS

3M has a broad line of synthetic casting tapes, splinting products and compatible casting accessories. Conformable, customizable and comfortable, 3M™ Scotchcast™ Splints provide strength and stability. Our ready-to-use 3M™ Scotchcast™ Conformable Splints are made from layered fiberglass with pre-padded, water-repellent felt to ensure quick drying. In addition, 3M™ Primacast™ Unpadded Splint technology offers the speed and strength users expect from synthetic splints, combined with conformability.



3M™ Scotchcast™ Plus Casting Tape


A lightweight, strong and durable casting tape that combines the benefits of a fiberglass casting tape with the handling ease of plaster. Available in 13 colors to meet patient preferences.

Product Benefits:

- Water-activated resin lets you control working time according to application needs
- High and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs or replacements
- Lightweight for easy mobility
- Stretches to conform well to contours for a more comfortable, better fitting cast
- Porous weave allows the skin to breathe

Suggested Applications:

- For use in the treatment of musculoskeletal injuries



Available Colors:

- Black (A)
- Green (G)
- Purple (U)
- Bright Green (V)
- Maroon (M)
- White
- Light Blue (L)
- Bright Orange (W)
- Glow-in-the-Dark (J)
- Blue (B)
- Red (R)
- Bright Pink (X)
- Yellow (Y)

Cat. No.	Description	Width	Length	Rolls/Case
82001	White	1" (2,5 cm)	2 yds. (1,8 m)	10
8200*	White	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*A	Black	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8200*B	Blue	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*G	Green	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*J	Glow-in-the-Dark	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8200*L	Light Blue	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*R	Red	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*U	Purple	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8200*V	Bright Green	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8200*W	Bright Orange	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8200*X	Bright Pink	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8200*Y	Yellow	2", 3" (5 cm, 7,5 cm)	4 yds. (3,6 m)	10
8200*M	Maroon	2", 3" (5 cm, 7,5 cm)	4 yds. (3,6 m)	10





3M™ Scotchcast™ Wet or Dry Cast Padding

For use with 3M™ Synthetic Cast Stockinet and 3M™ Scotchcast™ Plus Casting Tape.

A synthetic cast padding that allows patients to get their cast wet, when following the instructions for use.

Product Benefits:

- Allows patients to shower, bathe, swim or use a whirlpool
- Sheds water and feels comfortably dry to patient in approximately 1 to 3 hours
- Allows you to stock only one cast padding for all your patient needs
- Suitable for use with all synthetic casting materials such as Scotchcast Plus casting tape

Suggested Applications:

- For use in the treatment of musculoskeletal injuries



Cat. No.	Width x Length	Rolls/Bag	Bags/Case
WDP2	2" x 4 yds. (5 cm x 3,6 m)	20	4
WDP3	3" x 4 yds. (7,5 cm x 3,6 m)	20	4
WDP4	4" x 4 yds. (10 cm x 3,6 m)	20	4
WDP6	6" x 4 yds. (15 cm x 3,6 m)	20	4

3M™ Scotchcast™ Soft Cast Casting Tape

A fiberglass casting tape that does not set to a completely rigid state, but remains flexible and soft even when it is completely cured.

Product Benefits:

- Flexible, soft cast allows for mobility
- The soft, pliable cast edges improve patient comfort
- Can be used in combination with rigid fiberglass splints, as soft edges provide additional comfort

Suggested Applications:

- For use in the construction of casts where support is needed but rigid immobilization is not required, or a limited range of motion is desired

Applications Include:

- Sports Medicine: thumb, wrist and ankle sprains
- Pediatric Orthopedics: serial casting for club foot treatment
- General Orthopedics: secondary casting, functional stabilization casting
- Occupational Therapy: removable splints
- Diabetic Foot: total contact casting



Available Colors:

- White
- Blue (B)
- Purple (U)
- Red (R)

Cat. No.	Description	Width	Length	Rolls/Case
82101	White	1" (2,5 cm)	2 yds. (1,8 m)	10
8210*	White	2", 3", 4", 5" (5 cm, 7,5 cm, 10 cm, 12,5 cm)	4 yds. (3,6 m)	10
8210*B	Blue	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8210*U	Purple	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10
8210*R	Red	2", 3", 4" (5 cm, 7,5 cm, 10 cm)	4 yds. (3,6 m)	10



3M™ Synthetic Cast Padding

A non-woven polyester padding that is porous, highly conformable and strong, yet easy to tear. Suitable for use with all casting materials, plaster and synthetic.



Cat. No.	Size	Rolls/Bag	Bags/Case
CMW02	2" x 4 yds. (5 cm x 3,6 m)	20	4
CMW03	3" x 4 yds. (7,5 cm x 3,6 m)	20	4
CMW04	4" x 4 yds. (10 cm x 3,6 m)	20	4
CMW06	6" x 4 yds. (15 cm x 3,6 m)	20	4

3M™ Synthetic Cast Stockinet

Essential component of a high quality cast, this synthetic stockinet is intended for use as an underlayment for all standard casting applications. Suitable for use with all casting materials, plaster and synthetic.



Cat. No.	Size	Rolls/Case
MS01	1" x 25 yds. (2,5 cm x 22,8 m)	1
MS02	2" x 25 yds. (5 cm x 22,8 m)	1
MS03	3" x 25 yds. (7,5 cm x 22,8 m)	1
MS04	4" x 25 yds. (10,1 cm x 22,8 m)	1
MS06	6" x 25 yds. (15,2 cm x 22,8 m)	1
MS08	8" x 25 yds. (20,3 cm x 22,8 m)	1

3M™ Scotchcast™ Conformable Roll Splint

Cut-to-length fiberglass splinting material is packaged in an easy-to-use, easy-to-seal system. The splinting material is made of seven layers of fiberglass and is pre-padded with our proprietary water-repellent felt.

Product Benefits:

- Splint can be customized and cut to the exact length needed for less waste and cost-effective product use
- Easy-to-use delivery system helps ensure freshness of splinting material and reduced waste
- Provides strong, lightweight support in minutes to encourage early patient mobility
- Less messy than plaster splints for easy application and faster cleanup
- Water-repellent felt padding dries much more quickly than standard padding
- Water-repellent padding promotes ease of use and faster application

Suggested Applications:

- For use in the construction of common orthopedic splints
- Ideal for high volume users such as emergency departments



Seven layers of fiberglass have proven over time to provide the strength you need.

Cat. No.	Width x Length	Packaging
73002	2" x 15 yds. roll (5,1 cm x 13,7 m)	1 roll/box, 2 boxes/case
73003	3" x 15 yds. roll (7,6 cm x 13,7 m)	1 roll/box, 2 boxes/case
73004	4" x 15 yds. roll (10,1 cm x 13,7 m)	1 roll/box, 2 boxes/case
73005	5" x 15 yds. roll (12,7 cm x 13,7 m)	1 roll/box, 2 boxes/case

3M™ Primacast™ Unpadded Splint

A unique splint technology that combines the speed and strength users expect from synthetic splints with the conformability and feel of plaster. Available in a variety of sizes to meet patient treatment needs.

Product Benefits:

- Water-activated material sets in approximately 3–4 minutes with adequate functional strength
- Provides a plaster-like conformability for excellent molding and patient fit
- Provides superior radiolucency for easier diagnosis from x-rays
- Unpadded for easy customization

Suggested Applications:

- For use in the construction of common orthopedic splints
- Ideal for use in operating room applications



Cat. No.	Width x Length	Packaging
30210	2" x 10" (5,1 cm x 25,4 cm)	15 each/case
30312	3" x 12" (7,6 cm x 30,5 cm)	15 each/case
30335	3" x 35" (7,6 cm x 88,9 cm)	15 each/case
30415	4" x 15" (10,1 cm x 38,1 cm)	15 each/case
30430	4" x 30" (10,1 cm x 76,2 cm)	15 each/case
30530	5" x 30" (12,7 cm x 76,2 cm)	10 each/case
30545	5" x 45" (12,7 cm x 114,3 cm)	10 each/case
30630	6" x 30" (15,2 cm x 76,2 cm)	10 each/case
30645	6" x 45" (15,2 cm x 114,3 cm)	10 each/case
30830	8" x 30" (20,3 cm x 76,2 cm)	10 each/case
30845	8" x 45" (20,3 cm x 114,3 cm)	10 each/case



Pre-cut Splints

3M™ Scotchcast™ Conformable Splints

A ready-to-use, highly conformable splint made of layered fiberglass and pre-padded with water-repellent felt for quicker drying. Prepackaged in a variety of sizes to meet application needs.

Product Benefits:

- Low profile, one-sided felt padding is less bulky, so it molds more easily to the patient for a closer, more conformable fit
- All-in-one splint allows for easy application and time savings, resulting in faster application and patient turnaround time
- Water-repellent felt padding dries much more quickly than standard padding
- Provides strong, lightweight support in minutes
- Less messy than plaster splints for easy application and faster cleanup

Suggested Applications:

- For use in the construction of common orthopedic splints



Seven layers of fiberglass have proven over time to provide the strength you need.

Cat. No.	Width x Length	Packaging
72210	2" x 10" (5,1 cm x 25,4 cm)	10 each/case
72312	3" x 12" (7,6 cm x 30,5 cm)	10 each/case
72335	3" x 35" (7,6 cm x 88,9 cm)	10 each/case
72415	4" x 15" (10,1 cm x 38,1 cm)	10 each/case
72430	4" x 30" (10,1 cm x 76,2 cm)	10 each/case
72530	5" x 30" (12,7 cm x 76,2 cm)	10 each/case

3M™ Scotchcast™ One-Step Splint

A splint that incorporates foam padding, splint and backing for easy, one-step application. Prepackaged in a variety of sizes to meet application needs.

Product Benefits:

- All-in-one splint allows for easy application and time savings
- Conformable and moldable for better patient fit
- Provides strong, lightweight support to encourage early patient mobility
- Available in a variety of pre-cut sizes to accommodate a variety of patient sizes and reduce the need to cut and customize splints

Suggested Applications:

- For use in the construction of common orthopedic splints



Cat. No.	Width x Length	Packaging
76210A	2" x 10" (5 cm x 25 cm)	10 each/case
76312A	3" x 12" (7,5 cm x 30 cm)	10 each/case
76415A	4" x 15" (10 cm x 37,5 cm)	10 each/case
76430	4" x 30" (10 cm x 75 cm)	5 each/case
76530	5" x 30" (12,5 cm x 75 cm)	5 each/case
86335	3" x 35" (7,5 cm x 87,5 cm)	5 each/case
86545	5" x 45" (12,5 cm x 112,5 cm)	5 each/case

For more information

For additional information contact your 3M Health Care Sales Representative, call the 3M Health Care Customer Helpline at 1-800-228-3957, or visit 3M.com/Immobilization.



3M Infection Prevention Division 3M Health Care

2510 Conway Avenue
St. Paul, MN 55144-1000
USA
1-800-228-3597
www.3M.com/healthcare
www.futuro-usa.com

3M, Primacast and Scotchcast
are trademarks of 3M.
Please recycle. Printed in U.S.A.
© 3M 2014. All rights reserved.
70-2009-7709-1-A