

## Forsus™ Fatigue Resistant Device with Direct Push Rod





3M Unitek, an industry-leading innovator, is pushing the limitations of Class II Correction and patient satisfaction. Until now, your options were limited to fixed intraoral devices with breakage concerns, or other appliances, including headgear, requiring a high degree of patient compliance. That's why we created a more reliable appliance for Class II Correction—a 3-piece telescoping spring solution that fits perfectly with your current appliance system and is more attractive and easier for your patients to live with. The design incorporates an external spring assembly that won't compromise patient comfort and acceptance, a key to achieving successful clinical outcomes when using compliance-independent appliances. With the addition of the Direct Push Rod, installation can now be achieved in minutes. The Direct Push Rod attaches directly onto the archwire, eliminating the need for auxiliaries and simplifying chairside installation.

## Fatigue resistant by design.

The unique coaxial spring design of the Forsus<sup>™</sup> Fatigue Resistant Device addresses the issue of fatigue failure—a fracture caused by repeated application of stresses in the coil spring. To maximize fatigue life, the Forsus Fatigue Resistant Device is designed to keep these stresses to a minimum. And its superior strength and consistent forces have been shown in laboratory data.

Traditional solutions simply can't compare in terms of strength, reliability, ease of use and patient satisfaction. Put the Forsus Fatigue Resistant Device to work for you, and discover how it can force your productivity to new heights.










## Forsus™ Fatigue Resistant Device with Direct Push Rod

Description	Unit		Part No.
Forsus Fatigue Resistant Device 5-Patient Kit	Kit		885-121
Forsus Fatigue Resistant Device 20-Patient Kit	Kit		885-122
Universal Spring Module	pkg/5		885-100
Direct Push Rod Short (25mm)	pkg/5	L	885-111
Direct Push Rod Short (25mm)	pkg/5	R	885-112
Direct Push Rod Med (29mm)	pkg/5	L	885-113
Direct Push Rod Med (29mm)	pkg/5	R	885-114
Direct Push Rod Lg (32mm)	pkg/5	L	885-115
Direct Push Rod Lg (32mm)	pkg/5	R	885-116
Direct Push Rod XL (35mm)	pkg/5	L	885-117
Direct Push Rod XL (35mm)	pkg/5	R	885-118
Direct Push Rod XXL (38mm) no stop	pkg/5	L	885-107
Direct Push Rod XXL (38mm) no stop	pkg/5	R	885-108
Universal L Pin	pkg/10		885-110
Universal Split Crimp	pkg/50		885-203
Measurement Guide	pkg/5		807-014

Pushing the limitations of Class II Correctors.



A 3-piece telescoping coaxial spring designed to allow for enhanced Easy to install patient jaw opening L-Pin Connector



Direct Push Rod attaches directly onto the archwire



3M Unitek **Orthodontic Products** 

2724 South Peck Road Monrovia, CA 91016 USA

www.3MUnitek.com

In U.S. and Puerto Rico: 1-800-423-4588 • 626-574-4000

Technical Helpline: 1-800-265-1943 • 626-574-4577 CE Hotline: 1-800-852-1990 x 4649 • 626-574-4649 Outside these areas, contact your local representative.

In Canada: 1-800-443-1661



Printed on 50% recycled waste paper, including 10% post-consumer waste paper. Forsus and Unitek are trademarks of 3M Unitek.

© 2002 3M Unitek REF 016-922-1 0201