

# Unitek™ Miniature Twin Brackets

## Outstanding Performance – Smaller Size

Unitek™ Miniature Twin Brackets provide all the proven clinical characteristics of full size twin brackets, in a smaller size for greater comfort and improved aesthetics.

Manufactured from sturdy 17-4 stainless steel, metal injection molded Unitek Miniature Twin brackets incorporate a micro-etched base for bond strength equivalent to mesh. Also featured are easy-to-read permanent I.D. marks and tie-wing ink marks. In addition, Unitek Miniature Twin brackets include a smooth, bi-directional hook for superior elastic engagement. Unitek Miniature Twin brackets promise outstanding performance for patients and practitioners alike.



**APC™ II**  
Adhesive Coated Appliance System

**APC™ PLUS**  
Adhesive Coated Appliance System

**For more information, contact your  
3M representative today!**



A micro-etched base features strength comparable to mesh.

Torque-in-base design with compound base contour.



Ample under tie-wing area means smoother and more secure ligation.

# Unitek™ Miniature Twin Metal Brackets – MBT™ Appliance System



MAXILLARY			IN/OUT		Hook	D/R	Part Number – .018 in.		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm			L	R	L	R
Central	+17°	+4°	.031	0.79			020-601	020-602	020-701	020-702
Lateral	+8°	+8°	.044	1.12			020-603	020-604	020-703	020-704
Cuspid	0°	+8°	.026	0.66	D/Hk		020-605	020-606	020-705	020-706
Bicuspid	-7°	0°	.028	0.71	D/Hk		020-607	020-608	020-707	020-708
MANDIBULAR			IN/OUT		Hook	D/R	Part Number – .018 in.		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm			L	R	L	R
Anterior	-6°	0°	.044	1.12			020-609		020-709	
Cuspid	0°	+3°	.018	0.46	D/Hk	6°	020-610	020-611	020-710	020-711
1st Bicuspid	-12°	+2°	.030	0.76	D/Hk		020-612	020-613	020-712	020-713
2nd Bicuspid	-17°	+2°	.029	0.74	D/Hk		020-614	020-615	020-714	020-715

Adding Prefix "3" is APC™ II Adhesive System\*

## One Case Kits

One Case Kit – 5×5	Part Number – .018 in.	Part Number – .022 in.	One Case Kit – 3×3	Part Number – .018 in.	Part Number – .022 in.
U/L 5×5 Hooks (Cuspid/Bicuspid)	020-201	020-202	U/L 3×3 Hooks (Cuspid)	020-203	020-204

Adding Prefix "3" is APC™ II Adhesive System\*

Adding Prefix "3" is APC™ II Adhesive System\*

\*Ask for current availability.

# Unitek™ Miniature Twin Metal Brackets – Roth\* Prescription

MAXILLARY			IN/OUT		M/D	Hook	D/R	Part Number – .018 in.		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm	Width			L	R	L	R
Central	+12°	+5°	.031	0.79	3.5			017-201	017-202	017-301	017-302
Lateral	+8°	+9°	.044	1.12	3.0			017-203	017-204	017-303	017-304
Cuspid	0°	+11°	.026	0.66	3.5			017-205	017-206	017-305	017-306
	0°	+11°	.026	0.66	3.5	D/Hk		017-207	017-208	017-307	017-308
Bicuspid	-7°	0°	.028	0.71	3.0			017-209		017-309	
	-7°	0°	.028	0.71	3.0			7017-209		7017-309	
	-7°	0°	.028	0.71	3.0	D/Hk		017-211	017-212	017-311	017-312
	-7°	0°	.028	0.71	3.0	D/Hk		7017-211	7017-212	7017-311	7017-312
MANDIBULAR			IN/OUT		M/D	Hook	D/R	Part Number – .018 in.		Part Number – .022 in.	
Tooth	Torque	Angulation	in.	mm	Width			L	R	L	R
Anterior	-1°	0°	.044	1.12	2.4			017-213		017-313	
Cuspid	-11°	+5°	.018	0.46	3.0		6°	017-377	017-378	017-477	017-478
	-11°	+5°	.018	0.46	3.0	D/Hk	6°	017-379	017-380	017-479	017-480
	-11°	+5°	.028	.071	3.0			017-215	017-216	017-315	017-316
	-11°	+5°	.028	.071	3.0	D/Hk		017-217	017-218	017-317	017-318
1st Bicuspid	-17°	0°	.030	0.76	3.0			017-219		017-319	
	-17°	0°	.030	0.76	3.0			7017-219		7017-319	
	-17°	0°	.030	0.76	3.0	Hk		017-221	017-222	017-321	017-322
	-17°	0°	.030	0.76	3.0	D/Hk		7017-221	7017-222	7017-321	7017-322
2nd Bicuspid	-22°	0°	.029	0.74	3.0			017-220		017-320	
	-22°	0°	.029	0.74	3.0			7017-220		7017-320	
	-22°	0°	.029	0.74	3.0	D/Hk		017-223	017-224	017-323	017-324
	-22°	0°	.029	0.74	3.0	D/Hk		7017-223	7017-224	7017-323	7017-324

Adding Prefix "3" is APC™ II Adhesive System

Adding Prefix "5" is APC™ PLUS Adhesive System

## Ten Case Kits

Ten Case Kit – 5×5	Part Number – .022 in.
U/L 5×5	017-107
U/L 5×5 (Cuspid/Bicuspid HK)	017-105
U/L 5×5 (Cuspid HK)	017-115

\*3M version of this prescription. No endorsement by the Doctor is implied.

## One Case Kits

One Case Kit – 5×5	Part Number – .018 in.	Part Number – .022 in.
U/L 5×5 (Cuspid/Bicuspid HK)	017-102	017-104
U/L 5×5	–	017-106
U/L 5×5 (Cuspid HK)	017-112	017-114
U/L 5×5	017-1001	017-1061
U/L 5×5 (Cuspid/Bicuspid HK)	017-1021	017-1041
U/L 5×5 (Cuspid HK)	017-1121	–

Adding Prefix "3" is APC™ II Adhesive System

Adding Prefix "5" is APC™ PLUS Adhesive System



**3M Oral Care**  
 2510 Conway Avenue  
 St. Paul, MN 55144-1000 USA  
 Phone 1-800-423-4588  
 Web 3M.com/ortho

**3M Canada**  
 Health Care Division  
 300 Tartan Dr.  
 London, ON N5V 4M9 Canada  
 Phone 1-800-443-1661

3M, APC, MBT and Unitek are trademarks of 3M.  
 Used under license in Canada.  
 © 3M 2016. All rights reserved.  
 70-2021-4959-0 1601