

## Choose the Braces That Get the Looks.

Style and technology get together for Clarity™ SL Self-Ligating Braces. Advanced design. Fast and effective. And the smiles just keep on coming.



Clarity™ SL Self-Ligating Braces



Metal Braces With Ligatures

- Durable, translucent ceramic bracket
- Potential for fewer and faster appointments
- Potential for faster treatment
- Smooth, molded, comfortable edges
- Self-Ligating – no ligatures to attract plaque
- Open design helps you maintain good hygiene
- Engineered to move your teeth quickly

**C L A R I T Y™ | S L**  
SELF-LIGATING BRACES

## Clarity™ SL Self-Ligating Braces: A Very Attractive Package.

Every day, millions of people around the world rely on 3M products such as Scotch® Tape, Post-it® Notes and brands like Nexcare™, Scotch-Brite®, Scotchgard™ and Filtrete™ for exceptional quality and breakthrough technology. With Clarity™ SL self-ligating braces you and your orthodontist are on the leading edge of 3M technology.

Clarity SL braces are unlike any other braces available today. Unique design. Sleek good looks. First-rate comfort and excellent performance. Wear Clarity SL braces and you'll be well on your way to a winning smile.


**C L A R I T Y™ | S L**  
SELF-LIGATING BRACES

- Ask your orthodontist for instructions on the safe use of Clarity SL braces.
- The inner portion of Clarity SL braces contains nickel and chromium. Tell your orthodontist if you are allergic to nickel or chromium.

**3M Unitek**  
**Orthodontic Products**  
2724 South Peck Road  
Monrovia, CA 91016 USA  
www.3MUnitek.com

Please recycle. Printed in U.S.A.  
©3M 2007. All rights reserved.  
014-498 0705

The Smile Says It All.



**C L A R I T Y™ | S L**  
SELF-LIGATING BRACES

**3M** Unitek

**3M** Unitek

**Style.  
Comfort.  
Dazzling  
Results.**

You'll love wearing your Clarity SL braces. They are made from high-tech ceramic with smooth, rounded edges for a comfortable fit.

Clarity SL braces are translucent, so they are less visible than metal braces. And they are the coolest looking braces you'll find anywhere.


Without ligatures to stain or hold food particles, they're easier to keep clean. And, when the time comes to take off your Clarity SL braces, all that will remain is your dazzling new smile.

**C L A R I T Y™ | S L**  
SELF-LIGATING BRACES

Clarity™ SL Self-Ligating Braces are the choice you can feel good about. They use the latest technology, so they work fast and look good.

Self-ligating technology removes the need for ligatures, the elastic bands that were used to hold the wire in your braces. This innovation has major advantages. The Clarity SL braces design can help your teeth move faster, make your braces easier to clean, and eliminate the hassle of stained elastics. You can spend less time in the doctor's chair, and maybe even have fewer visits. And because your teeth can move into position more efficiently, your overall treatment time may be shortened too.

**C L A R I T Y™ | S L**  
SELF-LIGATING BRACES



**3M** Unitek