



**PRODUCT CLINICAL DATA SUMMARY**

**Product No. 1510**

**3M Hi Tack Conformable Double Coated Tape**

**Effective: April 2003**

3M Medical Specialties #1510 3M Hi Tack Conformable Double Coated Tape has been subjected to the following safety evaluations:

**In Vitro Cytotoxicity (Agar Overlay)**

Protocol reference: Guess, W. L. et al; "Agar Diffusion Method for Toxicity Screening of Plastics on Cultured Cell Monolayers" J. Pharm. Sci. 54:1545-1547 (1965).(HRFE 00915)

**Results:** 0.0/0.0

**Acute Primary Skin Irritation in Albino Rabbits**

Protocol reference: Draize: Appraisal of the Safety of Chemicals in Food, Drugs and Cosmetics (1965). Published by the Editorial Committee of the Association of Food and Drug Officials of the United States. (HRFE 00915.)

**Results:** 0.0/8.0

These tests are in accordance with the ISO 10993 Part-1 "Biological Evaluation of Medical Devices", as put forth by the FDA. All laboratory testing was conducted in accordance with the FDA Good Laboratory Practices Regulation of 1978.

It is the responsibility of our customers to determine the final suitability of our products for their application.