

Imprint™ 4 VPS Impression Material – Portfolio overview

Product	Dispensing system	Viscosity low  high	Setting version	Maximum working time (23 °C) 	Maximum intra-oral syringing time (37 °C) 	Constant intra-oral setting time (37 °C) 
TRAY MATERIALS						
Imprint™ 4 Penta™ Putty		Putty	Regular Set	1:30	–	2:30
Imprint™ 4 Penta™ Heavy		Heavy Body	Regular Set	2:00	–	2:00
Imprint™ 4 Penta™ Super Quick Heavy		Heavy Body	Fast Set	1:15	–	1:15
WASH MATERIALS						
Imprint™ 4 Super Quick Ultra-Light			Fast Set	1:15	or 0:35	1:15
Imprint™ 4 Light			Regular Set	2:00	or 1:00	2:00
Imprint™ 4 Super Quick Light			Fast Set	1:15	or 0:35	1:15
Imprint™ 4 Regular			Regular Set	2:00	or 1:00	2:00
Imprint™ 4 Super Quick Regular			Fast Set	1:15	or 0:35	1:15

Imprint™ 4 VPS Impression Material – Recommended material combinations

TRAY MATERIALS	RECOMMENDED WASH MATERIALS
1-step technique	
<p>Imprint™ 4 Penta™ Heavy Hydrophilic heavy body</p> 	<p>Imprint™ 4 Light </p> <p>Imprint™ 4 Regular </p>
<p>Imprint™ 4 Penta™ Super Quick Heavy Fast setting hydrophilic heavy body</p> 	<p>Imprint™ 4 Super Quick Light </p> <p>Imprint™ 4 Super Quick Regular </p>
<p>Imprint™ 4 Penta™ Putty Putty consistency</p> 	<p>Imprint™ 4 Regular </p>
2-step technique	
<p>Imprint™ 4 Penta™ Putty Putty consistency</p> 	<p>Imprint™ 4 Super Quick Ultra-Light </p> <p>Imprint™ 4 Light </p>