

A full line of restorative solutions to fit your unique needs.

Versatility & Esthetics



**Filtek™
Supreme Ultra**
Ultra Universal Restorative

- Exceptional esthetics
- Excellent polish retention
- Exceptional handling
- Outstanding strength for anterior and posterior
- Wear resistance similar to enamel*
- True nanotechnology
- 36 shades, 4 opacities



**Filtek™
Supreme Ultra**
Ultra Flowable Restorative

- Excellent esthetics
- Excellent cavity adaptation
- Excellent polish retention
- Low wear
- True nanotechnology
- 12 shades

Efficiency & Simplicity



**Filtek™
Bulk Fill**
Posterior Restorative

- One-step placement
- Fast and easy to use
- Excellent handling
- Excellent wear resistance
- Stress relief
- True nanotechnology
- 5 shades



**Filtek™
Bulk Fill**
Flowable Restorative

- Low shrinkage stress
- Good marginal adaptation
- Excellent wear resistance
- 4 shades

Value



**Paradigm™
Nano Hybrid**
Universal Restorative

- Good wear resistance
- Easy to polish
- Nano-hybrid technology
- 12 shades

*In vitro data on file.

3M ESPE has you covered for every class of direct restoration.

	Versatility & Esthetics		Efficiency & Simplicity		Value
	Filtek™ Supreme Ultra Ultra Universal Restorative (Incremental up to 2 mm)	Filtek™ Supreme Ultra Ultra Flowable Restorative (Incremental up to 2 mm)	Filtek™ Bulk Fill Posterior Restorative (Bulk up to 5 mm)	Filtek™ Bulk Fill Flowable Restorative (Bulk up to 4 mm)	Paradigm™ Nano Hybrid Universal Restorative (Incremental up to 2 mm)
Class I	✓	✓ Liner/Base	✓ Up to 4 mm	✓ Liner/Base	✓
Class II	✓	✓ Liner/Base	✓ Up to 5 mm	✓ Liner/Base	✓
Class III	✓	✓		✓	✓
Class IV	✓				✓
Class V	✓	✓		✓	✓
Minimally Invasive	✓	✓	✓	✓	✓
Stress Bearing	✓		✓		✓

✓ = Recommended uses

Refer to Instructions for Use for full list of indications.

3M ESPE

★★★★★ Review these products at www.3MESPE.com/review

Learn more at www.3MESPE.com/Filtek

Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

Please recycle. Printed in U.S.A.
70-2013-0646-4

3M, ESPE, Filtek and Paradigm are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2014. All rights reserved.