

Table of Contents

Section 2 – Table of Contents

Clarity™ Metal-Reinforced Ceramic Brackets	2.2
Transcend™ Series 6000 Ceramic Brackets	2.18
Transcend™ Series 6000 Debonding Instrument	2.26

non-metal brackets **2**

2

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

2

There is a reason why the Clarity brand of ceramic brackets is associated with high-performance aesthetic treatment. Let's look at the evidence:

- Over three million cases
- Worldwide acceptance
- 10 years of clinical experience
- Numerous reviews in major orthodontic journals worldwide



Clarity™ Ceramic Brackets are the #1 aesthetic bracket choice worldwide. Why? Clarity Brackets feature innovations from 3M that make an impact on practice efficiency, ease of use, patient satisfaction and the standard of care. The chart shows just a sampling of the technical features and benefits of Clarity Brackets.



Clarity™ Brackets

Experience metal-like sliding mechanics. Clarity Brackets feature an innovative metal slot liner.

Expect reliable squeeze debonding thanks to the Clarity Bracket's innovative microcrystalline bracket bonding surface and mesial-distal stress concentrator.

Choose both strength and beauty. Clarity features a strong and translucent monocrySTALLINE bracket material.

Benefit from efficiency on bonding day. 3M Unitek's family of ceramic brackets are the only ceramic brackets to offer the APC™ Adhesive pre-coat option.

Treat to your standards and preference. Clarity Brackets are available in a range of prescriptions, slot sizes, kit packages and hook configurations.


Start treatment with precise bracket placement. The metal slot liner of Clarity Brackets provide horizontal visual guidance for bracket positioning.

the #1 aesthetic bracket choice worldwide



Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets



Also available with your choice of APC™ II Adhesive or with the APC™ PLUS System with color changing adhesive for reliable and efficient bonding.



APC™ II

Available in APC™ II Adhesive Coated Appliance System

APC™ PLUS







Available in APC™ PLUS Adhesive Coated Appliance System

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection By Tooth – Maxillary .018 in.

2

	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Central 	+22°	5°	1.12	3.5		UL	6700-651	6800-651	6400-651
	+22°	5°	1.12	3.5		UR	6700-652	6800-652	6400-652
	+17°	4°	1.12	3.5		UL	6700-801	6800-801	6400-801
	+17°	4°	1.12	3.5		UR	6700-802	6800-802	6400-802
	+12°	5°	1.12	3.5		UL	6700-601	6800-601	6400-601
	+12°	5°	1.12	3.5		UR	6700-602	6800-602	6400-602
	0°	0°	.99	3.5		UL/UR	6700-820	6800-820	6400-820
Lateral 	+14°	8°	1.47	3.0		UL	6700-653	6800-653	6400-653
	+14°	8°	1.47	3.0		UR	6700-654	6800-654	6400-654
	+10°	8°	1.47	3.0		UL	6700-803	6800-803	6400-803
	+10°	8°	1.47	3.0		UR	6700-804	6800-804	6400-804
	+8°	9°	1.47	3.0		UL	6700-603	6800-603	6400-603
	+8°	9°	1.47	3.0		UR	6700-604	6800-604	6400-604
	0°	0°	.99	3.0		UL/UR	6700-821	6800-821	6400-821
Cuspid  	-7°	8°	1.04	3.5	D	UL	6700-683	6800-683	6400-683
	-7°	8°	1.04	3.5	D	UR	6700-682	6800-682	6400-682
	-7°	8°	1.04	3.5		UL	6700-805	6800-805	6400-805
	-7°	8°	1.04	3.5		UR	6700-806	6800-806	6400-806
	+7°	5°	1.04	3.5		UL	6700-671	6800-671	6400-671
	+7°	5°	1.04	3.5		UR	6700-670	6800-670	6400-670
	+7°	5°	1.04	3.5	D	UL	6700-655	6800-655	6400-655
	+7°	5°	1.04	3.5	D	UR	6700-656	6800-656	6400-656
	0°	11°	1.04	3.5		UL	6700-607	6800-607	6400-607
	0°	11°	1.04	3.5		UR	6700-608	6800-608	6400-608
	0°	11°	1.04	3.5	D	UL	6700-605	6800-605	6400-605
	0°	11°	1.04	3.5	D	UR	6700-606	6800-606	6400-606
	0°	8°	1.04	3.5	D	UL	6700-807	6800-807	6400-807
	0°	8°	1.04	3.5	D	UR	6700-808	6800-808	6400-808
	0°	8°	1.04	3.5		UL	6700-817	6800-817	6400-817
0°	8°	1.04	3.5		UR	6700-818	6800-818	6400-818	
0°	0°	.99	3.5		UNIV	6700-822	6800-822	6400-822	
Bicuspid  	-7°	0°	1.04	3.0		UL/UR	6700-611	6800-611	6400-611
	-7°	0°	1.04	3.5	D	UL	6700-609	6800-609	6400-609
	-7°	0°	1.04	3.5	D	UR	6700-610	6800-610	6400-610
	0°	0°	1.04	3.5	D	UL	6700-675	6800-675	6400-675
	0°	0°	1.04	3.5	D	UR	6700-674	6800-674	6400-674
	0°	0°	.99	3.0		UNIV	6700-823	6800-823	6400-823

6700


Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800

Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

.018 Kits – Page 2.16

Note






APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
See page 2.7 and 2.13 for Clarity™ Roth™ Prescription Brackets.	MBT™ Brackets Rx pictured in most cases. 

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection By Tooth – Mandibular .018 in.



	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Anterior 	-6°	0°	1.52	2.5		LL/LR	6700-810	6800-810	6400-810
	-1°	0°	1.32	2.5		LL/LR	6700-621	6800-621	6400-621
	0°	0°	.99	2.5		LL/LR	6700-824	6800-824	6400-824
Cuspid 	-11°	7°	1.04	3.5		LL	6700-627	6800-627	6400-627
	-11°	7°	1.04	3.5		LR	6700-628	6800-628	6400-628
	-11°	7°	1.04	3.5	D	LL	6700-625	6800-625	6400-625
	-11°	7°	1.04	3.5	D	LR	6700-626	6800-626	6400-626
	+7°	5°	1.04	3.5		LL	6700-673	6800-673	6400-673
	+7°	5°	1.04	3.5		LR	6700-672	6800-672	6400-672
	+7°	5°	1.04	3.5	D	LL	6700-657	6800-657	6400-657
	+7°	5°	1.04	3.5	D	LR	6700-658	6800-658	6400-658
	-6°	3°	.79	3.5		LL	6700-695	6800-695	6400-695
	-6°	3°	.79	3.5		LR	6700-696	6800-696	6400-696
	-6°	3°	.79	3.5	D	LL	6700-685	6800-685	6400-685
	-6°	3°	.79	3.5	D	LR	6700-684	6800-684	6400-684
	0°	3°	.79	3.5		LL	6700-697	6800-697	6400-697
	0°	3°	.79	3.5		LR	6700-698	6800-698	6400-698
	0°	3°	.79	3.5	D	LL	6700-699	6800-699	6400-699
	0°	3°	.79	3.5	D	LR	6700-700	6800-700	6400-700
	0°	0°	.99	3.5		UNIV	6700-822	6800-822	6400-822
First Bicusp 	-17°	0°	1.04	3.0		LL/LR	6700-636	6800-636	6400-636
	-17°	0°	1.04	3.5	D	LL	6700-629	6800-629	6400-629
	-17°	0°	1.04	3.5	D	LR	6700-630	6800-630	6400-630
	-12°	2°	1.04	3.0		LL	6700-687	6800-687	6400-687
	-12°	2°	1.04	3.0		LR	6700-686	6800-686	6400-686
	-12°	2°	1.04	3.5	D	LL	6700-689	6800-689	6400-689
	-12°	2°	1.04	3.5	D	LR	6700-688	6800-688	6400-688
	-7°	0°	1.04	3.0		LL/LR	6700-660	6800-660	6400-660
	-7°	0°	1.04	3.5	D	LL	6700-661	6800-661	6400-661
-7°	0°	1.04	3.5	D	LR	6700-662	6800-662	6400-662	
Bicuspid 	0°	0°	.99	3.0		UNIV	6700-823	6800-823	6400-823
Second Bicusp 	-22°	0°	1.04	3.0		LL/LR	6700-644	6800-644	6400-644
	-22°	0°	1.04	3.5	D	LL	6700-637	6800-637	6400-637
	-22°	0°	1.04	3.5	D	LR	6700-638	6800-638	6400-638
	-17°	2°	1.04	3.0		LL	6700-691	6800-691	6400-691
	-17°	2°	1.04	3.0		LR	6700-690	6800-690	6400-690
	-17°	2°	1.04	3.5	D	LL	6700-693	6800-693	6400-693
	-17°	2°	1.04	3.5	D	LR	6700-692	6800-692	6400-692
	-11°	0°	1.04	3.0		LL/LR	6700-659	6800-659	6400-659
	-11°	0°	1.04	3.5	D	LL	6700-677	6800-677	6400-677
	-11°	0°	1.04	3.5	D	LR	6700-676	6800-676	6400-676

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
See page 2.6 and 2.12 for the Clarity™ MBT™ Bracket Prescription.	MBT™ Brackets Rx pictured in most cases. MBT™

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.









6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection by MBT™ Appliance System Prescription .018 in.

2

	Torque	Ang.	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Maxillary								
Central 	+17°	4°	3.5		UL	6700-801	6800-801	6400-801
	+17°	4°	3.5		UR	6700-802	6800-802	6400-802
Lateral 	+10°	8°	3.0		UL	6700-803	6800-803	6400-803
	+10°	8°	3.0		UR	6700-804	6800-804	6400-804
Cuspid 	-7°	8°	3.5	D	UL	6700-683	6800-683	6400-683
	-7°	8°	3.5	D	UR	6700-682	6800-682	6400-682
	-7°	8°	3.5		UL	6700-805	6800-805	6400-805
	-7°	8°	3.5		UR	6700-806	6800-806	6400-806
	0°	8°	3.5	D	UL	6700-807	6800-807	6400-807
	0°	8°	3.5	D	UR	6700-808	6800-808	6400-808
	0°	8°	3.5		UL	6700-817	6800-817	6400-817
	0°	8°	3.5		UR	6700-818	6800-818	6400-818
Bicuspid 	-7°	0°	3.0		UL/UR	6700-611	6800-611	6400-611
	-7°	0°	3.5	D	UL	6700-609	6800-609	6400-609
	-7°	0°	3.5	D	UR	6700-610	6800-610	6400-610
Mandibular								
Anterior 	-6°	0°	2.5		LL/LR	6700-810	6800-810	6400-810
Cuspid 	-6°	3°	3.5		LL	6700-695	6800-695	6400-695
	-6°	3°	3.5		LR	6700-696	6800-696	6400-696
	-6°	3°	3.5	D	LL	6700-685	6800-685	6400-685
	-6°	3°	3.5	D	LR	6700-684	6800-684	6400-684
	0°	3°	3.5		LL	6700-697	6800-697	6400-697
	0°	3°	3.5		LR	6700-698	6800-698	6400-698
	0°	3°	3.5	D	LL	6700-699	6800-699	6400-699
	0°	3°	3.5	D	LR	6700-700	6800-700	6400-700
First Bicuspid 	-12°	2°	3.0		LL	6700-687	6800-687	6400-687
	-12°	2°	3.0		LR	6700-686	6800-686	6400-686
	-12°	2°	3.5	D	LL	6700-689	6800-689	6400-689
	-12°	2°	3.5	D	LR	6700-688	6800-688	6400-688
Second Bicuspid 	-17°	2°	3.0		LL	6700-691	6800-691	6400-691
	-17°	2°	3.0		LR	6700-690	6800-690	6400-690
	-17°	2°	3.5	D	LL	6700-693	6800-693	6400-693
	-17°	2°	3.5	D	LR	6700-692	6800-692	6400-692


6700

Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800

Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.









.018 Kits – Page 2.16

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
Measurements are in millimeters unless otherwise noted.	
For Clarity™ MBT™ Brackets Rx One Case Kits see page 2.16 and 2.17.	For detailed buccal tube information see the Buccal Tubes section. 

Non-Metal Brackets



Clarity™ Metal-Reinforced Ceramic Brackets

Selection by Roth* Prescription .018 in.

Maxillary	Torque	Ang.	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Central 	+12°	5°	3.5		UL	6700-601	6800-601	6400-601
	+12°	5°	3.5		UR	6700-602	6800-602	6400-602
Lateral 	+8°	9°	3.0		UL	6700-603	6800-603	6400-603
	+8°	9°	3.0		UR	6700-604	6800-604	6400-604
Cuspid 	0°	11°	3.5		UL	6700-607	6800-607	6400-607
	0°	11°	3.5		UR	6700-608	6800-608	6400-608
	0°	11°	3.5	D	UL	6700-605	6800-605	6400-605
	0°	11°	3.5	D	UR	6700-606	6800-606	6400-606
Bicuspid 	-7°	0°	3.0		UL/UR	6700-611	6800-611	6400-611
	-7°	0°	3.5	D	UL	6700-609	6800-609	6400-609
	-7°	0°	3.5	D	UR	6700-610	6800-610	6400-610
Mandibular								
Anterior 	-1°	0°	2.5		LL/LR	6700-621	6800-621	6400-621
Cuspid 	-11°	7°	3.5		LL	6700-627	6800-627	6400-627
	-11°	7°	3.5		LR	6700-628	6800-628	6400-628
	-11°	7°	3.5	D	LL	6700-625	6800-625	6400-625
	-11°	7°	3.5	D	LR	6700-626	6800-626	6400-626
First Bicuspid 	-17°	0°	3.0		LL/LR	6700-636	6800-636	6400-636
	-17°	0°	3.5	D	LL	6700-629	6800-629	6400-629
	-17°	0°	3.5	D	LR	6700-630	6800-629	6400-630
Second Bicuspid 	-22°	0°	3.0		LL/LR	6700-644	6800-644	6400-644
	-22°	0°	3.5	D	LL	6700-637	6800-637	6400-637
	-22°	0°	3.5	D	LR	6700-638	6800-638	6400-638

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.









Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
For Clarity™ Roth* One Case Kits please see page 2.16 and 2.17. 	For detailed buccal tube information see the Buccal Tubes section. 

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection by High Torque Prescription .018 in.



2

Maxillary	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Central 	+22°	5°	1.12	3.5		UL	6700-651	6800-651	6400-651
	+22°	5°	1.12	3.5		UR	6700-652	6800-652	6400-652
Lateral 	+14°	8°	1.47	3.0		UL	6700-653	6800-653	6400-653
	+14°	8°	1.47	3.0		UR	6700-654	6800-654	6400-654
Cuspid 	+7°	5°	1.04	3.5		UL	6700-671	6800-671	6400-671
	+7°	5°	1.04	3.5		UR	6700-670	6800-670	6400-670
Cuspid w/Hook	+7°	5°	1.04	3.5	D	UL	6700-655	6800-655	6400-655
	+7°	5°	1.04	3.5	D	UR	6700-656	6800-656	6400-656
Bicuspid 	0°	0°	1.04	3.0		UL/UR	6700-823	6800-823	6400-823
Bicuspid w/Hook	0°	0°	1.04	3.5	D	UL	6700-675	6800-675	6400-675
	0°	0°	1.04	3.5	D	UR	6700-674	6800-674	6400-674
Mandibular									
Anterior 	0°	0°	0.99	2.5		LL/LR	6700-824	6800-824	6400-824
Cuspid 	+7°	5°	1.04	3.5		LL	6700-673	6800-673	6400-673
	+7°	5°	1.04	3.5		LR	6700-672	6800-672	6400-672
Cuspid w/Hook	+7°	5°	1.04	3.5	D	LL	6700-657	6800-657	6400-657
	+7°	5°	1.04	3.5	D	LR	6700-658	6800-658	6400-658
1st Bicuspid 	-7°	0°	1.04	3.0		LL/LR	6700-660	6800-660	6400-660
2nd Bicuspid 	-11°	0°	1.04	3.0		LL/LR	6700-659	6800-659	6400-659
1st Bicuspid w/Hook	-7°	0°	1.04	3.5	D	LL	6700-661	6800-661	6400-661
	-7°	0°	1.04	3.5	D	LR	6700-662	6800-662	6400-662
2nd Bicuspid w/Hook	-11°	0°	1.04	3.5	D	LL	6700-677	6800-677	6400-677
	-11°	0°	1.04	3.5	D	LR	6700-676	6800-676	6400-676

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.






.018 Kits – Page 2.16

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted. 1 2 3 4 Buccal Tube is available with a bondable base.
For Clarity™ High Torque One Case Kits please see page 2.16 and 2.17. 	For detailed buccal tube information see the Buccal Tubes section. 


Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection by Standard Edgewise Prescription .018 in.

	Torque	Ang.	M/D Width	Hook	L/R	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Central 	0°	0°	3.5		UL/UR	6700-820	6800-820	6400-820
Lateral 	0°	0°	3.0		UL/UR	6700-821	6800-821	6400-821
Cuspid 	0°	0°	3.5		UNIV	6700-822	6800-822	6400-822
Bicuspid 	0°	0°	3.0		UNIV	6700-823	6800-823	6400-823
Anterior 	0°	0°	2.5		LL/LR	6700-824 6700 Use this prefix to order APC™ II Adhesive Coated Appliance System.	6800-824 6800 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.	6400-824









Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
For detailed buccal tube information see the Buccal Tubes section. 	

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection By Tooth – Maxillary .022 in.

2

	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Central 	+22°	5°	1.12	3.5		UL	6700-751	6800-751	6400-751
	+22°	5°	1.12	3.5		UR	6700-752	6800-752	6400-752
	+17°	4°	1.12	3.5		UL	6700-901	6800-901	6400-901
	+17°	4°	1.12	3.5		UR	6700-902	6800-902	6400-902
	+12°	5°	1.12	3.5		UL	6700-701	6800-701	6400-701
	+12°	5°	1.12	3.5		UR	6700-702	6800-702	6400-702
	0°	0°	.99	3.5		UL/UR	6700-920	6800-920	6400-920
Lateral 	+14°	8°	1.47	3.0		UL	6700-753	6800-753	6400-753
	+14°	8°	1.47	3.0		UR	6700-754	6800-754	6400-754
	+10°	8°	1.47	3.0		UL	6700-903	6800-903	6400-903
	+10°	8°	1.47	3.0		UR	6700-904	6800-904	6400-904
	+8°	9°	1.47	3.0		UL	6700-703	6800-703	6400-703
	+8°	9°	1.47	3.0		UR	6700-704	6800-704	6400-704
	0°	0°	.99	3.0		UL/UR	6700-921	6800-921	6400-921
Cuspid  	-7°	8°	1.04	3.5	D	UL	6700-783	6800-783	6400-783
	-7°	8°	1.04	3.5	D	UR	6700-782	6800-782	6400-782
	-7°	8°	1.04	3.5		UL	6700-905	6800-905	6400-905
	-7°	8°	1.04	3.5		UR	6700-906	6800-906	6400-906
	+7°	5°	1.04	3.5		UL	6700-771	6800-771	6400-771
	+7°	5°	1.04	3.5		UR	6700-770	6800-770	6400-770
	+7°	5°	1.04	3.5	D	UL	6700-755	6800-755	6400-755
	+7°	5°	1.04	3.5	D	UR	6700-756	6800-756	6400-756
	0°	11°	1.04	3.5		UL	6700-707	6800-707	6400-707
	0°	11°	1.04	3.5		UR	6700-708	6800-708	6400-708
	0°	11°	1.04	3.5	D	UL	6700-705	6800-705	6400-705
	0°	11°	1.04	3.5	D	UR	6700-706	6800-706	6400-706
	0°	8°	1.04	3.5	D	UL	6700-907	6800-907	6400-907
	0°	8°	1.04	3.5	D	UR	6700-908	6800-908	6400-908
	0°	8°	1.04	3.5		UL	6700-917	6800-917	6400-917
0°	8°	1.04	3.5		UR	6700-918	6800-918	6400-918	
0°	0°	.99	3.5		UNIV	6700-922	6800-922	6400-922	
Bicuspid  	-7°	0°	1.04	3.0		UL/UR	6700-711	6800-711	6400-711
	-7°	0°	1.04	3.5	D	UL	6700-709	6800-709	6400-709
	-7°	0°	1.04	3.5	D	UR	6700-710	6800-710	6400-710
	0°	0°	1.04	3.5	D	UL	6700-775	6800-775	6400-775
	0°	0°	1.04	3.5	D	UR	6700-774	6800-774	6400-774
	0°	0°	.99	3.0		UNIV	6700-923	6800-923	6400-923

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.






.022 Kits – Page 2.17

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
See page 2.7 and 2.13 for Clarity™ Roth® Prescription Brackets.	MBT™ Brackets Rx pictured in most cases. MBT™

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection By Tooth – Mandibular .022 in.

	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Anterior 	-6°	0°	1.52	2.5		LL/LR	6700-910	6800-910	6400-910
	-1°	0°	1.32	2.5		LL/LR	6700-721	6800-721	6400-721
	0°	0°	.99	2.5		LL/LR	6700-924	6800-924	6400-924
Cuspid 	-11°	7°	1.04	3.5		LL	6700-727	6800-727	6400-727
	-11°	7°	1.04	3.5		LR	6700-728	6800-728	6400-728
	-11°	7°	1.04	3.5	D	LL	6700-725	6800-725	6400-725
	-11°	7°	1.04	3.5	D	LR	6700-726	6800-726	6400-726
	+7°	5°	1.04	3.5		LL	6700-773	6800-773	6400-773
	+7°	5°	1.04	3.5		LR	6700-772	6800-772	6400-772
	+7°	5°	1.04	3.5	D	LL	6700-757	6800-757	6400-757
	+7°	5°	1.04	3.5	D	LR	6700-758	6800-758	6400-758
	-6°	3°	.79	3.5		LL	6700-795	6800-795	6400-795
	-6°	3°	.79	3.5		LR	6700-796	6800-796	6400-796
	-6°	3°	.79	3.5	D	LL	6700-785	6800-785	6400-785
	-6°	3°	.79	3.5	D	LR	6700-784	6800-784	6400-784
	0°	3°	.79	3.5		LL	6700-797	6800-797	6400-797
	0°	3°	.79	3.5		LR	6700-798	6800-798	6400-798
	0°	3°	.79	3.5	D	LL	6700-799	6800-799	6400-799
	0°	3°	.79	3.5	D	LR	6700-800	6800-800	6400-800
0°	0°	.99	3.5		UNIV	6700-922	6800-922	6400-922	
First Bicuspid 	-17°	0°	1.04	3.0		LL/LR	6700-736	6800-736	6400-736
	-17°	0°	1.04	3.5	D	LL	6700-729	6800-729	6400-729
	-17°	0°	1.04	3.5	D	LR	6700-730	6800-730	6400-730
	-12°	2°	1.04	3.0		LL	6700-787	6800-787	6400-787
	-12°	2°	1.04	3.0		LR	6700-786	6800-786	6400-786
	-12°	2°	1.04	3.5	D	LL	6700-789	6800-789	6400-789
	-12°	2°	1.04	3.5	D	LR	6700-788	6800-788	6400-788
	-7°	0°	1.04	3.0		LL/LR	6700-760	6800-760	6400-760
	-7°	0°	1.04	3.5	D	LL	6700-761	6800-761	6400-761
	-7°	0°	1.04	3.5	D	LR	6700-762	6800-762	6400-762
Bicuspid 	0°	0°	.99	3.0		UNIV	6700-923	6800-923	6400-923
Second Bicuspid 	-22°	0°	1.04	3.0		LL/LR	6700-744	6800-744	6400-744
	-22°	0°	1.04	3.5	D	LL	6700-737	6800-737	6400-737
	-22°	0°	1.04	3.5	D	LR	6700-738	6800-738	6400-738
	-17°	2°	1.04	3.0		LL	6700-791	6800-791	6400-791
	-17°	2°	1.04	3.0		LR	6700-790	6800-790	6400-790
	-17°	2°	1.04	3.5	D	LL	6700-793	6800-793	6400-793
	-17°	2°	1.04	3.5	D	LR	6700-792	6800-792	6400-792
	-11°	0°	1.04	3.0		LL/LR	6700-759	6800-759	6400-759
	-11°	0°	1.04	3.5	D	LL	6700-777	6800-777	6400-777
	-11°	0°	1.04	3.5	D	LR	6700-776	6800-776	6400-776

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
See page 2.6 and 2.12 for the Clarity™ MBT™ Brackets Prescription.	MBT™ Brackets Rx pictured in most cases. MBT™

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.












Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection by MBT™ Appliance System Prescription .022 in.


2

Maxillary	Torque	Ang.	M/D Width	Hook	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Central 	+17°	4°	3.5		UL	6700-901	6800-901	6400-901
	+17°	4°	3.5		UR	6700-902	6800-902	6400-902
Lateral 	+10°	8°	3.0		UL	6700-903	6800-903	6400-903
	+10°	8°	3.0		UR	6700-904	6800-904	6400-904
Cuspid 	-7°	8°	3.5	D	UL	6700-783	6800-783	6400-783
	-7°	8°	3.5	D	UR	6700-782	6800-782	6400-782
	-7°	8°	3.5		UL	6700-905	6800-782	6400-905
	-7°	8°	3.5		UR	6700-906	6800-906	6400-906
	0°	8°	3.5	D	UL	6700-907	6800-907	6400-907
	0°	8°	3.5	D	UR	6700-908	6800-908	6400-908
	0°	8°	3.5		UL	6700-917	6800-917	6400-917
	0°	8°	3.5		UR	6700-918	6800-918	6400-918
Bicuspid 	-7°	0°	3.0		UL/UR	6700-711	6800-711	6400-711
	-7°	0°	3.5	D	UL	6700-709	6800-709	6400-709
	-7°	0°	3.5	D	UR	6700-710	6800-710	6400-710
Mandibular								
Anterior 	-6°	0°	2.5		LL/LR	6700-910	6800-910	6400-910
Cuspid 	-6°	3°	3.5		LL	6700-795	6800-795	6400-795
	-6°	3°	3.5		LR	6700-796	6800-796	6400-796
	-6°	3°	3.5	D	LL	6700-785	6800-785	6400-785
	-6°	3°	3.5	D	LR	6700-784	6800-784	6400-784
	0°	3°	3.5		LL	6700-797	6800-797	6400-797
	0°	3°	3.5		LR	6700-798	6800-798	6400-798
	0°	3°	3.5	D	LL	6700-799	6800-799	6400-799
	0°	3°	3.5	D	LR	6700-800	6800-800	6400-800
First Bicuspid 	-12°	2°	3.0		LL	6700-787	6800-787	6400-787
Second Bicuspid 	-12°	2°	3.0		LR	6700-786	6800-786	6400-786
	-12°	2°	3.5	D	LL	6700-789	6800-789	6400-789
	-12°	2°	3.5	D	LR	6700-788	6800-788	6400-788
Second Bicuspid 	-17°	2°	3.0		LL	6700-791	6800-791	6400-791
	-17°	2°	3.0		LR	6700-790	6800-790	6400-790
	-17°	2°	3.5	D	LL	6700-793	6800-793	6400-793
	-17°	2°	3.5	D	LR	6700-792	6800-792	6400-792

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.










.022 Kits – Page 2.17

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
Measurements are in millimeters unless otherwise noted.	
For Clarity™ MBT™ Brackets Rx One Case Kits see page 2.16 and 2.17.	For detailed buccal tube information see the Buccal Tubes section. 

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets


Selection by Roth* Prescription .022 in.

Maxillary	Torque	Ang.	M/D Width	Hook	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Central 	+12°	5°	3.5		UL	6700-701	6800-701	6400-701
	+12°	5°	3.5		UR	6700-702	6800-702	6400-702
Lateral 	+8°	9°	3.0		UL	6700-703	6800-703	6400-703
	+8°	9°	3.0		UR	6700-704	6800-704	6400-704
Cuspid 	0°	11°	3.5		UL	6700-707	6800-707	6400-707
	0°	11°	3.5		UR	6700-708	6800-708	6400-708
	0°	11°	3.5	D	UL	6700-705	6800-705	6400-705
	0°	11°	3.5	D	UR	6700-706	6800-706	6400-706
Bicuspid 	-7°	0°	3.0		UL/UR	6700-711	6800-711	6400-711
	-7°	0°	3.5	D	UL	6700-709	6800-709	6400-709
	-7°	0°	3.5	D	UR	6700-710	6800-710	6400-710
Mandibular								
Anterior 	-1°	0°	2.5		LL/LR	6700-721	6800-721	6400-721
Cuspid 	-11°	7°	3.5		LL	6700-727	6800-727	6400-727
	-11°	7°	3.5		LR	6700-728	6800-728	6400-728
	-11°	7°	3.5	D	LL	6700-725	6800-725	6400-725
	-11°	7°	3.5	D	LR	6700-726	6800-726	6400-726
First Bicuspid 	-17°	0°	3.0		LL/LR	6700-736	6800-736	6400-736
Second Bicuspid 	-17°	0°	3.5	D	LL	6700-729	6800-729	6400-729
	-17°	0°	3.5	D	LR	6700-730	6800-730	6400-730
Second Bicuspid 	-22°	0°	3.0		LL/LR	6700-744	6800-744	6400-744
	-22°	0°	3.5	D	LL	6700-737	6800-737	6400-737
	-22°	0°	3.5	D	LR	6700-738	6800-738	6400-738

6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

.022 Kits – Page 2.17














Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
For Clarity™ Roth* One Case Kits please see page 2.16 and 2.17.	For detailed buccal tube information see the Buccal Tubes section. 

Non-Metal Brackets



Clarity™ Metal-Reinforced Ceramic Brackets

Selection by High Torque Prescription .022 in.

2

Maxillary	Torque	Ang.	In/Out	M/D Width	Hook	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
 Central	+22°	5°	1.12	3.5		UL	6700-751	6800-751	6400-751
	+22°	5°	1.12	3.5		UR	6700-752	6800-752	6400-752
 Lateral	+14°	8°	1.47	3.0		UL	6700-753	6800-753	6400-753
	+14°	8°	1.47	3.0		UR	6700-754	6800-754	6400-754
 Cuspid	+7°	5°	1.04	3.5		UL	6700-771	6800-771	6400-771
	+7°	5°	1.04	3.5		UR	6700-770	6800-770	6400-770
 Cuspid w/Hook	+7°	5°	1.04	3.5	D	UL	6700-755	6800-755	6400-755
	+7°	5°	1.04	3.5	D	UR	6700-756	6800-756	6400-756
 Bicuspid	0°	0°	1.04	3.0		UL/UR	6700-923	6800-923	6400-923
 Bicuspid w/Hook	0°	0°	1.04	3.5	D	UL	6700-775	6800-775	6400-775
	0°	0°	1.04	3.5	D	UR	6700-774	6800-774	6400-774
Mandibular									
 Anterior	0°	0°	.99	2.5		LL/LR	6700-924	6800-924	6400-924
 Cuspid	+7°	5°	1.04	3.5		LL	6700-773	6800-773	6400-773
	+7°	5°	1.04	3.5		LR	6700-772	6800-772	6400-772
 Cuspid w/Hook	+7°	5°	1.04	3.5	D	LL	6700-757	6800-757	6400-757
	+7°	5°	1.04	3.5	D	LR	6700-758	6800-758	6400-758
 1st Bicuspid	-7°	0°	1.04	3.0		LL/LR	6700-760	6800-760	6400-760
 2nd Bicuspid	-11°	0°	1.04	3.0		LL/LR	6700-759	6800-759	6400-759
 1st Bicuspid w/Hook	-7°	0°	1.04	3.5	D	LL	6700-761	6800-761	6400-761
	-7°	0°	1.04	3.5	D	LR	6700-762	6800-762	6400-762
 2nd Bicuspid w/Hook	-11°	0°	1.04	3.5	D	LL	6700-777	6800-777	6400-777
	-11°	0°	1.04	3.5	D	LR	6700-776	6800-776	6400-776

.022 Kits – Page 2.17

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted. 1 2 3 4 Buccal Tube is available with a bondable base.
For Clarity™ High Torque One Case Kits please see page 2.16 and 2.17. 	For detailed buccal tube information see the Buccal Tubes section. 






6700
Use this prefix to order APC™ II Adhesive Coated Appliance System.


6800
Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

Selection by Standard Edgewise Prescription .022 in.

	Torque	Ang.	M/D Width	L/R	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Central 	0°	0°	3.5	UL/UR	6700-920	6800-920	6400-920
Lateral 	0°	0°	3.0	UL/UR	6700-921	6800-921	6400-921
Cuspid 	0°	0°	3.5	UNIV	6700-922	6800-922	6400-922
Bicuspid 	0°	0°	3.0	UNIV	6700-923	6800-923	6400-923
Anterior 	0°	0°	2.5	LL/LR	6700 Use this prefix to order APC™ II Adhesive Coated Appliance System.	6800 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.	6400-924

Note	
APC™ II Clarity™ Brackets sold in units of five.	APC™ PLUS Clarity™ Brackets sold in units of five.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
For detailed buccal tube information see the Buccal Tubes section. 	

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

One Case Kits .018 in.

2

Clarity MBT™ Brackets One Case Kits

	Type	Hook	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Upper	3 x 3		6700-208	6800-208	6400-208
Upper	3 x 3	Cuspid	6700-212	6800-212	6400-212
Upper	5 x 5		6700-210	6800-210	6400-210
Upper	5 x 5	Cuspid	6700-214	6800-214	6400-214
U/L	3 x 3		6700-209	6800-209	6400-209
U/L	3 x 3	Cuspid	6700-213	6800-213	6400-213
U/L	5 x 5		6700-211	6800-211	6400-211
U/L	5 x 5	Cuspid	6700-215	6800-215	6400-215

Clarity Roth* One Case Kits

	Type	Hook	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
Upper	3 x 3		6700-102	6800-102	6400-102
Upper	3 x 3	Cuspid	6700-103	6800-103	6400-103
Upper	5 x 5		6700-104	6800-104	6400-104
Upper	5 x 5	Cuspid & Bicuspid	6700-105	6800-105	6400-105
U/L	3 x 3		6700-106	6800-106	6400-106
U/L	3 x 3	Cuspid	6700-107	6800-107	6400-107
U/L	5 x 5		6700-100	6800-100	6400-100
U/L	5 x 5	Cuspid & Bicuspid	6700-101	6800-101	6400-101
U/L	5 x 5	Cuspid	6700-132	6800-132	6400-132
Upper	5 x 5	Cuspid	6700-134	6800-134	6400-134

Clarity High Torque One Case Kits

	Type	Hook	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
U	5 x 5	Cuspid	6700-400	6800-400	6400-400
U/L	5 x 5	Cuspid	6700-401	6800-401	6400-401
U/L	3 x 3	Cuspid	6700-402	6800-402	6400-402
U	3 x 3	Cuspid	6700-403	6800-403	6400-403
U	5 x 5		6700-404	6800-404	6400-404
U/L	5 x 5		6700-405	6800-405	6400-405
U/L	3 x 3		6700-406	6800-406	6400-406
U	3 x 3		6700-407	6800-407	6400-407

Clarity Standard Edgewise One Case Kits

	Type	Hook	APC™ II .018 in.	APC™ PLUS .018 in.	Non-coated .018 in.
U/L	5 x 5		6700-300	6800-300	6400-300

6700

Use this prefix to order APC™ II Adhesive Coated Appliance System.

6800

Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Note

* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.

Non-Metal Brackets

Clarity™ Metal-Reinforced Ceramic Brackets

One Case Kits .022 in.



Clarity MBT™ Brackets One Case Kits					
	Type	Hook	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Upper	3 x 3		6700-200	6800-200	6400-200
Upper	3 x 3	Cuspid	6700-204	6800-204	6400-204
Upper	5 x 5		6700-202	6800-202	6400-202
Upper	5 x 5	Cuspid	6700-206	6800-206	6400-206
U/L	3 x 3		6700-201	6800-201	6400-201
U/L	3 x 3	Cuspid	6700-205	6800-205	6400-205
U/L	5 x 5		6700-203	6800-203	6400-203
U/L	5 x 5	Cuspid	6700-207	6800-207	6400-207

Clarity Roth* One Case Kits					
	Type	Hook	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
Upper	3 x 3		6700-122	6800-122	6400-122
Upper	3 x 3	Cuspid	6700-123	6800-123	6400-123
Upper	5 x 5		6700-124	6800-124	6400-124
Upper	5 x 5	Cuspid & Bicuspид	6700-125	6800-125	6400-125
U/L	3 x 3		6700-126	6800-126	6400-126
U/L	3 x 3	Cuspid & Bicuspид	6700-127	6800-127	6400-127
U/L	5 x 5		6700-120	6800-120	6400-120
U/L	5 x 5	Cuspid & Bicuspид	6700-121	6800-121	6400-121
U/L	5 x 5	Cuspid	6700-133	6800-133	6400-133
Upper	5 x 5	Cuspid	6700-135	6800-135	6400-135

Clarity High Torque One Case Kits					
	Type	Hook	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
U	5 x 5	Cuspid	6700-420	6800-420	6400-420
U/L	5 x 5	Cuspid	6700-421	6800-421	6400-421
U/L	3 x 3	Cuspid	6700-422	6800-422	6400-422
U	3 x 3	Cuspid	6700-423	6800-423	6400-423
U	5 x 5		6700-424	6800-424	6400-424
U/L	5 x 5		6700-425	6800-425	6400-425
U/L	3 x 3		6700-426	6800-426	6400-426
U	3 x 3		6700-427	6800-427	6400-427

Clarity Standard Edgewise One Case Kits					
	Type	Hook	APC™ II .022 in.	APC™ PLUS .022 in.	Non-coated .022 in.
U/L	5 x 5		6700-320	6800-320	6400-320

6700	6800
Use this prefix to order APC™ II Adhesive Coated Appliance System.	Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Note
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

2




Radiused corners and edges enhance patient comfort and aesthetics.

Unique micro-crystalline bonding surface provides the consistency, strength and reliability of mechanical retention.

Underwire tie-wing protection ensures the labial tooth surface will never come in contact with AlastiK™ Ligatures.



Mechanical lock base provides reliable bond strength with all chemical and light cure bonding adhesive systems.



- Removable **color coded** indicators provide positive bracket identification and allow easier, accurate placement.
- Millimeter ledges aid in obtaining proper occlusal- gingival height and eliminate the need for Unitek™ Boone Gauge or other placement instruments.







APC™ II

Available in APC™ II
Adhesive Coated
Appliance System

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Selection By Tooth – Maxillary

	Torque	Ang.	In/Out	M/D Width	Hook	L/R	.018 in.	.022 in.	.018 in.	.022 in.
Central 	+22°	0°	.91	3.5		UL/UR	3300-650	-	6001-650	-
	+12°	5°	1.04	3.5		UL	3300-601	3300-701	6001-601	6001-701
	+12°	5°	1.04	3.5		UR	3300-602	3300-702	6001-602	6001-702
	0°	0°	.79	3.5		UL/UR	3300-801	3300-901	6001-801	6001-901
Lateral 	+14°	8°	.86	3.5		UL	3300-651	-	6001-651	-
	+14°	8°	.86	3.5		UR	3300-652	-	6001-652	-
	+8°	9°	1.42	3.5		UL	3300-603	3300-703	6001-603	6001-703
	+8°	9°	1.42	3.5		UR	3300-604	3300-704	6001-604	6001-704
	0°	0°	.79	3.5		UL/UR	3300-801	3300-901	6001-801	6001-901
Cuspid  	+7°	7°	.84	3.5		UL	3300-653	-	6001-653	-
	+7°	7°	.84	3.5		UR	3300-654	-	6001-654	-
	+7°	7°	.84	3.5	Hk	UL	3300-655	-	6001-655	-
	+7°	7°	.84	3.5	Hk	UR	3300-656	-	6001-656	-
	0°	11°	.94	3.5		UL	3300-607	3300-707	6001-607	6001-707
	0°	11°	.94	3.5		UR	3300-608	3300-708	6001-608	6001-708
	0°	0°	.79	3.5		UL/UR	3300-805	3300-905	6001-805	6001-905
	0°	11°	.94	3.5	Hk	UL	3300-605	3300-705	6001-605	6001-705
	0°	11°	.94	3.5	Hk	UR	3300-606	3300-706	6001-606	6001-706
Bicuspid  	-7°	0°	.94	3.5		UL/UR	3300-616	3300-716	6001-616	6001-716
	-7°	0°	.94	3.5	Hk	UL/UR	3300-615	3300-715	6001-615	6001-715
	-7°	0°	.94	2.5		UL/UR	3300-618	3300-718	6001-618	6001-718
	-7°	0°	.94	2.5	Hk	UL/UR	3300-617	3300-717	6001-617	6001-717
	0°	0°	.79	3.5		UL/UR	3300-805	3300-905	6001-805	6001-905

APC™ II Adhesive Coated Appliance System

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Note	
APC™ II Transcend™ Brackets sold in units of four.	3300-650 Appliance is available in APC™ II Adhesive Coated Appliance System.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
Roth* prescription brackets pictured in most cases.	













Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Selection By Tooth – Mandibular

2

	Torque	Ang.	In/Out	M/D Width	Hook	L/R	.018 in.	.022 in.	.018 in.	.022 in.
Anterior 	-1°	0°	1.42	2.5		LL/LR	3300-621	3300-721	6001-621	6001-721
	0°	0°	.79	2.5		LL/LR	3300-821	3300-921	6001-821	6001-921
Cuspid   	-11°	7°	.81	3.5		LL	3300-627	3300-727	6001-627	6001-727
	-11°	7°	.81	3.5		LR	3300-628	3300-728	6001-628	6001-728
	-11°	7°	.81	3.5	Hk	LL	3300-625	3300-725	6001-625	6001-725
	-11°	7°	.81	3.5	Hk	LR	3300-626	3300-726	6001-626	6001-726
	+7°	7°	.84	3.5		LL	3300-654	-	6001-654	-
	+7°	7°	.84	3.5		LR	3300-653	-	6001-653	-
	+7°	7°	.84	3.5	Hk	LL	3300-656	-	6001-656	-
	+7°	7°	.84	3.5	Hk	LR	3300-655	-	6001-655	-
	0°	0°	.79	3.5		LL/LR	3300-805	3300-905	6001-805	6001-905
Bicuspid 	0°	0°	.79	3.5		LL/LR	3300-805	3300-905	6001-805	6001-905
First Bicuspid   	-17°	0°	.86	3.5		LL/LR	3300-636	3300-736	6001-636	6001-736
	-17°	0°	.86	3.5	Hk	LL/LR	3300-635	3300-735	6001-635	6001-735
	-17°	0°	.86	2.5		LL/LR	3300-638	3300-738	6001-638	6001-738
	-17°	0°	.86	2.5	Hk	LL/LR	3300-637	3300-737	6001-637	6001-737
	0°	0°	.91	3.5		LL/LR	3300-805	-	6001-805	-
Second Bicuspid  	-22°	0°	.86	3.5		LL/LR	3300-644	3300-744	6001-644	6001-744
	-22°	0°	.86	3.5	Hk	LL/LR	3300-643	3300-743	6001-643	6001-743
	-22°	0°	.86	2.5		LL/LR	3300-646	3300-746	6001-646	6001-746
	-22°	0°	.86	2.5	Hk	LL/LR	3300-645	3300-745	6001-645	6001-745
	-17°	0°	.86	3.5		LL/LR	3300-636	-	6001-636	-
	APC™ II Adhesive Coated Appliance System									

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Note	
APC™ II Transcend™ Brackets sold in units of four.	3300-650 Appliance is available in APC™ II Adhesive Coated Appliance System.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
Roth* prescription brackets pictured in most cases.	

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Selection By Roth* Prescription

Maxillary	Torque	Ang.	Hook	L/R	.018 in.	.022 in.	.018 in.	.022 in.
Central	+12°	5°		UL	3300-601	3300-701	6001-601	6001-701
	+12°	5°		UR	3300-602	3300-702	6001-602	6001-702
Lateral	+8°	9°		UL	3300-603	3300-703	6001-603	6001-703
	+8°	9°		UR	3300-604	3300-704	6001-604	6001-704
Cuspid	0°	11°		UL	3300-607	3300-707	6001-607	6001-707
	0°	11°		UR	3300-608	3300-708	6001-608	6001-708
	0°	11°	Hk	UL	3300-605	3300-705	6001-605	6001-705
	0°	11°	Hk	UR	3300-606	3300-706	6001-606	6001-706
Bicuspid	-7°	0°		UL/UR	3300-616	3300-716	6001-616	6001-716
	-7°	0°	Hk	UL/UR	3300-615	3300-715	6001-615	6001-715
	-7°	0°		UL/UR	3300-618	3300-718	6001-618	6001-718
	-7°	0°	Hk	UL/UR	3300-617	3300-717	6001-617	6001-717
Mandibular								
Anterior	-1°	0°		LL/LR	3300-621	3300-721	6001-621	6001-721
Cuspid	-11°	7°		LL	3300-627	3300-727	6001-627	6001-727
	-11°	7°		LR	3300-628	3300-728	6001-628	6001-728
	-11°	7°	Hk	LL	3300-625	3300-725	6001-625	6001-725
	-11°	7°	Hk	LR	3300-626	3300-726	6001-626	6001-726
First Bicuspid	-17°	0°		LL/LR	3300-636	3300-736	6001-636	6001-736
	-17°	0°	Hk	LL/LR	3300-635	3300-735	6001-635	6001-735
	-17°	0°		LL/LR	3300-638	3300-738	6001-638	6001-738
	-17°	0°	Hk	LL/LR	3300-637	3300-737	6001-637	6001-737
Second Bicuspid	-22°	0°		LL/LR	3300-644	3300-744	6001-644	6001-744
	-22°	0°	Hk	LL/LR	3300-643	3300-743	6001-643	6001-743
	-22°	0°	NAR	LL/LR	3300-646	3300-746	6001-646	6001-746
	-22°	0°	Hk/NAR	LL/LR	3300-645	3300-745	6001-645	6001-745

APC™ II Adhesive Coated Appliance System


Victory Series™ Buccal Tubes Roth* Prescription

	Torque	Offset	L/R	APC™ II Adhesive Coated Appliance System		APC™ PLUS Adhesive Coated Appliance System		Page
				.018 in.	.022 in.	.018 in.	.022 in.	
Upper 1 ST	-10°	7°	L	068-805	068-905	068-805	068-905	4.4
	-10°	7°	R	068-806	068-906	068-806	068-906	
Upper 2 ND	-10°	6°	L	3 067-829 1	3 067-929 1	5 067-829 1	5 067-929 1	4.11
	-10°	6°	R	3 067-830 2	3 067-930 2	5 067-830 2	5 067-930 2	
Lower 1 ST	-25°	6°	L	3 067-801 3	3 067-901 3	5 067-801 3	5 067-901 3	4.8
	-25°	6°	R	3 067-802 4	3 067-902 4	5 067-802 4	5 067-902 4	
Lower 2 ND	-30°	6°	L	3 067-835 3	3 067-935 3	5 067-835 3	5 067-935 3	4.12
	-30°	6°	R	3 067-836 4	3 067-936 4	5 067-836 4	5 067-936 4	

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Note

APC™ II Transcend™ Brackets sold in units of four.	3300-650	Appliance is available in APC™ II Adhesive Coated Appliance System.
APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.		
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.	
For Transcend Series 6000 Roth* One/Ten Case Kits see page 2.24.	For detailed buccal tube information see the Buccal Tubes section.	
1 2 3 4 Buccal Tube is available with a bondable base.		

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Selection By High Torque/Ricketts* Prescription

2

Maxillary	Torque	Ang.	Hook	L/R	.018 in.	.022 in.	.018 in.	.022 in.
Central	+22°	0°		UL/UR	3300-650	-	6001-650	-
Lateral	+14°	8°		UL	3300-651	-	6001-651	-
	+14°	8°		UR	3300-652	-	6001-652	-
Cuspid	+7°	7°		UL	3300-653	-	6001-653	-
	+7°	7°		UR	3300-654	-	6001-654	-
	+7°	7°	Hk	UL	3300-655	-	6001-655	-
	+7°	7°	Hk	UR	3300-656	-	6001-656	-
Bicuspid	0°	0°		UL/UR	3300-805	-	6001-805	-
Mandibular								
Anterior	0°	0°		LL/LR	3300-821	-	6001-821	-
Cuspid	+7°	7°		LL	3300-654	-	6001-654	-
	+7°	7°		LR	3300-653	-	6001-653	-
	+7°	7°	Hk	LL	3300-656	-	6001-656	-
	+7°	7°	Hk	LR	3300-655	-	6001-655	-
First Bicuspid	0°	0°		LL/LR	3300-805	-	6001-805	-
Second Bicuspid	-17°	0°		LL/LR	3300-636	-	6001-636	-

APC™ II Adhesive Coated Appliance System

Victory Series™ Buccal Tubes Ricketts* Prescription

Torque	Offset	L/R	APC™ II Adhesive Coated Appliance System		APC™ PLUS Adhesive Coated Appliance System		Page
			.018 in.	.022 in.	.018 in.	.022 in.	
Upper 1 ST	0°	L	068-601	-	068-601	-	4.8
	0°	R	068-602	-	068-602	-	
Upper 2 ND	0°	L	3 067-825 1	-	5 067-825 1	-	4.11
Lower 1 ST	0°	R	3 067-826 2	-	5 067-826 2	-	
Lower 2 ND	0°	L/R	067-843	-	067-843	-	4.10
	0°	L/R	067-843	-	067-843	-	
	0°	L	3 067-826 3	-	5 067-826 3	-	4.11
	0°	R	3 067-825 4	-	5 067-825 4	-	

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.





Note	
APC™ II Transcend™ Brackets sold in units of four.	3300-650 Appliance is available in APC™ II Adhesive Coated Appliance System.
APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.	
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.
For Transcend Series 6000 High Torque/Ricketts* One Case Kits see page 2.25.	For detailed buccal tube information see the Buccal Tubes section.

1 2 3 4 Buccal Tube is available with a bondable base.

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets Selection By Standard Edgewise Prescription

	Torque	Ang.	Hook	L/R	.018 in.	.022 in.	.018 in.	.022 in.
Central Laterals 	0°	0°		UL/UR	3300-801	3300-901	6001-801	6001-901
Cuspids Bicuspids 	0°	0°		UL/UR	3300-805	3300-905	6001-805	6001-905
Anteriors 	0°	0°		LL/LR	3300-821	3300-921	6001-821	6001-921
Cuspids Bicuspids 	0°	0°		LL/LR	3300-805	3300-905	6001-805	6001-905

APC™ II Adhesive Coated Appliance System


Victory Series™ Buccal Tubes Roth* Prescription

	Torque	Offset	L/R	APC™ II Adhesive Coated Appliance System		APC™ PLUS Adhesive Coated Appliance System		Page
				.018 in.	.022 in.	.018 in.	.022 in.	
Upper 1 ST	0°	7°	L	068-801	068-901	068-801	068-901	4.4
	0°	7°	R	068-802	068-902	068-802	068-902	
Upper 2 ND	0°	6°	L	3 067-827 1	3 067-927 1	5 067-827 1	5 067-927 1	4.11
	0°	6°	R	3 067-828 2	3 067-928 2	5 067-828 2	5 067-928 2	
Lower 1 ST	0°	6°	L	3 067-809 3	3 067-909 3	5 067-809 3	5 067-909 3	4.8
	0°	6°	R	3 067-810 4	3 067-910 4	5 067-810 4	5 067-910 4	
Lower 2 ND	0°	6°	L	3 067-828 3	3 067-928 3	5 067-828 3	5 067-928 3	4.10
	0°	6°	R	3 067-827 4	3 067-927 4	5 067-827 4	5 067-927 4	

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Note		
APC™ II Transcend™ Brackets sold in units of four.	3300-650	Appliance is available in APC™ II Adhesive Coated Appliance System.
APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.		
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.	1 2 3 4 Buccal Tube is available with a bondable base.
For detailed buccal tube information see the Buccal Tubes section. 		

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Case Kits

2

Transcend™ Series 6000 Roth* One Case Kits						
	Type	Hook	.018 in.	.022 in.	.018 in.	.022 in.
Upper	3 x 3		3300-972	3300-992	-	-
Upper	3 x 3	Cuspid	3300-973	3300-993	-	-
U/L	3 x 3		3300-976	3300-996	-	-
U/L	3 x 3	Cuspid	3300-977	3300-997	-	-
U/L	5 x 5	Cuspid/Bicuspid	-	-	6010-601	6010-701

APC™ II Adhesive Coated Appliance System

Transcend Series 6000 Roth* One Case Kits (Narrow Bicuspid Brackets Included)						
	Type	Hook	.018 in.	.022 in.	.018 in.	.022 in.
Upper	5 x 5		3300-974	3300-994	-	-
Upper	5 x 5	Cuspid/Bicuspid	3300-975	3300-995	-	-
U/L	5 x 5		3300-978	3300-998	6010-602	6010-702
U/L	5 x 5	Cuspid/Bicuspid	3300-979	3300-999	-	-
U/L	5 x 5	Bicuspid	-	-	6010-603	6010-703
U/L	5 x 5		-	-	6010-604	6010-704

APC™ II Adhesive Coated Appliance System

Transcend Series 6000 Roth* Ten Case Kits						
	Type	Hook	.018 in.	.022 in.	.018 in.	.022 in.
U/L	5 x 5	Cuspid/Bicuspid	-	-	6010-610	6010-710
U/L	5 x 5		-	-	6010-612	6010-712

Note	
APC™ II Transcend™ Brackets sold in units of four.	3300-972 Appliance is available in APC™ II Adhesive Coated Appliance System.
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets


Case Kits



Transcend™ Series 6000 Standard Edgewise One Case Kits						
	Type	Hook	.018 in.	.022 in.	.018 in.	.022 in.
U/L	5 x 5		3300-969	3300-989	6010-801	6010-901
APC™ II Adhesive Coated Appliance System						

Transcend Series 6000 High Torque/Ricketts* One Case Kits						
	Type	Hook	.018 in.	.022 in.	.018 in.	.022 in.
U/L	5 x 5		3300-959	-	6010-802	-
APC™ II Adhesive Coated Appliance System						

Debonding Instrument:
See page 2.26 for proper
debonding instrument.

Note			
APC™ II Transcend™ Brackets sold in units of four.	3300-650	Appliance is available in APC™ II Adhesive Coated Appliance System.	
* 3M Unitek version of this prescription. No endorsement by the Doctor is implied.	Measurements are in millimeters unless otherwise noted.	1 2 3 4	Buccal Tube is available with a bondable base.
For detailed buccal tube information see the Buccal Tubes section. 			

Non-Metal Brackets

Transcend™ Series 6000 Ceramic Brackets

Transcend™ Series 6000 Debonding Instrument

2



Transcend™ Series 6000 Debonding Instrument

The Transcend Debonding instrument is designed exclusively for the debonding of Transcend™ Series 6000 Ceramic Brackets.

For exclusive use with Transcend ceramic brackets only.

444-770

CAUTION: Use of other instruments or techniques may increase the risk of enamel damage or other patient injury.