

Victory Series™ First Molar Buccal Tubes.....4.4
Victory Series™ Second Molar Buccal Tubes4.11
Victory Series™ Second Molar Mini Tubes.....4.13
Victory Series™ 2nd Bicuspid Tubes4.13
Light Wire Buccal Tubes.....4.14
Bondable Buccal Tubes4.15
Buccal Tubes on Bonding Bases4.16



Buccal Tubes

Victory Series™ Buccal Tubes

4



Our compact, technically advanced Victory Series™ Buccal Tubes include a patented, trouble-free integral convertible cap that converts easily and consistently, every time. All you have to do is lift and peel. This exclusive cap also offers solid, one piece construction, eliminating the possibility of braze or solder materials that can leave a residue in the archwire slot.

Victory Series Buccal Tubes are designed with funneled and radiused entries on both the headgear tube and archwire slot, resulting in easy archwire insertion which can reduce patient chair time. Smooth, patient-friendly hooks improve comfort for your patients. You'll also appreciate the strong construction, low-profile design and distal notch that secures tie backs.



Victory Series Buccal Tubes. Reducing chair time, increasing patient comfort and providing you with a well-designed, reliable and predictable orthodontic solution.



Buccal Tubes

Victory Series™ Molar Bonding Bases



Now you can achieve a stronger, more consistent bond in every case with the new Victory Series™ Molar Bonding Base. Look closely at its profile and you'll see a curved anatomical design for optimal tube-to-tooth fit and a buccal indent that naturally conforms to the surface of the molar, making it easier to correctly position the buccal tube in place.

The micro-etched 80 gauge mesh, along with the buccal indent and the correct anatomical fit, are designed for increased bond strength. The buccal indent assists in more accurate tube positioning; a key advantage that can eliminate time consuming archwire adjustments. And with the Victory Series molar bonding base, you'll never have to worry about tooth separation for band placement.

For added convenience, two different sizes accommodate a wide array of buccal attachments. The Victory Series molar bonding base from 3M Unitek. Strong molar bonding is our strength.



APC™ II

Available in APC™ II
Adhesive Coated
Appliance System

APC™ PLUS

Available in APC™ PLUS
Adhesive Coated
Appliance System

Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Upper Convertible Double Combination



4

Slot Torque Arch	Distal Offset	Headgear	Length	L/R	.018 Slot	.022 Slot
0°	0°	.045	3.6 mm	UL	068-841	068-941
0°	0°	.045	3.6 mm	UR	068-842	068-942
0°	0°	.051	3.6 mm	UL	068-843	068-943
0°	0°	.051	3.6 mm	UR	068-844	068-944
0°	7°	.045	4.3 mm	UL	068-801	068-901
0°	7°	.045	4.3 mm	UR	068-802	068-902
0°	7°	.051	4.3 mm	UL	068-803	068-903
0°	7°	.051	4.3 mm	UR	068-804	068-904
0°	7°	.045	3.6 mm	UL	068-875	068-975
0°	7°	.045	3.6 mm	UR	068-874	068-974
-10°	7°	.045	4.3 mm	UL	068-805	068-905
-10°	7°	.045	4.3 mm	UR	068-806	068-906
-10°	7°	.051	4.3 mm	UL	068-807	068-907
-10°	7°	.051	4.3 mm	UR	068-808	068-908
-10°	7°	.045	3.6 mm	UL	068-871	068-971
-10°	7°	.045	3.6 mm	UR	068-870	068-970
-10°	7°	.051	3.6 mm	UL	068-873	068-973
-10°	7°	.051	3.6 mm	UR	068-872	068-972
-10°	14°	.045	4.3 mm	UL	068-809	068-909
-10°	14°	.045	4.3 mm	UR	068-810	068-910
-10°	14°	.051	4.3 mm	UL	068-811	068-911
-10°	14°	.051	4.3 mm	UR	068-812	068-912
-10°	14°	.045	3.6 mm	UL	068-879	068-979
-10°	14°	.045	3.6 mm	UR	068-878	068-978
-10°	14°	.051	3.6 mm	UL	068-881	068-981
-10°	14°	.051	3.6 mm	UR	068-880	068-980

1 2 3 4
UL UR LL LR

Suffix indicates Buccal Tubes are available with bonding bases.

3 Prefix indicates Buccal Tube, with bonding base, is available in APC™ II Adhesive Coated Appliance System.

5 Prefix indicates Buccal Tube, with bonding base, is available in APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

Buccal Tube is not available with these options when blank.

Compatible with the Forsus™ Fatigue Resistant Device EZ Module



Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Upper Convertible Single



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot			MBT™		
Arch					3	5	1	3	5	1		2	
-14°		10°	3.6 mm	UL	3	5	068-897	1	3	5	068-997	1	MBT™
-14°		10°	3.6 mm	UR	3	5	068-898	2	3	5	068-998	2	MBT™

APC™ II
 APC™ PLUS
 Base
 APC™ II
 APC™ PLUS
 Base

Upper Convertible Double Combination With Gingival Headgear Tube



Slot Torque		Distal Offset	Headgear	Length	L/R	.018 Slot		.022 Slot		MBT™
Arch						068-821	068-822	068-921	068-922	
-10°		7°	.045	4.3 mm	UL	068-821	068-822	068-921	068-922	
-10°		7°	.045	4.3 mm	UR	068-821	068-822	068-921	068-922	
-14°		10°	.045	4.3 mm	UL	068-885	068-886	068-985	068-986	MBT™
-14°		10°	.045	4.3 mm	UR	068-885	068-886	068-985	068-986	MBT™
-14°		10°	.051	4.3 mm	UL	068-887	068-888	068-987	068-988	MBT™
-14°		10°	.051	4.3 mm	UR	068-887	068-888	068-987	068-988	MBT™

Upper Convertible Double Combination - Occlusal Headgear



Slot Torque		Distal Offset	Headgear	Length	L/R	.018 Slot		.022 Slot		MBT™
Arch						068-893	068-894	068-993	068-994	
-14°		10°	.045	4.6 mm	UL	068-893	068-894	068-993	068-994	MBT™
-14°		10°	.045	4.6 mm	UR	068-893	068-894	068-993	068-994	MBT™
-14°		10°	.051	4.6 mm	UL	068-895	068-896	068-995	068-996	MBT™
-14°		10°	.051	4.6 mm	UR	068-895	068-896	068-995	068-996	MBT™

Compatible with the Forsus™ Fatigue Resistant Device EZ Module



Note	
3 Use this prefix to order APC™ II Adhesive Coated Appliance System.	5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.
APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.	Measurements are in millimeters unless otherwise noted.
1 2 3 4 Buccal Tube is available with a bondable base.	

Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Upper Convertible Triple Combination - Occlusal Headgear



Slot Torque		Distal Offset	Headgear	Length	L/R	.018 Slot	.022 Slot
Arch	Aux						
0°	0°	0°	.045	3.6 mm	UL	068-845	068-945
0°	0°	0°	.045	3.6 mm	UR	068-846	068-946
0°	0°	0°	.051	3.6 mm	UL	068-847	068-947
0°	0°	0°	.051	3.6 mm	UR	068-848	068-948
0°	0°	15°	.045	4.0 mm	UL	068-825	-
0°	0°	15°	.045	4.0 mm	UR	068-826	-
-10°	0°	7°	.045	3.6 mm	UL	068-813	068-913
-10°	0°	7°	.045	3.6 mm	UR	068-814	068-914
-10°	0°	7°	.051	3.6 mm	UL	068-815	068-915
-10°	0°	7°	.051	3.6 mm	UR	068-816	068-916
-10°	0°	14°	.045	3.6 mm	UL	068-817	068-917
-10°	0°	14°	.045	3.6 mm	UR	068-818	068-918
-10°	0°	14°	.051	3.6 mm	UL	068-819	068-919
-10°	0°	14°	.051	3.6 mm	UR	068-820	068-920
-14°	0°	10°	.045	4.3 mm	UL	068-889	068-989
-14°	0°	10°	.045	4.3 mm	UR	068-890	068-990
-14°	0°	10°	.051	4.3 mm	UL	068-891	068-991
-14°	0°	10°	.051	4.3 mm	UR	068-892	068-992

MBT™
MBT™
MBT™
MBT™

Compatible with the
Forsus™ Fatigue Resistant Device Module



Upper Convertible Triple Combination Differentiated Slot



Slot Torque		Distal Offset	Headgear			Length	L/R	.018 Slot	.022 Slot
Arch	Aux		Main	Aux	Aux				
-10°	0°	14°	.045	.018	.022	3.6 mm	UL	068-853	-
-10°	0°	14°	.045	.018	.022	3.6 mm	UR	068-854	-
-10°	0°	7°	.045	.022	.018	3.6 mm	UL	-	068-959
-10°	0°	7°	.045	.022	.018	3.6 mm	UR	-	068-960
0°	0°	0°	.045	.022	.018	3.6 mm	UL	-	068-957
0°	0°	0°	.045	.022	.018	3.6 mm	UR	-	068-958

Compatible with the
Forsus™ Fatigue Resistant Device EZ Module




Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Upper Convertible Single



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch					3	5	1	2	3	5	1	2
0°		0°	3.6 mm	UL	3	5	067-808	1	3	5	067-908	1
0°		0°	3.6 mm	UR	3	5	067-807	2	3	5	067-907	2
0°		6°	3.6 mm	UL	3	5	067-810	1	3	5	067-910	1
0°		6°	3.6 mm	UR	3	5	067-809	2	3	5	067-909	2
-10°		7°	4.3 mm	UL	3	5	068-823	1	3	5	068-923	1
-10°		7°	4.3 mm	UR	3	5	068-824	2	3	5	068-924	2



Upper 1ST Non-Convertible Single - Funneled

Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch	Aux				3	5	1	2	3	5	1	2
-14°		10°	4.6 mm	UL	3	5	067-879	1	3	5	067-979	1
-14°		10°	4.6 mm	UR	3	5	067-878	2	3	5	067-978	2

MBT™
MBT™

Upper Convertible Double



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch	Aux				3	5	1	2	3	5	1	2
0°	0°	0°	3.6 mm	UL	3	5	068-849	1	3	5	068-949	1
0°	0°	0°	3.6 mm	UR	3	5	068-850	2	3	5	068-950	2
-10°	0°	7°	3.6 mm	UL	3	5	068-851	1	3	5	068-951	1
-10°	0°	7°	3.6 mm	UR	3	5	068-852	2	3	5	068-952	2

Upper Convertible Double Differentiated Slot



Slot Torque		Distal Offset	Main	Aux	Length	L/R	.018 Slot	.022 Slot
Arch	Aux						1	2
0°	0°	0°	.022	.018	3.6 mm	UL/LR	-	068-953
0°	0°	0°	.022	.018	3.6 mm	UR/LL	-	068-954

Note	3 Use this prefix to order APC™ II Adhesive Coated Appliance System.	5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.
APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.	Measurements are in millimeters unless otherwise noted.	1 2 3 4 Buccal Tube is available with a bondable base.



Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Upper Non-Convertible Double Combination



Slot Torque		Distal Offset	Headgear	Length	L/R	.018 Slot	.022 Slot
Arch							
0°		0°	.045	5.1 mm	UL	068-861	068-961
0°		0°	.045	5.1 mm	UR	068-862	068-962
0°		0°	.051	5.1 mm	UL	068-863	068-963
0°		0°	.051	5.1 mm	UR	068-864	068-964
0°		6°	.045	5.1 mm	UL	068-865	068-965
0°		6°	.045	5.1 mm	UR	068-866	068-966
0°		6°	.051	5.1 mm	UL	068-867	068-967
0°		6°	.051	5.1 mm	UR	068-868	068-968

Upper Non-Convertible Triple Combination



Slot Torque		Distal Offset	Headgear	Hook	Length	L/R	.018 Slot	.022 Slot
Arch	Aux							
0°	0°	0°	.045	CTR	4.3 mm	UNIV	068-882	068-982
0°	0°	0°	.051	CTR	4.3 mm	UNIV	068-883	068-983
0°	0°	0°	.045	MES	4.3 mm	UL	068-601	**
0°	0°	0°	.045	MES	4.3 mm	UR	068-602	**



Compatible with the Forsus™ Fatigue Resistant Device EZ Module

Lower Convertible Single



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch												
0°		0°	3.6 mm	LL	3	5	067-807	3	3	5	067-907	3
0°		0°	3.6 mm	LR	3	5	067-808	4	3	5	067-908	4
0°		6°	3.6 mm	LL	3	5	067-809	3	3	5	067-909	3
0°		6°	3.6 mm	LR	3	5	067-810	4	3	5	067-910	4
-10°		0°	4.3 mm	LR	3	5	067-813	4	3	5	067-913	4
-10°		0°	4.3 mm	LL	3	5	067-814	3	3	5	067-914	3
-25°		6°	3.6 mm	LL	3	5	067-851	3	3	5	067-951	3
-25°		6°	3.6 mm	LR	3	5	067-850	4	3	5	067-950	4
-25°		6°	4.3 mm	LL	3	5	067-801	3	3	5	067-901	3
-25°		6°	4.3 mm	LR	3	5	067-802	4	3	5	067-902	4
-20°		0°	4.3 mm	LL	3	5	067-819	3	3	5	067-919	3
-20°		0°	4.3 mm	LR	3	5	067-820	4	3	5	067-920	4
-25°		0°	4.3 mm	LL	3	5	067-817	3	3	5	067-917	3
-25°		0°	4.3 mm	LR	3	5	067-818	4	3	5	067-918	4

MBT™
MBT™

Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Lower 1ST Non-Convertible Single - Funneled



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot		
Arch					APC™ II	Part No.	Base	APC™ II	Part No.	Base
-20°		0°	4.6 mm	LL	3	067-877	3	3	067-977	3
-20°		0°	4.6 mm	LR	3	067-876	4	3	067-976	4

Lower Non-Convertible Single



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot			MBT™		
Arch					APC™ II	APC™ PLUS	Base	APC™ II	APC™ PLUS	Base		MBT™	
-20°		0°	4.3 mm	LL	3	5	067-867	3	3	5	067-967	3	MBT™
-20°		0°	4.3 mm	LR	3	5	067-868	4	3	5	067-968	4	MBT™

Lower Convertible Double



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot			MBT™		
Arch	Aux				APC™ II	APC™ PLUS	Base	APC™ II	APC™ PLUS	Base		MBT™	
0°	0°	0°	3.6 mm	LL	3	5	068-850	3	3	5	068-950	3	MBT™
0°	0°	0°	3.6 mm	LR	3	5	068-849	4	3	5	068-949	4	MBT™
-25°	0°	7°	3.6 mm	LL	3	5	067-803	3	3	5	067-903	3	MBT™
-25°	0°	7°	3.6 mm	LR	3	5	067-804	4	3	5	067-904	4	MBT™
-20°	0°	0°	4.3 mm	LL	3	5	067-861	3	3	5	067-961	3	MBT™
-20°	0°	0°	4.3 mm	LR	3	5	067-860	4	3	5	067-960	4	MBT™

Note

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.

Measurements are in millimeters unless otherwise noted.

1 2 3 4 Buccal Tube is available with a bondable base.



Buccal Tubes

Victory Series™ First Molar Buccal Tubes

Lower Convertible Double Differentiated Slot

Slot Torque		Distal Offset	Main	Aux	Length	L/R	.018 Slot	.022 Slot
Arch	Aux							
-25°	0°	7°	.022	.018	3.6 mm	LL	-	067-955
-25°	0°	7°	.022	.018	3.6 mm	LR	-	067-956

4

Lower Lip Bumper Convertible Double



Slot Torque		Distal Offset	Lip Bumper	Length	L/R	.018 Slot	.022 Slot
Arch							
-25°		5°	.045	4.3 mm	LL	061-251	061-351
-25°		5°	.045	4.3 mm	LR	061-252	061-352

Lower Lip Bumper Convertible Triple



Slot Torque		Distal Offset	Lip Bumper	Length	L/R	.018 Slot	.022 Slot
Arch	Aux						
-25°	0°	5°	.045	4.3 mm	LL	067-811	067-911
-25°	0°	5°	.045	4.3 mm	LR	067-812	067-912

Upper/Lower Non-Convertible Double



Slot Torque		Distal Offset	Length	L/R	.018 Slot	.022 Slot
Arch	Aux					
0°	0°	0°	4.3 mm	UNIV	067-843	067-943

Buccal Tubes

Victory Series™ Second Molar Buccal Tubes

Upper Second Non-Convertible Single with Distal Notch



Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch					3	5	1	2	3	5	1	2
-10°		6°	3.6 mm	UL	3	5	067-829	1	3	5	067-929	1
-10°		6°	3.6 mm	UR	3	5	067-830	2	3	5	067-930	2
-10°		10°	3.6 mm	UL	3	5	067-831	1	3	5	067-931	1
-10°		10°	3.6 mm	UR	3	5	067-832	2	3	5	067-932	2
-14°		10°	4.3 mm	UL	3	5	067-857	1	3	5	067-957	1
-14°		10°	4.3 mm	UR	3	5	067-858	2	3	5	067-958	2

MBT™
MBT™



Upper Second Non-Convertible Single - Funneled



Slot Torque		Distal Offset	Length	L/R	.018 Slot		.022 Slot			
Arch					3	5	3	5		
-19°		8°	3.5 mm	UL	3	067-873	1	3	067-973	1
-19°		8°	3.5 mm	UR	3	067-872	2	3	067-972	2

MBT™/Sondhi™
MBT™/Sondhi™

Upper/Lower Second Non-Convertible Single with Distal Notch

Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot				
Arch					3	5	1	2	3	5	1	2
0°		0°	3.6 mm	UL	3	5	067-825	1	3	5	067-925	1
0°		0°	3.6 mm	UR	3	5	067-826	2	3	5	067-926	2
0°		0°	3.6 mm	LL	3	5	067-826	3	3	5	067-926	3
0°		0°	3.6 mm	LR	3	5	067-825	4	3	5	067-925	4
0°		6°	3.6 mm	UL	3	5	067-827	1	3	5	067-927	1
0°		6°	3.6 mm	UR	3	5	067-828	2	3	5	067-928	2
0°		6°	3.6 mm	LL	3	5	067-828	3	3	5	067-928	3
0°		6°	3.6 mm	LR	3	5	067-827	4	3	5	067-927	4

Note

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.

Measurements are in millimeters unless otherwise noted.

1 2 3 4 Buccal Tube is available with a bondable base.

Buccal Tubes

Victory Series™ Second Molar Buccal Tubes

Lower Second Non-Convertible Single with Distal Notch

Slot Torque		Distal Offset	Length	L/R	.018 Slot			.022 Slot			MBT™	
Arch					3	5		3	5			3
-10°		0°	4.3 mm	LL	3	5	067-816	3	3	5	067-916	3
-10°		0°	4.3 mm	LR	3	5	067-815	4	3	5	067-915	4
-25°		6°	3.6 mm	LL	3	5	067-837	3	3	5	067-937	3
-25°		6°	3.6 mm	LR	3	5	067-838	4	3	5	067-938	4
-25°		10°	3.6 mm	LL	3	5	067-839	3	3	5	067-939	3
-25°		10°	3.6 mm	LR	3	5	067-840	4	3	5	067-940	4
-30°		6°	3.6 mm	LL	3	5	067-835	3	3	5	067-935	3
-30°		6°	3.6 mm	LR	3	5	067-836	4	3	5	067-936	4

APC™ II
 APC™ PLUS
 Base
 APC™ II
 APC™ PLUS
 Base

Lower Second Non-Convertible Single with Distal Notch and Hook

Slot Torque		Distal Offset	Length	L/R	.018 Slot	.022 Slot
Arch					067-844	067-944
0°		0°	5.1 mm	UL/LR	067-844	067-944
0°		0°	5.1 mm	UR/LL	067-845	067-945
0°		6°	5.1 mm	UL/LR	067-846	067-946
0°		6°	5.1 mm	UR/LL	067-847	067-947

Note

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.

Measurements are in millimeters unless otherwise noted.

1 2 3 4 Buccal Tube is available with a bondable base.

Buccal Tubes

Victory Series™ Second Molar Mini Tubes

Victory Series™ 2nd Bicuspid Tubes

Buccal Tubes For Partially Erupted Second Molars



Upper Victory Series™ Mini Tube

Slot Torque Arch		Distal Offset	Length	L/R			.018 Slot		.022 Slot	
0°					0°	2.5 mm	UR	3	5	066-4092
0°		0°	2.5 mm	UL	3	5	066-4101	3	5	066-5101
-10°		0°	2.5 mm	UR	3	5	066-4112	3	5	066-5112
-10°		0°	2.5 mm	UL	3	5	066-4121	3	5	066-5121
-14°		0°	2.5 mm	UR	3	5	066-4082	3	5	066-5082
-14°		0°	2.5 mm	UL	3	5	066-4071	3	5	066-5071

MBT™
MBT™



Lower Victory Series™ Mini Tube

Slot Torque Arch		Distal Offset	Length	L/R			.018 Slot		.022 Slot	
0°					0°	2.5 mm	LR	3	5	066-4024
0°		0°	2.5 mm	LL	3	5	066-4013	3	5	066-5013
-10°		0°	2.5 mm	LR	3	5	066-4044	3	5	066-5044
-10°		0°	2.5 mm	LL	3	5	066-4033	3	5	066-5033
-25°		0°	2.5 mm	LR	3	5	066-4064	3	5	066-5064
-25°		0°	2.5 mm	LL	3	5	066-4053	3	5	066-5053

MBT™
MBT™

These tubes are direct bond only, not available welded to bands.

Victory Series™ 2nd Bicuspid Tubes



3	5	058-801	0	Lower Left 2nd Bicuspid Tube .018 Slot	(-17° torque, 2° angulation)	MBT™
3	5	058-802	0	Lower Right 2nd Bicuspid Tube .018 Slot	(-17° torque, 2° angulation)	MBT™
3	5	058-803	0	Upper Left / Right 2nd Bicuspid Tube .018 Slot	(-7° torque, 0° angulation)	MBT™
3	5	058-901	0	Lower Left 2nd Bicuspid Tube .022 Slot	(-17° torque, 2° angulation)	MBT™
3	5	058-902	0	Lower Right 2nd Bicuspid Tube .022 Slot	(-17° torque, 2° angulation)	MBT™
3	5	058-903	0	Upper Left / Right 2nd Bicuspid Tube .022 Slot	(-7° torque, 0° angulation)	MBT™

APC™ II
APC™ PLUS

Offset Erase

For more information about Victory Series 2nd Bicuspid Tubes – see page 3.50.

Note

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

APC™ II Buccal Tubes sold in units of five, APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.


Measurements are in millimeters unless otherwise noted.

1 2 3 4 Buccal Tube is available with a bondable base.

Buccal Tubes


Light Wire Buccal Tubes

Upper First Molars with Vertical Slot




V-Slot Size	Distal Offset	Light Wire Tube I.D.	Headgear I.D.	Headgear Length	Light Wire Length	L/R	Order Number
.021 x .021	0°	.036	.045	4.5 mm	6.5 mm	UL	046-001
.021 x .021	0°	.036	.045	4.5 mm	6.5 mm	UR	046-002

Upper First Molars




V-Slot Size	Distal Offset	Light Wire Tube I.D.	Headgear I.D.	Headgear Length	Light Wire Length	L/R	Order Number
-	0°	.036	.045	4.5 mm	4.5 mm	UNIV	047-012 No Hook
.021 x .021	6°	.036	.045	4.5 mm	6.5 mm	UL	046-011
.021 x .021	6°	.036	.045	4.5 mm	6.5 mm	UR	046-012
-	6°	.036	.045	4.5 mm	6.5 mm	UL	046-013
-	6°	.036	.045	4.5 mm	6.5 mm	UR	046-014

Upper/Lower First Molars with Vertical Slot



V-Slot Size	Distal Offset	Light Wire Tube I.D.	Headgear I.D.	Headgear Length	Light Wire Length	L/R	Order Number
.021 x .021	0°	.036			6.5 mm	UL/LR	046-146
.021 x .021	0°	.036			6.5 mm	LL/UR	046-147
.021 x .021	6°	.036			4.5 mm	UL/LR	046-174
.021 x .021	6°	.036			4.5 mm	LL/UR	046-175
.021 x .021	6°	.036			6.5 mm	UL/LR	046-176
.021 x .021	6°	.036			6.5 mm	LL/UR	046-177

Upper/Lower First Molars





Distal Offset	Light Wire Tube I.D.	Headgear I.D.	Headgear Length	Light Wire Length	L/R	Order Number
0°	.036			4.5 mm	UL/LR	046-114
					LL/UR	046-115
0°	.036			6.5 mm	UL/LR	046-116
					LL/UR	046-117

All Light Wire Buccal Tubes have hooks unless specified.

Buccal Tubes

Bondable Buccal Tubes

First Molar Bonding Bases (#80 Gauge Mesh)




Buccal Tubes


Buccal Tubes on Bonding Bases

4.3 mm Single Convertible Buccal Tube


Bonding Base
019-434



Bonding Base
019-325




Bonding Base
019-326




3.6 mm Single Convertible Buccal Tube

Bonding Base
019-435



Bonding Base
019-325



Bonding Base
019-326

