



Transbond™ Plus Light Cure Band Adhesive

*No mixing
No contamination.
No waste.*

*A convenient, single paste, no mix
band adhesive from 3M Unitek.*

It should come as no surprise that the company that revolutionized orthodontic adhesives with APC™ Adhesive Coated Appliances also brings you a single paste band adhesive. Once again, we've simplified the process. Transbond™ Plus Light Cure Band Adhesive can maximize your productivity by helping you minimize chair time, reduce adhesive waste and by providing superior band retention.

3M Unitek

Products that make your life easier.

Make the smart move to Transbond™ Plus Light Cure Band Adhesive.



A faster, safer, more predictable solution.

As a single paste, the major advantage of Transbond™ Plus light cure band adhesive is the elimination of mixing, resulting in reduced clean-up, more consistent adhesive performance, increased office efficiency, and less chair time. The adhesive kit comes with five convenient syringes and 50 disposable angulated tips. We've eliminated the mixing pads and concerns about cross-contamination, because you apply the adhesive directly onto the band. Plus, you and your patients will rest easier knowing that the adhesive releases fluoride.



Function you can see.

Unlike other designs, angulated Transbond™ Plus Light Cure Band Adhesive disposable tips apply adhesive directly to the band, reducing waste and the chance for cross-contamination.

The distinctive blue color of Transbond Plus light cure band adhesive is another smart feature. It allows you to easily see the adhesive during banding, making clean-up of excess adhesive easy. During de-banding, you can see any adhesive left on the teeth, which is clearly easier on your staff as well as your patients. Used with a high performance light cure unit, such as the Ortholux™ LED Curing Light, cure time is a mere 20 seconds.

Contact 3M Unitek today and discover more about Transbond Plus Light Cure Band Adhesive. It's the smartest call you can make to enhance efficiency and patient satisfaction.



3M

**3M Unitek
Orthodontic Products**
2724 South Peck Road
Monrovia, CA 91016 USA
www.3MUnitek.com

In U.S. and Puerto Rico: 1-800-423-4588 • 626-574-4000
In Canada: 1-800-443-1661
Technical Helpline: 1-800-265-1943 • 626-574-4577
CE Hotline: 1-800-852-1990 x 4649 • 626-574-4649
Outside these areas, contact your local representative.

Please recycle. Printed in U.S.A.
© 1999-2008 3M. All rights reserved.
REF 016-910-3 0801