

Victory Series™ Molar Bonding Base

*Secure a stronger
bond, on contact.*



Now you can achieve
a stronger, more
consistent bond in every
case with the new Victory

Series™ Molar Bonding Base. Look closely
at its profile and you'll see a curved anatomical design for
optimal tube-to-tooth fit and a buccal indent that naturally
conforms to the surface of the molar, making it easier to
correctly position the buccal tube in place.

The micro-etched 80 gauge mesh, along with the buccal
indent and the correct anatomical fit, are designed for
significantly increased bond strength. The buccal indent
assists in more accurate tube positioning; a key advantage
that can eliminate the need for time consuming archwire
adjustments. And with the Victory Series molar bonding
base, you'll never have to worry about tooth separation for
band placement.

For added convenience, two different sizes accommodate
a wide array of buccal attachments.

The new Victory Series molar bonding base from
3M Unitek. Superior molar bonding is our strength.



*Optimal base fit for increased bond
strength and reduced occlusal interference.
Color coded ID dots make identification
quick and easy.*

3M Unitek
Products that make your life easier.

Victory Series™ Buccal Tubes

Our compact, technically advanced Victory Series Buccal Tubes include a patented, trouble-free integral convertible cap that converts easily and consistently, every time. All you have to do is lift and peel. This exclusive cap also offers solid, one piece construction, eliminating the possibility of braze or solder materials that can leave a residue in the archwire slot. Victory Series Buccal Tubes are designed with funneled and radiused entries on both the headgear tube and archwire slot, resulting in easy archwire insertion which can reduce patient chairtime. Smooth, patient-friendly hooks improve comfort for your patients. You'll also appreciate the strong construction, low-profile design and distal notch that secures tie backs.



Victory Series Buccal Tubes. Reducing chair time, increasing patient comfort and providing you with a well-designed, reliable and predictable orthodontic solution.



3M Unitek's patented convertible buccal tube cap removing pliers for quick, smooth cap removal.



To peel the convertible cap away, close the pliers handles slowly and steadily, driving the slide bar tip into the archwire slot. The buccal tube is now a molar bracket.



For partial cap conversion, securely grasp the cap with utility pliers and rotate distally to break away the attachment.

3M

**3M Unitek
Dental Products Division**

2724 South Peck Road
Monrovia, CA 91016

In U.S. and Puerto Rico: 1-800-423-4588
In Canada: 1-800-443-1661
In Mexico: 95-800-833-3215
Technical Hotline: 1-800-265-1943
Outside these areas, contact your local representative.



Printed on 50% recycled waste paper, including 10% post-consumer waste paper.

Unitek and Victory Series are trademarks of 3M Unitek Corporation
©1995 3M Unitek Corporation 16-881 9505