

ORTHO IMPLANT™

Exclusively Distributed by
3M Unitek



3.5mm
LT035



10.5mm
LT105



17.5mm
LT175

ORTHO Locking Closed Coil Springs*

- **Unique locking eyelet** allows for **ligature-free**, locking attachment of the spring to the head of the implant.
- Available in three lengths: 3mm (LCC3-10), 6mm (LCC6-10) & 8mm (LCC8-10).

*Patent Pending

LT-Driver

- **Access** in areas where the straight handle of the O-Driver is not ideal, i.e. the lateral palatal wall, depth of the palatal vault, retromolar area and external oblique ridge.
- **Manual contra angle** design makes placement comfortable and simple.
- **ORTHO LT Driver Kit** includes ORTHO LT Driver and 3.5mm ORTHO LT Driver Adapter (Item LT-KIT).

LT-Adapters

- LT Adapters allow for direct engagement of the ORTHO Implant and are available in lengths of **3.5mm, 10.5mm and 17.5mm** for easy access to multiple areas.
- **Laser markings** on all LT Adapters correspond with the location of the attachment holes on the implant head, allowing the doctor to know where the attachment holes are positioned without removing the instrument from the implant.

Slides over implant head...



locks into place



3mm
LCC3-10



6mm
LCC6-10



8mm
LCC8-10