



Littmann®
Stethoscopes

Classic

Love your stethoscope.

The compact and sensitive **3M™ Littmann® Classic III™ Monitoring Stethoscope** is the perfect instrument for clinicians who monitor and assess children or adults in non-critical care environments. Manufactured with strong yet lightweight materials for hours of comfortable use.

Monitor and assess a wide range of patients.

- Useful in non-critical care environments such as a medical office, general ward, OB/GYN, ambulatory clinic or urgent care
- Fun and vibrant; match color and finish to your personality
- Designed for use with adult and pediatric patients
- Tunable, dual-sided stainless steel chestpiece with open or closed bell
- Non-chill bell sleeve for greater patient comfort
- Soft-sealing eartips provide optimal comfort and excellent sound occlusion
- Made in the USA without natural rubber latex



Product Features

| | | |
|-------------------------------------|---|--|
| Acoustic Performance | 8 (1–10 scale) | |
| Chestpiece Design/Technology | Machined stainless steel, two-sided | |
| Headset Material | Lightweight, aircraft-grade, aluminum alloy | |
| Approximate Length | 27"/69 cm | |
| Approximate Weight | 150 grams | |
| Diaphragm Details/Diameter | Two tunable diaphragms 1.7"/4.3 cm, 1.4"/3.4cm | |
| Warranty | 5 years | |

3M Health Care
3M Medical Solutions Division
2510 Conway Ave., St. Paul, MN 55144 USA
1-800-228-3957 | Littmann.com

3M, Littmann, the L Littmann logo, Classic III, are trademarks of 3M. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2019. All rights reserved. 70-2011-7198-3