

Meet the next generation  
of a clinical icon.



### Open bell versatility.

The pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

## 3M™ Littmann® Classic III™ Stethoscope

**Right out of the box, it lives up to its name.**

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The Classic III stethoscope offers you:



**High acoustic sensitivity** for performing general physical assessments.



**Tunable diaphragms** on both the adult and pediatric sides of the chestpiece.



**Single-piece diaphragm** that's easy to attach and easy to clean because the surface is smooth without crevices.



**Next-generation tubing** with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.

# Here's what makes the 3M™ Littmann® Classic III™ Stethoscope stand out.

## Features

**Tunable technology** is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

**Soft-sealing eartips** give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

**Open bell stays clear** of dirt and debris by keeping it covered with the small diaphragm.

**Adjustable headset** lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

**No natural rubber latex or phthalate plasticizers** are used in the tubing or any other component—to help protect human health and the environment.










**Made in the USA** for consistent quality. We apply rigorous process and quality controls to every Littmann stethoscope. You'll always get a durable, dependable, clinical tool that will serve you for years to come.

Classic III stethoscope specifications	
Length	69cm / 27.2 inches
Weight	150g
Acoustic performance	7
Warranty	5 years
<b>CHESTPIECE</b>	
Construction	Two-sided, adult and pediatric
Material	Stainless steel
Diaphragm type	Tunable, single-piece
Diaphragm diameter	Adult: 4.3cm / 1.7 in Pediatric: 3.3cm / 1.3 in
<b>TUBING</b>	
Construction	Single lumen (single acoustic channel within tube)
Material	Not made with natural rubber latex or phthalate plasticizers
<b>HEADSET</b>	
Material	Wide diameter aerospace alloy / anodized aluminum

## Available tubing colors and finish combinations

Black	Burgundy	Navy Blue	Plum	Chocolate	Gray	Ceil Blue	Raspberry	Caribbean Blue	Pearl Pink	Dark Olive Green	Lavender	Turquoise	Lemon-Lime	Rose Pink
5620 <sup>1</sup> 5803 <sup>9</sup> 5811 <sup>5</sup> 5861 <sup>15</sup> 5870 <sup>12</sup>	5627 <sup>1</sup> 5864 <sup>15</sup> 5868 <sup>9</sup>	5622 <sup>1</sup> 5863 <sup>13</sup> 5867 <sup>9</sup>	5831 <sup>1</sup>	5809 <sup>4</sup>	5621 <sup>1</sup>	5630 <sup>1</sup>	5806 <sup>7</sup> 5862 <sup>13</sup> 5871 <sup>6</sup>	5623 <sup>1</sup> 5807 <sup>7</sup> 5869 <sup>9</sup>	5633 <sup>1</sup>	5812 <sup>5</sup>	5832 <sup>1</sup> 5865 <sup>13</sup>	5835 <sup>1</sup>	5839 <sup>1</sup>	5631 <sup>1</sup> *

\* Only available in the U.S.

 <b>1 Standard Finish</b> – stainless steel chestpiece and eartubes.	 <b>5 Smoke Edition</b> – smoke finish chestpiece and eartubes.	 <b>9 Black Edition</b> – black finish chestpiece, stem and eartubes.	 <b>13 Mirror Edition</b> – mirror finish chestpiece and smoke finish stem and eartubes.
 <b>2 Black Edition</b> – black finish chestpiece and eartubes.	 <b>6 Smoke Edition</b> – smoke finish chestpiece and black finish eartubes.	 <b>10 Smoke Edition</b> – smoke finish chestpiece, stem and eartubes.	 <b>14 Mirror Edition</b> – mirror finish chestpiece and stem and stainless steel finish eartubes.
 <b>3 Brass Edition</b> – brass finish chestpiece and eartubes.	 <b>7 Rainbow Edition</b> – rainbow finish chestpiece and brass finish eartubes.	 <b>11 Brass Edition</b> – brass finish chestpiece and smoke finish stem and eartubes.	 <b>15 Champagne Edition</b> – champagne finish chestpiece and smoke finish stem and eartubes.
 <b>4 Copper Edition</b> – copper finish chestpiece and eartubes.	 <b>8 Rainbow Edition</b> – rainbow finish chestpiece and black finish eartubes.	 <b>12 Rainbow Edition</b> – rainbow finish chestpiece and black finish stem and eartubes.	

*Due to limitations of the printing process, the colors shown above are approximations of the actual tubing colors.*

**3M**  
**Littmann®**  
Stethoscopes

3M Health Care  
3M Infection Prevention Division  
2510 Conway Ave., St. Paul, MN 55144 USA  
1-800-228-3957 | www.littmann.com

3M, Littmann, the Littmann Quality seal, the L, and Classic III are trademarks of 3M, used under license in Canada.

Please recycle. Printed in USA. © 3M 2018.  
All rights reserved. Made in USA. Issued: 3/18 14014HB  
70-2011-5600-0