



Littmann®
Stethoscopes

Classic

Love your stethoscope.

The compact and sensitive **3M™ Littmann® Classic III™ Monitoring Stethoscope** is the perfect instrument for clinicians who monitor and assess children or adults in non-critical care environments. Manufactured with strong yet lightweight materials for hours of comfortable use.



Monitor and assess a wide range of patients.

- Useful in non-critical care environments such as a medical office, general ward, OB/GYN, ambulatory clinic or urgent care
- Fun and vibrant; match color and finish to your personality
- Designed for use with adult and pediatric patients
- Tunable, dual-sided stainless steel chestpiece with open or closed bell
- Non-chill bell sleeve for greater patient comfort
- Soft-sealing eartips provide optimal comfort and excellent sound occlusion
- Made in the USA without natural rubber latex

Product Features

Acoustic Performance	8 (1–10 scale)	
Chestpiece Design/Technology	Machined stainless steel, two-sided	
Headset Material	Lightweight, aircraft-grade, aluminum alloy	
Approximate Length	27"/69 cm	
Approximate Weight	150 grams	
Diaphragm Details/Diameter	Two tunable diaphragms 1.7"/4.3 cm, 1.4"/3.4cm	
Warranty	5 years	