



















	Procedure	Typical Incision Size	Overall Island Dressing Size*	Medipore™ +Pad Soft Cloth Dressing	Tegaderm™ +Pad Film Dressing
CARDIOVASCULAR/ THORACIC	 Coronary Artery Bypass Graft (CABG)	Traditional sternal: Avg. 6–8" (15–20cm) up to 12" (31 cm)	3½" x 10"	3571	3591
		Minimally invasive: 3–4" (8–10 cm)	3½" x 13¾"	3573	3593
			3½" x 6"	3569	3589
		Vein harvesting: Several 1–½" (4 cm) incisions	3½" x 8"	3570	3590
 Peripheral Stent	Very small (approx. 1–3cm), sometimes in groin	2" x 2¾"	3562	3582	
		2¾" x 4"	3564	3584	
GENERAL SURGERY	 Appendectomy	Lap: 1" (2.5 cm)	2" x 2¾"	3562	3582
			2¾" x 4"	3564	3584
		Open incision: 3–6" (8–15 cm)	3½" x 4"	3566	3587 (oval)
			3½" x 6"	3569	3589
	 Breast Biopsy	2" (5 cm)	2¾" x 4"	3564	3584
			3½" x 4"	3566	3586 3587 (oval)
	 Cholecystectomy	Lap: 4 to 5 incisions ½" (1 cm) in size	2" x 2¾"	3562	3582
			2¾" x 4"	3564	3584
		Open incision: 3–6" (8–15 cm)	2¾" x 4"	3564	3584
			3½" x 4"	3566	3586
	 Gastric Bypass	Lap procedure: Four 1" (2.5 cm)	2" x 2¾"	3562	3582
			2¾" x 4"	3564	3584
		Midline: 6–10" (15–25 cm)	3½" x 4"	3566	3586
			3½" x 8"	3570	3590
 Inguinal Hernia	Lap: 1" (2.5 cm)	2" x 2¾"	3562	3582	
		2¾" x 4"	3564	3584	
	Open: up to 4" (10 cm)	3½" x 4"	3566	3586	
		3½" x 6"	3569	3589	
 Mastectomy	Avg. 6" (15 cm), depending on the procedure	3½" x 8"	3570	3590	
		3½" x 10"	3571	3591	
 Partial Colectomy	4" (10 cm) transverse abdominal or 7–8" (18–20 cm) midline	6" x 6"	3568	3588	
		3½" x 6"	3569	3589	
		3½" x 8"	3570	3590	
 Cesarean Section	6–8" (15–20 cm) bikini cut, pubic; lower abdominal	3½" x 10"	3571	3591	
		3½" x 13¾"	3573	3593	
OB-GYN	 Hysterectomy	Abdominal: 5–7" (13–18 cm)	3½" x 8"	3570	3590
			3½" x 10"	3571	3591
		Laparoscopic: Three ½" (1 cm) LAVH: Two ½" (1 cm)	2" x 2¾"	3562	3582
			2¾" x 4"	3564	3584
			3½" x 4"	3566	3587 (oval)
			3½" x 6"	3569	3589
 Carpal Tunnel	1–3" (2.5–8 cm)	2¾" x 4"	3564	3584	
		3½" x 4"	3566	3587 (oval)	
		3½" x 6"	3569	3589	
ORTHOPEDIC	 Fracture Fixation	Varies, depending on location of fracture. Often approx. 4" (10 cm) long. May have 2 incisions.	3½" x 6"	3569	3589
			3½" x 8"	3570	3590
	 Hip Replacement	8" (20 cm) over incision; some 8–12" (20–31 cm), or as small as 4–6" (10–15 cm)	3½" x 8"	3570	3590
			3½" x 10"	3571	3591
			3½" x 13¾"	3573	3593
	 Knee Replacement	4–6" (10–15 cm) total knee vertical incision, arthroplasty	3½" x 6"	3569	3589
3½" x 8"			3570	3590	
3½" x 10"			3571	3591	
 Spine	Laminectomy spinal infusion 2–6" (5–15 cm)	2¾" x 4"	3564	3584	
		3½" x 4"	3566	3586	
		3½" x 6"	3569	3589	
		3½" x 8"	3570	3590	
PLASTIC SURGERY	 Breast Augmentation or Reduction	3" (8cm) vertical and/or 3" 'smile' under breast reduction	2¾" x 4"	3564	3584
			3½" x 4"	3566	3586
			3½" x 6"	3569	3589
	 Liposuction	Laparoscopic procedure: Can be up to 12 incisions, ¼"–½" (½–1 cm)	2" x 2¾"	3562	3582
		2¾" x 4"	3564	3584	

dress for SUCCESS

Suggested Island Dressings for Surgical Procedures

3M™ Tegaderm™ and 3M™ Medipore™ +Pad dressings are all-in-one sterile dressings designed:

- for ease of use
- to reduce application time
- to conform to body contours, providing comfort
- with breathable film to allow skin to function normally with good exchange of moisture vapor and oxygen
- to conform to body contours, and to flex with movement or swelling

Both feature non-adherent, absorbent pads:

- to promote normal healing
- with less pain and trauma
- for light to moderately draining wounds
- won't stick to the wound

3M™ Tegaderm™ +Pad Film Dressing with Non-Adherent Pad

Film Barrier – provides both a viral** and a bacterial barrier. Requires fewer dressing changes, while helping to reduce the risk of surgical site infections

Transparent Border – allows monitoring of peri-wound skin without removing dressing

Absorbent, Non-Adherent Pad – for light to moderately draining wounds; won't stick to the wound, allowing normal wound healing process

Waterproof – patient may sponge bathe or shower without removing dressing

Versatile – available in 9 different sizes

3M™ Medipore™ +Pad Soft Cloth Adhesive Wound Dressing

Gentle to skin – provides excellent adhesion, yet is gentle to the skin

All-in-one design – reduces application time

Absorbent, Non-Adherent Pad – for light to moderately draining wounds; won't stick to the wound, allowing normal wound healing process

Breathable – better exchange of moisture vapor and oxygen, and reduced risk of skin maceration

Versatile – available in 8 different sizes

*Choose a dressing size large enough to provide a margin that adheres to dry healthy skin around the incision site.

**In vitro testing shows that the transparent film dressing provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.