

RelyX™ Cements



RelyX™ Luting Plus Automix Tips

Item No.	Product Information
3540	RelyX™ Luting Plus Mixing Tips. Includes: 12 Regular Mixing Tips Pink

RelyX™ Unicem 2 Automix Tips

Item No.	Product Information
56918	RelyX™ Unicem 2 Automix Tips. Includes: 15 Mixing Tips Wide Yellow and 15 Intraoral Tips.
56919	RelyX™ Unicem 2 Automix Tips. Includes: 15 Mixing Tips Wide Yellow and 15 Endo Tips.
56851	RelyX™ Unicem 2 Automix Mixing Tips. Includes: 30 Regular Mixing Tips Yellow.

RelyX™ Ultimate Tips

Item No.	Product Information
56883	RelyX™ Ultimate Adhesive Resin Cement Mixing Tips. Includes: 30 Regular Mixing Tips Green
56885	RelyX™ Ultimate Adhesive Resin Cement Tips. Includes: 10 Wide Mixing Tips Green plus 10 Endo Tips.
56886	RelyX™ Ultimate Adhesive Resin Cement Tips. Includes: 10 Wide Mixing Tips Green plus 10 Intraoral Tips.

3M Australia Pty Ltd

Building A, 1 Rivett Road
North Ryde NSW 2113
Ph: 1300 363 454
www.3MESPE.com.au

3M New Zealand Ltd

94 Apollo Drive
Rosedale Auckland 0632
Ph: 0800 80 81 82
www.3MESPE.co.nz



3M, ESPE, Lava, RelyX and Scotchbond are trademarks of 3M or 3M ESPE.

All other trademarks are owned by other companies.

© 2015, 3M. All rights reserved.

RelyX™ Cements

For each indication –
the ideal cement



	RelyX™ Ultimate Adhesive Resin Cement	RelyX™ Unicem RelyX™ Unicem 2 Self-Adhesive Resin Cement	RelyX™ Luting Plus Resin Modified Glass Ionomer Cement	RelyX™ Veneer Veneer Cement
Metal / Metal based				
Inlays / Onlays	+	++	+	-
Crowns / Bridges	+	++	++	-
Endodontic Posts	+	++	++	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Glass Ceramics (incl. Li. Disilicate) (e.g. e.max, VITA Mark II, IPS Empress 2)				
Inlays / Onlays / Table Tops	++	+	-	-
Crowns / Bridges	++	++	-	-
Veneers	+	-	-	++
On Implant Abutments	+	+	-	-
Oxide Ceramics (e.g. Lava™ Plus, Brux Zir, Procera)				
Inlays / Onlays	++	++	+	-
Crowns / Bridges	+	++	+	-
Endodontic Posts	+	++	+	-
Maryland Bridges	++	+	-	-
On Implant Abutments	+	++	++	-
Resin Nano Ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)				
Inlays / Onlays	++	-	-	-
Veneers	+	-	-	++
Resin Composites				
Inlays / Onlays	++	+	-	-
Crowns	++	++	-	-
Endodontic Posts	+	++	-	-
Veneers	+	-	-	++
On Implant Abutments	++	++	-	-

++ Highly recommended*

+ Recommended

- Not indicated

*Either better performance for this indication or easier handling with equal performance.

3M ESPE