


**Forsus™**  
**Fatigue Resistant Device**  
**with Direct Push Rod**




3M, an industry-leading innovator, is pushing the limitations of Class II Correction and patient satisfaction. Until now, your options were limited to fixed intraoral devices with breakage concerns, or other appliances, including headgear, requiring a high degree of patient compliance. That's why we created a more reliable appliance for Class II Correction—a 3-piece telescoping spring solution that fits perfectly with your current appliance system and is more attractive and easier for your patients to live with. The design incorporates an external spring assembly that won't compromise patient comfort and acceptance, a key to achieving successful clinical outcomes when using compliance-independent appliances. ***With the addition of the Direct Push Rod, installation can now be achieved in minutes. The Direct Push Rod attaches directly onto the archwire, eliminating the need for auxiliaries and simplifying chairside installation.***

**Fatigue resistant by design.**

The unique coaxial spring design of the Forsus™ Fatigue Resistant Device addresses the issue of fatigue failure—a fracture caused by repeated application of stresses in the coil spring. To maximize fatigue life, the Forsus Fatigue Resistant Device is designed to keep these stresses to a minimum. And its superior strength and consistent forces have been shown in laboratory data.

Traditional solutions simply can't compare in terms of strength, reliability, ease of use and patient satisfaction. Put the Forsus Fatigue Resistant Device to work for you, and discover how it can force your productivity to new heights.



## Forsus™ Fatigue Resistant Device with Direct Push Rod


Description	Unit	Part No.
Forsus Fatigue Resistant Device 5-Patient Kit	Kit	885-121
Forsus Fatigue Resistant Device 20-Patient Kit	Kit	885-122
Universal Spring Module	pkg/5	885-100
Direct Push Rod Short (25mm)	pkg/5 L	885-111
Direct Push Rod Short (25mm)	pkg/5 R	885-112
Direct Push Rod Med (29mm)	pkg/5 L	885-113
Direct Push Rod Med (29mm)	pkg/5 R	885-114
Direct Push Rod Lg (32mm)	pkg/5 L	885-115
Direct Push Rod Lg (32mm)	pkg/5 R	885-116
Direct Push Rod XL (35mm)	pkg/5 L	885-117
Direct Push Rod XL (35mm)	pkg/5 R	885-118
Direct Push Rod XXL (38mm) no stop	pkg/5 L	885-107
Direct Push Rod XXL (38mm) no stop	pkg/5 R	885-108
Universal L Pin	pkg/10	885-110
Universal Split Crimp	pkg/50	885-203
Measurement Guide	pkg/5	807-014

**Pushing the limitations  
of Class II Correctors.**



*Easy to install  
L-Pin Connector*

*A 3-piece telescoping coaxial  
spring designed to allow for  
enhanced patient jaw opening*



*Direct Push Rod attaches directly  
onto the archwire*



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA  
Phone 1-800-423-4588  
Web 3M.com/ortho

3M Canada  
Health Care Division  
300 Tartan Dr.  
London, ON N5V 4M9 Canada  
Phone 1-800-443-1661

3M and Forsus are trademarks of 3M.  
Used under license in Canada.  
© 3M 2016. All rights reserved.  
70-2021-3922-9 1601