

Dear Doctor _____ ,

today you are receiving this high quality restoration:

- Metal/Metal based**
- Glass Ceramics incl. Li. Disilicate**
(e.g. IPS e.max®, VITABLOCS® Mark II, IPS Empress®)
- Zirconia**
(e.g. Lava™ Plus, BruxZir®, Procera®)
- Composites**
- Resin Nano Ceramics**
(e.g. Lava™ Ultimate CAD/CAM Restorative)
- _____

We recommend using the following cement:

- Adhesive Resin Cement**
(e.g. RelyX™ Ultimate with Scotchbond™ Universal Adhesive or: RelyX™ Veneer)
- Self-adhesive Resin Cement**
(e.g. RelyX™ Unicem, RelyX™ Unicem 2)
- Conventional Luting Cement**
(e.g. RelyX™ Luting Plus, Ketac™ Cem)



Pretreatment recommendations:

	HF etch	Silane treatment (incl. Rocatec™ / CoJet™)	Sand blast ≤ 50 µm
Done by lab			
Still needed			
<input type="checkbox"/> No pretreatment required			

Best regards,

Lab stamp / date / signature

3M, ESPE, CoJet, Ketac, Lava, RelyX, Rocatec and Scotchbond are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by other companies. © 2014, 3M. All rights reserved. 70-2013-0685-2