

3M[™] PICC/CVC Securement Device + Tegaderm[™] I.V. Advanced Dressing

With 3M's new PICC/CVC securement system, you don't have to sacrifice patient comfort to get reliable securement. The Tegaderm[®] I.V. Advanced Dressing was specifically designed to enhance securement of the device, helping to ensure the best possible experience and outcomes. The system is designed to stay in place for up to 7 days, yet removes gently, delivering the securement and comfort your patients deserve.



3M[™] PICC/CVC Securement

Device + Tegaderm[™] I.V. Advanced Dressing

Securement without Sacrifice


- Sutureless securement balances reliable adhesion with gentleness to skin
- Designed to be worn for up to 7 days
- Removes cleanly, without causing patients undue pain or distress
- Skin-friendly materials help improve comfort and mobility

Helps make your life easier

- · Quick and easy to apply and remove
- Accommodates the majority of PICC and short-term CVC catheters up to and including 12 French with up to 3 lumens
- Can be repositioned upon initial application without losing adhesion
- Single package with device and dressing simplifies product selection
- Supports CDC guidelines for sutureless securement and INS standards recommending the use of a catheter stabilization device

Help Reduce Overall Cost of Care

• Eliminates the risks and costs of suture-related needlestick injuries





Silicone adhesive

- Holds securely, yet removes gently
- Can be repositioned upon initial application
- Does not require alcohol to remove

Easy application securement base

- No mechanical doors or wings
- Visual cues to aid application

Comfortable, soft-cloth material

- Conformable and water-resistant
- Perforated for additional breathability

Documentation tape strip

- Preprinted for documenting dressing changes
 Provides additional securement
- Provides additional securement
 Covered with water-resistant transparent film
- Covered with water-resistant transparent in

, Tegaderm[™] brand transparent film

- Provides continuous site observation
- Is a waterproof, sterile barrier to external contaminants*

Reinforced stabilization border

- Maximizes securement, breathability and wear time
- Dressing and tape strips are water-resistant

Securement tape strip

- Perforated slit in large tape strip to ease removal
- Adhesive-free frame minimizes potential to stick
 to gloves or itself
- Allows for one-handed delivery
- Allows for one-nanded delivery
 Covered with water-resistant transparent film

or larger while the dressing remains intact without leakage

* In vitro testing shows that the transparent film of 3M[°] Tegaderm[°] I.V. Advanced Dressing provides a viral barrier for viruses 27 mm in diameter

Ordering Information

Product Code	Material ID	Overall Device Size	Overall Dressing Size	Device + Dressing/Box	Boxes/Case
1837-2100	7100033117	2 in x 2 1/8 in 5.1 cm x 5.4 cm	3 1/2 in x 4 1/2 in 8.5 cm x 11.5 cm	20	4
1839-2100	7100054100	2 in x 2 1/8 in 5.1 cm x 5.4 cm	4 in x 6 1/8 in 10 cm x 15.5 cm	20	4

Learn More

To learn about 3M[®] PICC/CVC Securement Device + Tegaderm[®] I.V. Advanced Dressing visit us at go.3M.com/3MSecurement, contact your 3M Critical & Chronic Care Solutions representative or call the 3M Health Care Customer Helpline at 1-800-364-3577.

3M Critical & Chronic Care Solutions Division collaborates to deliver proven, innovative solutions to help prevent and treat critical and chronic conditions across the continuum of care.



Available in Canada from:

Critical & Chronic Care Solutions Division 3M Canada Company P.O. Box 5757 London, Ontario N6A 4T1 Canada 1-800-563-2921 www.3M.ca/skinwound

3M Health Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA 1-800-228-3957 www.3M.com/C3SD

Please recycle. Printed in Canada. © 2015, 3M. All rights reserved. 3M and Tegaderm are trademarks of 3M. Used under license in Canada. 1502-00339 E