

# 3M™ Tegaderm™ Roll





## Transparent Film Roll

A non-sterile version of Tegaderm™ Transparent Film



*The roll version allows for maximum versatility as it can easily be cut to size.*

### Application Technique (Tegaderm™ Roll over primary dressing)



### Primary Uses

- Protection of at risk skin from friction and body fluids
- Waterproof fixation of primary dressings such as gauzes, alginates, foams, etc.
- Fixation of medical devices (eg. drainage tubes, etc.)
- Provides a bacterial and viral barrier to protect patients and health care workers

Roll Style	Cat. No.	Size	Roll/Box	Boxes/Case
	16002	2 in x 11 yd 5 cm x 10 m	1	4
	16004	4 in x 11 yd 10 cm x 10 m	1	4
	16006	6 in x 11 yd 15 cm x 10 m	1	4

**Tegaderm™**  
Dressings

SIMPLE. DEPENDABLE. TRUSTED.

**3M Canada Company**  
P.O. Box 5757  
London, Ontario N6A 4T1  
1 800 563-2921 (Ontario)  
1 800 361-4488 (Quebec and Atlantic  
Canada)  
1 800 661-8126 (Western Canada)  
[www.3M.com/ca/healthcare](http://www.3M.com/ca/healthcare)

**3M Health Care**  
3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000

3M and Tegaderm are trademarks of 3M.  
Used under license in Canada.  
©3M 2005 0509-MS-23348-E

For more information visit us  
at contact your 3M Health Care  
Sales Representative, or call  
the 3M Health Care Customer  
Helpline at 1 800 364-3577.

**3M** Health Care