

3M Health Care

3M™ Tegaderm™ Absorbent Clear Acrylic Dressing

# Tegaderm™ Absorbent Clear

## Acrylic Dressings

a clear advance  
in wound management

### Allows wound monitoring without changing the dressing

- Versatile – ideal for skin tears, pressure ulcers, abrasions, post-operative wounds and donor sites
- Provides a moist wound environment, which promotes autolytic wound debridement
- Very conformable and comfortable for patients
- Easy to apply and remove

### Rated significantly better\* than the leading hydrocolloid dressing on:


- Ability to assess the wound before and after absorption
- Barrier properties
- Patient comfort
- Conformability before and after absorption
- Ease of removal
- Non-adherence to wound bed
- Residue in wound bed and skin
- Odour
- Wear time

\*Multi-site, randomized, controlled study data on file at 3M.

**3M**

# The benefits of 3M™ Tegaderm™ Transparent Dressings with the added advantage of absorbency

## Tegaderm™ Absorbent Clear Acrylic Dressing - Advanced 3M technology



### Easy Wound Monitoring

- Monitor wounds without changing the dressing
- May eliminate unnecessary dressing changes
- May save supply costs and nursing time

### Comfortable for Patients

- Gentle removal from skin
- Very conformable
- Long wear time

### Industry Best Practice

- Recommended industry best practice in the treatment of skin tears\*

### Indications for Use

Tegaderm™ Absorbent Clear Acrylic Dressings are ideal for:

- Pressure ulcers
- Skin tears
- Abrasions
- Donor sites
- Superficial burns
- Post-operative Wounds

### Convenient to Use

- Large, non-adhesive tabs for easy handling with gloves
- Easy to centre on the wound and apply
- User-friendly packaging


### Safe and Effective

- Tested for safety in accordance with ISO 10993 Part 1
- Low-friction surface minimizes potential for friction and shear
- Low potential for periwound maceration
- Adherent film border minimizes likelihood of dressing failure
- Hypoallergenic and latex-free

This dressing can also be used as a secondary (cover) dressing over wound fillers (such as alginate dressings).

\* Leblanc Kimberly, et al, "Best Practice Recommendations for the Prevention and Treatment of Skin Tears," Wound Care Canada Volume 6, Number 1 (May 2008): p.14-30

# Where to use 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing



As a primary dressing over Stage II/Stage III minimal to moderately draining pressure ulcers



Over skin tears and superficial traumatic wounds



As a secondary dressing cover over wounds requiring an absorbent wound filler









As a primary dressing over skin graft donor sites



Over pressure-related injuries

# 3M™ Tegaderm™ Absorbent Clear Acrylic Dressing

## Ordering Information

Catalogue No.	Dressing Pad Size	Overall Dressing Size	Dressings/Box	Boxes/Case
 90800 Oval	1½ in. x 2½ in. 3.8 x 5.7 cm	3 in. x 3½ in. 7.6 x 9.5 cm	5	6
 90801 Oval	2½ in. x 3 in. 6.0 x 7.6 cm	4½ in. x 5 in. 11.1 x 12.7 cm	5	6
 90802 Square	3½ in. x 4 in. 9.8 x 10.1 cm	5½ in. x 6 in. 14.9 x 15.2 cm	5	6
 90803 Oval	3½ in. x 4½ in. 8.5 x 10.7 cm	5½ in. x 6½ in. 14.2 x 15.8 cm	5	6
 90805 Square	5½ in. x 6 in. 14.9 x 15.2 cm	7½ in. x 8 in. 20.0 x 20.3 cm	5	4
 90807 Sacral	4½ in. x 5 in. 11.75 x 14.3 cm	6½ in. x 7½ in. 16.8 x 19 cm	5	4



### A clear solution from the makers of 3M™ Tegaderm™ Transparent Dressings

Tegaderm™ Absorbent Clear Acrylic dressing, a transparent, absorbent dressing for optimal clinical performance and wear time.

### Compatible 3M Skin Health Products

- 3M™ Cavilon™ No Sting Barrier Film
- 3M™ Tegaderm™ High Gelling Alginate Dressing
- 3M™ Tegaderm™ High Integrity Alginate Dressing
- 3M™ Tegaderm™ Non-Adherent Contact Layer

**3M™ Tegaderm™ Dressings** SIMPLE. DEPENDABLE. TRUSTED.

For more information, contact your 3M Skin Health representative or call the 3M Health Care Customer Helpline at 1 800 364-3577. These products can be ordered from your local distributor.



**3M Canada**  
P.O. Box 5757  
London, ON N6A 4T1  
Canada

1 800 364-3577

[www.3M.com/ca/healthcare](http://www.3M.com/ca/healthcare)

**3M Health Care**  
3M Centre, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA

3M, Cavilon and Tegaderm are trademarks of 3M.  
Used under license in Canada.  
Please recycle. Printed in Canada.  
© 2008, 3M. All rights reserved.  
0806-2466