


3M™ Tegaderm™ Alginate Ag Silver Dressing

# Tegaderm™ Alginate Ag

Silver Dressing



Delivering the  
Power of Silver

A new dressing in the trusted  
Tegaderm™ brand family of products


**3M**



# Versatile, powerful, easy to use

## 3M™ Tegaderm™ Alginate Ag Silver Dressing


New 3M™ Tegaderm™ Alginate Ag Silver Dressing provides the antimicrobial silver, soft-gelling absorbency and fibre strength that professionals need to manage critically colonized wounds and wounds at increased risk of infection. In addition to the dressing's proven silver compound, its carboxymethylcellulose (CMC) and alginate fibre formula creates a dressing that conforms to the wound surface and has high wet-strength for easy removal. Available in both pads and ropes, these soft dressings are easy to use and compliment our broad range of trusted and reliable Tegaderm™ brand dressings.



## Antimicrobial effectiveness

A highly soluble form of silver that:

- Reduces a broad spectrum of pathogens, including MRSA and VRE
- Helps control wound odour
- Provides a sustained antimicrobial barrier for up to 14 days



## Soft-gelling absorbency

The unique combination of CMC and alginate creates a soft-gelling formula that:

- Conforms to the wound and is soft and pliable
- Assists in controlling minor bleeding
- Provides excellent exudate management with greater absorbency than Aquacel® Ag Hydrofiber® Wound Dressing with Silver




## One-piece removal

The CMC and alginate formula delivers high wet-strength, which:

- Allows for use in wounds with undermining or tunnels
- Facilitates easy removal of moist dressing
- Ensures a moist environment for optimal wound healing

### Absorbency of dressings\*






Statistical significant differences (P < 0.001)

\*Absorbency of "flat" dressings only. Data on file at 3M.

# 3M™ Tegaderm™ Alginate Ag Silver Dressing

## Ordering Information

	Catalogue number	Description	Dressings/ Box	Boxes/ Case
	90301	Alginate Ag Dressing 5.08 cm x 5.08 cm (2" x 2")	10	4
	90303	Alginate Ag Dressing 10.16 cm x 12.7 cm (4" x 5")	10	4
	90307	Alginate Ag Rope 2.54 cm x 30.48 cm (1" x 12")	5	4

## Compatible with:

3M™ Cavilon™ No Sting Barrier Film

3M™ Tegaderm™ Absorbent Clear Acrylic Dressings

3M™ Tegaderm™ Contact Layer

3M™ Tegaderm™ Foam Dressings

3M™ Tegaderm™ Hydrocolloid Dressings

3M™ Coban™ 2 Layer Compression System

**3M™ Tegaderm™** Dressings SIMPLE. DEPENDABLE. TRUSTED.

To learn more about Tegaderm™ Alginate Ag silver dressing, visit us at [www.3m.com/Tegaderm-AlginateAg](http://www.3m.com/Tegaderm-AlginateAg). For more information about the Skin & Wound Care family of products, visit us at [www.3M.com/skinhealth](http://www.3M.com/skinhealth), contact your 3M Skin & Wound Care representative or call the 3M Health Care Customer Helpline at 1-800-364-3577.



**Skin & Wound Care Division**  
**3M Canada Company**  
 P.O. Box 5757  
 London, Ontario N6A 4T1  
 Canada  
 1 800-364-3577  
[www.3m.com/HealthCare](http://www.3m.com/HealthCare)

**3M Health Care**  
 3M Center Building 275-4W-02  
 St. Paul, MN 55144-1000  
 U.S.A.  
 1 800 228-3957

Please recycle. Printed in Canada.  
 © 2009, 3M. All rights reserved.  
 3M, Tegaderm, Cavilon and Coban are trademarks of 3M.  
 Aquacel Ag is a trademark of Convatec. Acticoat is a trademark of Smith & Nephew.  
 Used under license in Canada.  
 0910-07478E