

For a Life in Motion



Conformable Compression that Works

Ideal Compression Therapy for Clients of all Sizes, Shapes and Lifestyles

Three convenient kits, making compression therapy easier for both the clinician and the client

3M™ Coban™ 2 Layer Compression System has been clinically proven in two randomized controlled trials^{1,2}.

- Choose moderate or high compression based on client needs
- Provide sustained, therapeutic compression for up to 7 days.²
- Application is quick and easy to teach and learn.⁴
- Creates a thin, lightweight, breathable sleeve, allowing clients to wear their own shoes and return to their daily activities

Therapeutic Compression Designed for Everyday Comfort



3M™ Coban™ 2 Layer Compression System

3M™ Coban™ 2 Layer Compression System, with the purple packaging, provides high compression (35-40 mm Hg). It is ideal for the majority of clients with venous leg ulcers, lymphedema and other conditions where compression therapy is appropriate.

Indicated for clients with ABPI ≥ 0.8 .

The 10 cm (4") kit is most commonly used to bandage the leg from the foot to below the knee. The 15 cm (6") kit is most commonly used to bandage the knee and the thigh. It can also be used below the knee on a larger limb.

2094 4 in x 2.5 yd (10 cm x 2,7 m)

20096 6 in x 3.8 yd (15 cm x 3,5 m)

3M™ Coban™ 2 Layer Lite Compression System

3M™ Coban™ 2 Layer Lite Compression System, with the green packaging, provides moderate compression (30 -35 mm Hg). It is designed for clients less tolerant of compression therapy³, including those who have mixed etiology with ABPI ≥ 0.5 , are new to compression or where tolerance is not known, are frail, or are less mobile.

2794 4 in x 2.5 yd (10 cm x 2,7 m)





To view application **videos** and download
reference material, visit www.3m.ca/Coban2Layer.

Ordering Information: 3M™ Coban™ 2 Layer Compression Therapy			
	3M™ Coban™ 2 Layer Lite Compression System	3M™ Coban™ 2 Layer Compression System	
Catalogue Number	2794	2094	20096
Dimensions	4 in x 2.5 yd (10 cm x 2,7 m)	4 in x 2.5 yd (10 cm x 2,7 m)	6 in x 3.8 yd (15 cm x 3,5 m)
Rolls/ box	2 rolls, 1 of each layer	2 rolls, 1 of each layer	2 rolls, 1 of each layer
Boxes/ case	8	8	8

References

1. Moffatt C, Edwards L, Collier M, Treadwell T, Miller M, Shafer L, Sibbald RG, Brassard A, McIntosh A, Reyzelman A, Price P, Kraus SM, Walters SA, Harding K. A randomized controlled 8-week crossover clinical evaluation of the 3M™ Coban™ 2 Layer Compression System versus Profore™ to evaluate the product performance in patients with venous leg ulcers. *Int Wound J* 2008; 5 (2): 267-279.
2. Mosti G, Crespi A, Mattaliano V. Comparison Between a New, Two-component Compression System With Zinc Paste Bandages for Leg Ulcer Healing: A Prospective, Multicenter, Randomized, Controlled Trial Monitoring Sub-bandage Pressures. *Wounds* 2011; 23 (5): 126-134.
3. Junger, M, Haase H, Ladwig A, Schwenke L, Bichel J, Schuren J. Compression therapy in patients with peripheral arterial occlusive disease: A prospective clinical study with the 3M™ Coban™ 2 Layer Lite Compression System for ABPI \geq 0.5. Data on file at 3M. 2010.
4. Collier M, Schuren J. Ease of use and reproducibility of five compression systems. *J Wound Care* 2007; 3M Supplement: 8-10.

To learn more, visit us at www.3M.ca/coban2layer, contact your 3M Critical & Chronic Care Solutions representative or call the 3M Health Care Customer Helpline at **1-800-228-3957**.



Available in Canada from:

**3M Critical & Chronic Care
Solutions Division 3M Health Care**
P.O. Box 5757
London, Ontario N6A 4T1
Canada
1-800-364-3577
www.3m.ca/coban2layer

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str.1
41453 Neuss Germany

Profore is a trademark of Smith & Nephew.
3M and Coban are trademarks of 3M.
Used under license in Canada.
Please recycle. Printed in Canada.
© 2015, 3M. All rights reserved.
1506-01601 E