

Efficient Posterior Restorative Solution

**Takes less time,
so you have
more time.**



Simplify your day and reduce chair time with faster, easier posterior restorations.

We all know that the more complicated a procedure is, the greater chance there is for something to go wrong and disrupt your day. That's why 3M simplified posterior restorations, one of the most common and complicated procedures.

By using four of our innovative technologies together, you can move through a posterior restoration with speed and simplicity.

Essential to this simplified procedure is Filtek™ Bulk Fill Posterior Restorative, designed to improve productivity by enabling fast and easy posterior restorations. You can place it in one convenient increment up to 5 mm—no need for layering or additional dispensing devices that can slow you down.

Let our simple Posterior Restorative Procedure help you stay on schedule ... because your time outside of work matters.

Efficient and consistent procedure for Class II restorations

- ▶ Simple, one-step placement up to 5 mm with Filtek™ Bulk Fill Posterior Restorative
- ▶ High bond strength and virtually no post-op sensitivity with Scotchbond™ Universal Adhesive
- ▶ Deep, uniform cure with Elipar™ DeepCure-S LED Curing Light
- ▶ Quickly polish restorations to a high-gloss finish with Sof-Lex™ Diamond Polishing System



Fast and easy

Designed to improve productivity from start to finish.




Scotchbond™ Universal Adhesive

Versatile. High bond strength. One simple technique and quick application. Virtually no post-operative sensitivity.



Filtek™ Bulk Fill Posterior Restorative

Simplicity of bulk placement up to 5 mm deep. Low wear for durability and stress relief.



Elipar™ DeepCure-S LED Curing Light

Predictable, reliable cure. Optimized light guide is more comfortable for patients and provides easier access to posterior teeth.



Sof-Lex™ Diamond Polishing System

Easily adapts to irregular, convex and concave tooth surfaces and imparts paste-like gloss with the convenience of a rubberized system. Fast and easy to use.



Ordering Information



Filtek™ Bulk Fill Posterior Restorative

Refill SHADE	Item #	Item #
A1	4863A1	4864A1
A2	4863A2	4864A2
A3	4863A3	4864A3
B1	4863B1	4864B1
C2	4863C2	4864C2

SYRINGE SHADE REFILLS include 1 – 4g syringe; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS include 20 – 0.2g capsules; Technique Guide; Instructions for Use.



Scotchbond™ Universal Adhesive

Item #	Product Information
41254	Scotchbond™ Universal Adhesive Intro Kit Kit includes: 1 – 5ml adhesive vial; 50 disposable mixing wells; 50 disposable microbrushes; 1 – 3ml etchant syringe
41258	Scotchbond™ Universal Adhesive Refill: 1 – 5ml adhesive vial
41255	Scotchbond™ Universal Adhesive Intro Kit: Unit Dose 40: Includes 40 unit dose; 1 – 3ml etchant syringe
41256	Scotchbond™ Universal Adhesive Intro Kit: Unit Dose 100: Includes 100 unit dose; 1 – 3ml etchant syringe
41257	Scotchbond™ Universal Adhesive Intro Kit: Unit Dose 200: Includes 200 unit dose; 1 – 3ml etchant syringe
41263	Scotchbond™ Universal Etchant Etching Gel: 2 – 3ml etchant syringes; 50 dispensing tips
41264	Scotchbond™ Universal DCA Dual Cure Activator: 1 – 5ml vial



Elipar™ DeepCure-S LED Curing Light

Item #	Product Information
76976	Elipar™ DeepCure-S LED Curing Light Contains: Handpiece (Cordless), Charging Base (230V); Li-ion Battery; 10mm Light Guide; Eye Shield
76981	Elipar™ DeepCure-S Light Guide, 10mm
76984	Elipar™ DeepCure Eye Shield
76985	Elipar™ DeepCure-S Rechargeable Li-ion Battery



Sof-Lex™ Diamond Polishing System

Item #	Product Information
5092-IM	Sof-Lex™ Diamond Polishing System—Introductory Kit Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); 1 Sof-Lex™ RA Mandrel; Instructions for Use; Technique Card
5090	Sof-Lex™ Pre-Polishing Spirals—Refill Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
5091	Sof-Lex™ Diamond Polishing Spirals—Refill Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)

www.3M.com/EfficientRestorations



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, ESPE, Elipar, Filtek, Scotchbond and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA.
© 3M 2016. All rights reserved.
70-2013-0712-4