

Portfolio Overview

Impregum™ Polyether Impression Material

PRODUCT	DISPENSING SYSTEM	VISCOSITY LOW  HIGH	SETTING VERSION	WORKING TIME* AT 23°C/74°F MIN:SEC	TOTAL SETTING TIME** MIN:SEC
TRAY AND MONOPHASE MATERIALS					
Impregum™ Penta™ Soft Heavy Body			Regular Set	2:30	6:00
Impregum™ Penta™ Soft Medium Body			Regular Set	2:45	6:00
Impregum™ Penta™ Medium Body			Regular Set	2:45	6:00
Impregum™ Soft Medium Body (Tray)			Regular Set	1:45	6:00
Impregum™ Penta™ Soft Quick Step Heavy Body			Fast Set	1:00	4:00
Impregum™ Penta™ Soft Quick Step Medium Body			Fast Set	1:00	4:00
Impregum™ Soft Quick Step Medium Body (Tray)			Fast Set	1:00	4:00
WASH MATERIALS					
Impregum™ Penta™ Soft Light Body			Regular Set	3:15	6:30
Impregum™ Soft Light Body			Regular Set	2:00	5:30
Impregum™ Soft Quick Step Light Body			Fast Set	1:00	4:00

*Working time includes mixing time.

**includes working time.

Material Combinations per Technique Impregum™ Polyether Impression Material

TRAY MATERIAL		RECOMMENDED WASH MATERIAL	
MONOPHASE TECHNIQUE – PENTA™			
Impregum™ Penta™ Soft Medium Body			
Impregum™ Penta™ Medium Body			
Impregum™ Penta™ Soft Quick Step Medium Body			
1-STEP TECHNIQUE – PENTA™			
Impregum™ Penta™ Soft Heavy Body		Impregum™ Soft Light Body	
Impregum™ Penta™ Soft Heavy Body		Impregum™ Penta™ Soft Light Body	
Impregum™ Penta™ Soft Quick Step Heavy Body		Impregum™ Soft Quick Step Light Body	
1-STEP TECHNIQUE – CARTRIDGE			
Impregum™ Soft Medium Body (Tray)		Impregum™ Soft Light Body	
Impregum™ Soft Quick Step Medium Body (Tray)		Impregum™ Soft Quick Step Light Body	

Customer Care Center: 1-800-634-2249 www.3MESPE.com



3M ESPE Dental
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1 Canada
1-888-363-3685

3M, ESPE, Impregum and Penta are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
Please recycle. Printed in U.S.A.
© 3M 2015. All rights reserved. 70-2013-0154-9 (Rev. B)