

Portfolio Overview

Imprint™ 4 VPS Impression Material

PRODUCT	DISPENSING SYSTEM	VISCOSITY	SETTING VERSION	MAXIMUM WORKING TIME (23°C/73°F)	MAXIMUM INTRA-ORAL SYRINGING TIME (37°C/98°F)	INTRA-ORAL SETTING TIME (37°C/98°F)
		LOW  HIGH				
TRAY MATERIALS						
Imprint™ 4 Penta™ Putty		Putty	Regular Set	1:30	–	2:30
Imprint™ 4 Penta™ Heavy Imprint™ 4 Heavy		Heavy Body	Regular Set	2:00	–	2:00
Imprint™ 4 Penta™ Super Quick Heavy Imprint™ 4 Super Quick Heavy		Heavy Body	Fast Set	1:15	–	1:15
WASH MATERIALS						
Imprint™ 4 Light			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Light			Fast Set	–	0:35	1:15
Imprint™ 4 Regular			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Regular			Fast Set	–	0:35	1:15

Material Combinations per Technique

Imprint™ 4 VPS Impression Material

TRAY MATERIAL		RECOMMENDED WASH MATERIALS	
1-STEP TECHNIQUE—PENTA™			
Imprint™ 4 Penta™ Heavy Hydrophilic heavy body		Imprint™ 4 Light	
		Imprint™ 4 Regular	
Imprint™ 4 Penta™ Super Quick Heavy Fast setting hydrophilic heavy body		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Super Quick Regular	
Imprint™ 4 Penta™ Putty Putty consistency		Imprint™ 4 Regular	
1-STEP TECHNIQUE—CARTRIDGE			
Imprint™ 4 Super Quick Heavy Fast setting hydrophilic heavy body		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Super Quick Regular	
Imprint™ 4 Heavy Hydrophilic heavy body		Imprint™ 4 Light	
		Imprint™ 4 Regular	
2-STEP TECHNIQUE			
Imprint™ 4 Penta™ Putty Putty consistency		Imprint™ 4 Super Quick Light	
		Imprint™ 4 Light	

Customer Care Center: 1-800-634-2249 www.3MESPE.com/Imprint4



3M ESPE Dental
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1 Canada
1-888-363-3685

3M, ESPE, Imprint and Penta are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
Please recycle. Printed in U.S.A.
© 3M 2015. All rights reserved. 70-2013-0583-9 (Rev. A)