

Procedure Guide

Sof-Lex™ Finishing and Polishing System

A full line of finishing and polishing solutions ... from beginning to end.

The Sof-Lex system is color-coded from dark (coarse) to light (superfine) for an easily followed step-by-step process. Easy, pop-on mounting lets you quickly change discs or spirals—so you can move efficiently through your finishing and polishing procedure.

Gross Reduction

High-Gloss Polishing



Sof-Lex™ Diamond Polishing System:
Forget the messy paste. Our pre-polishing spiral prepares the restoration for final gloss, while our diamond-impregnated polishing spiral gives your restorations that gorgeous paste-like gloss.

Pre-Polishing



PP

Diamond Polishing



HGP

Gross Reduction

High-Gloss Polishing



Sof-Lex™ Discs: Flexible discs conform to contours of teeth

Gross Reduction



C

Contouring



M

Finishing



F

Polishing



SF

Diamond Polishing



HGP



Sof-Lex™ XT Discs: Extra-thin discs for easier access to proximal surfaces

Gross Reduction



C

Contouring



M

Finishing



F

Polishing



SF


Gross Reduction

High-Gloss Polishing

Sof-Lex™ Strips:

Different grits on each end with uncoated center for easy interproximal insertion


Gross Reduction & Contouring



C

M

Finishing & Polishing



F

SF

Sof-Lex™ Finishing and Polishing System







Sof-Lex™ Discs



Sof-Lex™ Diamond
Polishing System



Sof-Lex™ Strips



3M.com/SofLex

3M

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3m.com/SofLex

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, ESPE and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2015. All rights reserved. 70-2013-0721-5