

3M Science.
Applied to Life.™

Esthetic Restoration Solution

**A beautiful smile
she can't wait to share.**



Because you want every patient to leave with a beautiful smile.

Patients love the way highly polished restorations look and feel. But creating a high-end, natural-looking gloss can take more of your valuable time. So what do you do ... sacrifice higher esthetics for higher productivity? Or spend the extra time to create a beautiful smile ... because, "that's just the way it is?"

The days of time-consuming esthetic restorations are gone. With our Esthetic Restoration Solution, you'll love how simple it is to create natural-looking restorations with a gorgeous paste-like gloss.

Oh, don't be surprised if word of your great work gets around ... because a beautiful smile is worth sharing.



The truth is, the composite you use matters. You can spend a lot of time polishing most composites, only to find the initial gloss doesn't last long.

Or, you can create smiles your patients will love by using our Esthetic Restoration Solution. Thanks to 3M's nanotechnology, Filtek™ Ultimate Universal Restorative is easy to polish and also offers unsurpassed polish retention.

Plus, nanotechnology makes Filtek Ultimate restorative beautifully strong, giving you the option of using it in any place and for any Class.

Having one restorative that can do it all is a beautiful thing.



Class I



Class II



Class III



Class IV



Class V

Nanotechnology. For beautifully strong restorations ... anywhere.

Diamond particles. For a gorgeous, paste-like gloss.

The esthetic appearance of a restoration depends on the shade, shape and the gloss achieved by final polishing.

With the 3M™ ESPE™ Sof-Lex™ Diamond Polishing System, you can get a gorgeous, diamond paste-like gloss—with the convenience of a rubberized system!

The Sof-Lex™ Pre-Polishing Spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing. The Sof-Lex™ Diamond Polishing Spiral (pink) delivers the final high-gloss polish.

You'll find it reassuring to know that you no longer need to stop and switch between traditional points, cups, discs and brushes. The flexible spirals adapt to all tooth surfaces—anterior and posterior, convex and concave. They can be sterilized and reused, making them an economical choice for your practice.

A difference you can see!



Filtek™ Ultimate
Universal Restorative polished
with Sof-Lex™ Diamond
Polishing System



TPH Spectra® Universal
Composite polished with
Enhance® Finishing System
and PoGo® Polishing System

Notice a clearer reflection when
Filtek™ Ultimate Universal Restorative is
polished with the Sof-Lex™ Diamond Polishing System.

Now simple esthetics can be as easy as 1, 2, 3 ...

A factor in the success of this solution is 3M's unprecedented history in abrasive technologies. We leveraged that expertise to create a polishing solution that can help you maximize the esthetic benefits of Filtek™ Ultimate Universal Restorative with ease and simplicity.

Filtek™ Ultimate
Universal Restorative



Sof-Lex™
Pre-Polishing Spiral



Sof-Lex™ Diamond
Polishing Spiral



Ordering Information



Simple Esthetics

While Filtek™ Ultimate Universal Restorative is known for its wide range of shades and opacities, most single-shade restorations are made with one of 8 Body shades.

A1B A2B A3B A3.5B

B1B B2B C1B C2B

Filtek™ Ultimate Universal Restorative

Item # Product Information

3920S-R Filtek™ Ultimate Universal Restorative Single Shade Intro Kit—Syringe

Kit includes: 4 – 4g syringes: 2 each of Body shades A2 and A3; 1 – 5ml Single Bond Universal Adhesive Vial; 1 – 3ml Scotchbond™ Universal Etchant Syringe; Technique Guide; Instructions for Use

3920D Filtek™ Ultimate Universal Restorative Dual Shade Intro Kit—Syringe

Kit includes: 4 – 4g syringes: 1 each of Dentin shades A3, A4 and Enamel shades A2, A3; 1 Shade Selector Wheel; Technique Guide; Instructions for Use

DENTIN SHADE REFILLS

Refill	Item #
SHADE	SYRINGE
A1D	3920A1D
A2D	3920A2D
A3D	3920A3D
A4D	3920A4D
B3D	3920B3D
C4D	3920C4D
White	3920WD

ENAMEL SHADE REFILLS

Refill	Item #
SHADE	SYRINGE
A1E	3920A1E
A2E	3920A2E
A3E	3920A3E
B1E	3920B1E
B2E	3920B2E
D2E	3920D2E
White	3920WE
Extra White	3920XWE

BODY SHADE REFILLS

Refill	Item #
SHADE	SYRINGE
A1B	3920A1B
A2B	3920A2B
A3B	3920A3B
A3.5B	3920A3.5B
B1B	3920B1B
B2B	3920B2B
C1B	3920C1B
C2B	3920C2B
A4B	3920A4B
A6B	3920A6B
B3B	3920B3B
B5B	3920B5B
C3B	3920C3B
D2B	3920D2B
D3B	3920D3B
White	3920WB
Extra White	3920XWB

TRANSLUCENT SHADE REFILLS

Refill	Item #
SHADE	SYRINGE
Clear	3920CT
Blue	3920BT
Grey	3920GT
Amber	3920AT

SYRINGE SHADE REFILLS: Include 1 – 4g syringe; Instructions for Use

Sof-Lex™ Diamond Polishing System

Item # Product Information

5092-1 Sof-Lex™ Diamond Polishing System—Introductory Kit

Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); Instructions for Use; Technique Card

5090 Sof-Lex™ Pre-Polishing Spirals—Refill

Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)

5091 Sof-Lex™ Diamond Polishing Spirals—Refill

Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)



www.3MESPE.com



3M ESPE Dental
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dental

3M Deutschland GmbH
Location Seefeld
3M ESPE • ESPE Platz
82229 Seefeld • Germany
Info3MESPE@mmm.com
www.3MESPE.com

3M, ESPE, Filtek and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Printed in USA. © 3M 2016. All rights reserved.